

# OCULUS

Pentacam® | Pentacam® HR

Anterior Segment Tomography



INSTRUCTIONS FOR USE



## Preface

The Pentacam® / Pentacam® HR (High Resolution) has been manufactured and tested according to strict quality criteria.

- The operating concept of the Oculus Pentacam® and Pentacam® HR is basically identical.
- Additional functions pertaining to the Pentacam® HR (high resolution) are indicated accordingly.

To ensure safe operation, it is essential that you use the device correctly. For this reason you should familiarise yourself thoroughly with the contents of this Instructions for Use before operating the device. In particular, pay attention to the safety instructions in the instructions for use.

The following user information is available for your device:

- **Instructions for Use** describes how to manage patient data, the default settings of the Pentacam® program and the measuring process. The design of the unit is described in detail in this document. The instructions for Use also contain all safety-related instructions of the Pentacam® / Pentacam® HR.
- **User Guide** contains information supplementing the description of the operating concept. All features of the examination and analysis software are described in the User Guide, along with detailed information about the Patient Data Management system.
- **Interpretation Guide Ophthalmology** serves as a support for interpreting measurement results and graphical representations generated with a model from the Pentacam® family.
- **Software Installation** describes how to install the Pentacam® / Pentacam® HR software and the associated drivers.
- **Manual Floating License Key**: information on the use of the Pentacam® / Pentacam® HR within networks.

Due to ongoing development, the figures shown may depict minor changes to the actual device delivered.

If you have any queries or would like additional information about your device, do not hesitate to call or send us an email or a fax. Our service team will gladly assist.

OCULUS Optikgeräte GmbH

Item number: 10028075 / EN

Revision: 05

Release: 26.03.2026



## Table of Contents

<b>1</b>	<b>Scope of Delivery .....</b>	<b>9</b>
<b>2</b>	<b>Safety.....</b>	<b>10</b>
2.1	Symbols.....	10
2.1.1	On the device / name plate .....	10
2.1.2	On the packaging .....	11
2.1.3	In this Manual.....	11
2.2	Safety Instructions for Use.....	12
2.2.1	Instructions for the operation of a ME system.....	12
2.2.2	Instructions on Electrical Safety .....	13
2.3	Instructions on Cyber Security .....	14
2.3.1	Precautionary Measures Against Unauthorized Access....	14
2.3.2	Precautionary Measures if the Computer is Connect- ed to a Local or Internet Network.....	15
2.3.3	Device Security.....	15
2.3.4	Data Responsibility.....	16
2.3.5	Reporting and Handling of Security Incidents .....	16
<b>3</b>	<b>Device Description .....</b>	<b>17</b>
3.1	Overview of the device components .....	17
3.2	Applied parts.....	18
3.3	Functionality of the Pentacam® und Pentacam® HR.....	18
3.4	Intended Purpose.....	19
3.4.1	Intended Use .....	19
3.4.2	Intended Medical Indication.....	19
3.4.3	Contraindication.....	19
3.4.4	Possible side effects.....	19
3.4.5	Intended users .....	20
3.4.6	Patient group.....	20
<b>4</b>	<b>Setup and Installation.....</b>	<b>21</b>
4.1	Installation and Operating Conditions.....	21
4.2	Patient environment information.....	21
4.3	Setup .....	22
4.4	Electrical Connection.....	23
4.5	Switching On .....	24
4.6	Switching Off.....	24
4.7	Update Software or Install on separate PCs.....	24
<b>5</b>	<b>Pentacam® Program .....</b>	<b>26</b>
5.1	Menu Bar in the Pentacam® Programm.....	26
5.2	Scan Screen.....	27
5.2.1	Scheimpflug image settings .....	28
5.2.2	Parameters applicable only to the Pentacam® HR.....	28
5.3	Loading previous examinations.....	29

5.4	Online help.....	29
5.5	Information for recording Scheimpflug images.....	29
<b>6</b>	<b>Measuring Procedure.....</b>	<b>31</b>
6.1	Default settings .....	32
6.2	Darken the room.....	32
6.3	Rough adjustment.....	32
6.4	Fine Adjustment.....	35
6.5	Quality Specifications.....	36
6.6	Measuring Procedure for Tomography .....	40
6.6.1	Quality Specifications for the Tomography.....	42
6.7	CSP measurement (Pentacam® only).....	44
6.8	CSP Pro measurement (Pentacam® HR only).....	47
6.8.1	Quality specification for CSP Pro Measurement.....	49
6.8.2	Repeat or delete measurement.....	50
<b>7</b>	<b>Patient Data Management .....</b>	<b>53</b>
7.1	Starting Patient Data Management.....	53
7.2	Entering a New Patient.....	54
7.3	Selecting an Existing Patient .....	54
7.4	Extended Patient Search: [Extended] Checkbox.....	54
7.5	Rename Patient Data .....	55
7.6	Exporting Patient Data .....	55
7.7	Importing Patient Data.....	57
7.8	Data Backup.....	58
7.8.1	Backup Data .....	58
7.8.2	Reconstructing Data .....	59
7.8.3	Automatic Backup.....	59
<b>8</b>	<b>Cleaning, Disinfection and Maintenance.....</b>	<b>60</b>
8.1	Intervals for Cleaning, Disinfection and Maintenance .....	60
8.2	Consumables.....	60
8.3	Cleaning .....	60
8.3.1	Cleaning the Housing.....	61
8.3.2	Clean chin rest and forehead rest .....	61
8.3.3	Cleaning the illuminated slit.....	61
8.4	Disinfection .....	61
8.5	Test measurements.....	62
8.6	Maintenance.....	62
8.7	Attaching Paper to the Chin Rest.....	63
<b>9</b>	<b>Troubleshooting.....</b>	<b>64</b>
<b>10</b>	<b>Transport and Storage.....</b>	<b>65</b>
10.1	Storage conditions.....	65
10.2	Transport conditions .....	65
10.3	Disassembly .....	66

<b>11 Disposal .....</b>	<b>67</b>
<b>12 Terms of Warranty and Servicing.....</b>	<b>68</b>
<b>13 Technical Data .....</b>	<b>69</b>
13.1 Description of the Connection .....	76
13.2 Data sheet HEMG 49-S240210-7 (05150150).....	77
13.3 Instructions for integration into an IT-Network .....	79
<b>Annex</b>	
Annex A) Electromagnetic Compatibility .....	63
Annex B) Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity .....	64
Annex C) Description of the Connection .....	68
Annex D) Data sheet HEMG 49-S240210-7 (05150150).....	69
Annex E) Instructions for integration into an IT-Network.....	71



# 1 Scope of Delivery

Product and accessories
<b>Pentacam®   Pentacam® HR</b>
<ul style="list-style-type: none"> <li>■ x-y base</li> <li>■ Cograil</li> <li>■ Cover</li> <li>■ Support plate</li> <li>■ Sliding plate with chin rest paper</li> <li>■ Head and chin rest</li> <li>■ User information</li> <li>■ Pentacam® Basic Software</li> <li>■ Y-cable with Med. secure isolator + USB connection</li> <li>■ Country specific power supply</li> <li>■ Dark sheet with washing manual</li> </ul>
Optional software packages
<ul style="list-style-type: none"> <li>■ Screening Package</li> <li>■ Refractive Package</li> <li>■ Cataract Package</li> <li>■ Contact Lens Fitting Package with CSP Pro</li> </ul>
Single Licenses
<ul style="list-style-type: none"> <li>■ Holladay Report and Holladay EKR Detail Report</li> <li>■ 3D pIOL Simulation and Aging Prediction (only Pentacam® HR)</li> <li>■ IOL Calculator (only Pentacam® HR)</li> <li>■ DICOM</li> </ul>

We reserve the right to make changes to the scope of delivery in the course of technical development.

- ➔ If you find transport damage upon delivery, immediately file a claim with the transport company.
- ➔ Have the damages noted on the bill of lading, so that your claim for damages can be handled properly.

For more information regarding shipping and handling, see → [chap. 10 Transport and Storage, page 58](#).

## 2 Safety

All safety-related instructions of the device are given in the Instructions for Use only.

- ➔ Carefully read through the Instructions for Use.
- ➔ Keep the Instructions for Use in good condition near the device.
- ➔ Observe the legal regulations with regard to accident prevention.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer ([vigilance@oculus.de](mailto:vigilance@oculus.de)) and the competent authority of the region in which the user and/or patient is established.

### 2.1 Symbols

#### 2.1.1 On the device / name plate




Fig. 2-1: Name plate Pentacam® (example)











Fig. 2-2: Name plate Pentacam® HR (example)

No.	Description	No.	Description
1	Name and address of manufacturer	10	Applied part Type B
2	CE conformity and number of the notified body	11	Medical device
3	Date of manufacture	12	Power supply
4	Follow instruction for use	13	Model number
5	UDI-PI (Production Identifier)	14	Article number
6	machine-readable matrix code	15	Device type
7	UDI-DI (Device Identifier)		
8	Safety class		
9	Disposal in household trash is prohibited		

### 2.1.2 On the packaging

Symbol	Description
	Keep dry
	This way up
	Fragile
<b>Transport</b> 	Limit of temperature for transport
<b>Storage</b> 	Limit of temperature for storage
	Limit of humidity
	Limit of air pressure

### 2.1.3 In this Manual



#### Warning

Identifies a potentially dangerous situation which may cause serious injury.



#### Caution

Identifies a potentially dangerous situation which may cause minor injury.



#### Attention

Identifies situations which may cause damage to the device or incorrect test results.



Indicates important instructions for use and important information about the device.

- > This symbol denotes menu paths and screen shots. Example for starting a new examination:  
Pentacam<sup>®</sup> und Pentacam<sup>®</sup> HR > Examination > Scan  
which means:
  - ➔ Select the "Examination" menu from the menu bar.
  - ➔ Select the menu item "Scan".
- [ ] Marks buttons.
- Cross reference

## 2.2 Safety Instructions for Use



### Caution

Personal injury or property damage due to improper operation

- ➔ Observe the following safety instructions.



### Caution

Personal injury or property damage due to equipment modifications that could jeopardize safety

- ➔ No modifications may be made to this device without the permission of the manufacturer. Only the OCULUS service and authorized dealers are allowed
  - to modify the device or the associated lifting table
  - to install software and software updates.

### 2.2.1 Instructions for the operation of a ME system

The device and a connected computer form together a medical electrical system (ME system) according to IEC 60601-1. If you connect additional devices, such as, for example a printer, those devices become part of the ME system.

- ➔ All devices of the ME system must comply with the requirements of IEC 60601-1 or IEC 62368-1.

## 2.2.2 Instructions on Electrical Safety


**Caution**
**Risk of personal injury or damage to property due to an incorrect level of safety**

Connecting the Pentacam® und Pentacm® HR with its non-medical electrical equipment (e.g. data processing equipment) to a medical electrical system must not result in a patient safety level below that prescribed by IEC 60601-1. If making this connection leads to the leakage current threshold being exceeded, protective measures including a circuit breaker must be in place.

- Ensure that connections with non-medical devices are made correctly.
- Only use the power adapter listed in the packing list.
- Use only a computer that meets the specifications given in this Instructions for Use, → [chap. 12 Technical Data, page 61](#).


**Caution**
**Risk of personal injury or material damage caused by unsafe multiple socket extension cord**

If you use a multiple socket extension cord to connect the device to the power supply, you must heed the following information:

- Use an extension cord that complies with the requirements of IEC 60601-1, section 16.
- Do not place the multiple socket extension cord on the floor.
- Do not use more than one multiple socket extension cord.
- Plug only the device and the computer that is being used with the unit (if applicable) into the multiple socket extension cord.
- If you are using a multiple socket extension chord it has to be supplied with a isolation transformer.
- If you are using a new computer for the device, you must have the electrical safety checked. Call OCULUS Service for this purpose.


**Caution**
**Risk of personal injury or damage to property due to electromagnetic interference**

Portable and mobile RF communications equipment can affect medical electrical equipment see chap. "Annex A Electromagnetic Compatibility" on page 63.

- Make sure that portable and mobile RF communications equipment do not cause interference.
- Recommendation: Maintain a minimum distance (refer "Annex A Electromagnetic Compatibility" on page 63). If the distance is shorter, you must ensure that the Pentacam® und Pentacm® HR functions correctly.

## 2.3 Instructions on Cyber Security



### Note

Please also follow the regulations, guidelines, and recommendations of the relevant authorities responsible for information security and the protection of critical infrastructure in your country.



The device is designed to operate independently without requiring an internet connection, network access, or portable media. It functions solely through a connected computer.

If you choose to connect the computer to the internet or another network for unrelated purposes, you are responsible for ensuring that the connection is secure and controlled.

### 2.3.1 Precautionary Measures Against Unauthorized Access

To increase the cybersecurity of the device:

➔ Secure the device against unauthorized access by unauthorized persons.

Observe all precautionary measures:

- Secure the computer with a strong password (e.g. during Windows startup).
- Choose a complex password of at least twelve characters that includes letters, numbers and special characters. Avoid dictionary words.
- Do not select a name or a device name as a password (such as "Pentacam").
- Change the default password after the first authentication.
- Change the password regularly.
- Do not write down the password at an accessible location.
- Use unique passwords for different user accounts.
- Do not share user names or passwords with colleagues or anyone else, even if authorized by law and employer policy to view the same type of information (e.g. two users reviewing the same patient samples).
- Enable a screen saver that requires re-entry of the password upon deactivation.
- Set an appropriate screen saver timeout (e.g. 10 minutes) based on operational conditions such as examination duration and patient flow.
- Ensure that the device is locked (keyboard shortcut: Windows logo button + 'L') or otherwise secured when not in use to prevent unauthorized access to electronic Protected Health Information (ePHI).
- Operators are required to be trained according to privacy awareness and handling personal data.
- If necessary, contact your Healthcare Organization's IT Department.

### 2.3.2 Precautionary Measures if the Computer is Connected to a Local or Internet Network

- Do not establish an internet connection while using the device. This is considered misuse!
- If you connect the computer to the internet for any other purpose, you are responsible for ensuring data security.

If you connect your computer to a local network, you are responsible for ensuring data security and must observe at least the following precautionary measures:

- Preferably connect the computer to the network using a cable connection over wireless connection.
- Use robust security methods incl. advanced Encryption Standards with a strong network key also for wired connections (not only for Wi-Fi connections).
- Use of a firewall (software or hardware) is recommended.
- Observe the instructions for integration into an IT network refer "Annex E Instructions for integration into an IT-Network" on page 71.



#### Note

The Healthcare Organization's IT Department should implement a risk management framework in alignment with IEC 80001-1 to support the secure and safe integration of medical IT networks. This includes assessing risks, enforcing access control, securing networks, applying software updates, monitoring incidents, protecting data, managing device life cycles, and training staff to help safeguard patient safety and data integrity.

The Manufacturer Disclosure Statement for Medical Device Security (MDS2) is available upon request for detailed security information.

---

### 2.3.3 Device Security

Ensure that the device is secured against unauthorized access refer "2.3.1 Precautionary Measures Against Unauthorized Access" on page 14.

- ➔ Protect the device and attached systems from malicious software.
- ➔ Implement new software versions when available.
- ➔ Implement operators access only on a need-to-know-basis.

The Healthcare Organization's IT Department is responsible for implementing controls for the handling and disposal of media and assets.

### 2.3.4 Data Responsibility

Operators should avoid entering unnecessary identifying data. Whenever possible, data should be de-identified and linked to the sample ID instead of the patient. Use only the input data essential for the intended purpose.

Operators have access to sensitive patient data (ePHI).

- Do not take any snapshots, screenshots, or images (e.g. using another device) of the information displayed on the device.

Data shall be deleted on a regular basis according to the Healthcare Organization's deletion policy, if respective data are processed by the device.

The Healthcare Organization's IT Department is responsible for deleting unused user accounts.

Only authorized personal is allowed to take backups. The Healthcare Organization's IT Department shall manage the location of each backup in order to respond to potential data subject requests. Backups and archive files are required to be transmitted and stored securely.

### 2.3.5 Reporting and Handling of Security Incidents

Operators must inform their Healthcare Organization's IT Department about any suspected or confirmed privacy or security breaches, including suspected or compromised user accounts, and report any service outages or access issues.

- If accounts are deemed compromised, devices are lost, or unauthorized access has been discovered or assumed, the Healthcare Organization's IT Department locks or changes the user login criteria and issues new login information so that the user can safely access his or her account.

## 3 Device Description

### 3.1 Overview of the device components




Fig. 3-1: Device components

No.	Description
1	Ventilation opening
2	Measuring window
3	Camera opening
4	On/off switch with indicator light
5	Y-cable connector
6	Cross slide
7	Sliding plate
8	Circular markings on sliding plate
9	Joystick
10	Nameplate

## 3.2 Applied parts

The chin and head rest are applied parts of type B.



Fig. 3-1: Anwendungsteile

No.	Description
1	Head rest
2	Chin rest

## 3.3 Functionality of the Pentacam® und Pentacam® HR

While rotating around the eye, the Pentacam® captures Scheimpflug images of the anterior eye segment through varying axes. The Scheimpflug images created during an examination are transmitted to the connected PC.

Up to 50 (Pentacam) or 100 (Pentacam HR) Scheimpflug images can be captured within maximum two seconds.

Up to 25,000 (HR: 138,000) genuine height values are measured and analysed from the Scheimpflug images.

The Scheimpflug images are the basis for the height data which are used to calculate a mathematical 3D model of the anterior eye segment. At same time any eye movements are recorded and taken into account.

A quality specification (QS) lets you gauge the quality of the measurement taken.

The mathematical 3D model, corrected for eye movements, provides the basis for all subsequent analysis.

The topography of the front and rear surfaces of the cornea and the pachymetry are calculated and displayed from limbus to limbus for the entire surface of the cornea.

An analysis of the anterior chamber provides the basis for calculating the chamber angle, chamber volume and chamber depth.

Densitometry of the cornea and lens automatically produces quantified values.

Colour images displayed on the screen show the results of the measurement.

A rotatable 3D model displays the front and rear surfaces of the cornea, the iris and the lens.

### 3.4 Intended Purpose

---



OCULUS Optikgeräte GmbH is not liable for the evaluation and interpretation of measurements taken with the Pentacam® / Pentacam® HR. The user manual and interpretation guide can provide assistance in this regard.

---



The device is intended exclusively for the use specified in this manual and in compliance with the safety instructions.

---

#### 3.4.1 Intended Use

The OCULUS Pentacam® und Pentacam® HR is designed to take photos of the anterior segment of the eye which includes the cornea, pupil, anterior chamber and lens of the eye to evaluate:

- corneal shape,
- analyse condition of the lens (opaque crystalline lens),
- analyse the anterior chamber angle,
- analyse anterior chamber depth,
- analyse the volume of the anterior chamber,
- analyse anterior or posterior cortical opacity,
- analyse the location of cataracts (nuclear, subcapsular and or cortical), using cross
- slit images with densitometry and
- corneal thickness.

#### 3.4.2 Intended Medical Indication

The Pentacam® und Pentacam® HR is indicated as an aid to screen several eye diseases for example, but not limited to:

- Keratoconus classification and progression
- Early ectatic diseases
- Quantification of optical corneal density
- Quantification of optical lens density
- Close angle glaucoma
- Planning support for IOLs.

#### 3.4.3 Contraindication

None known.

#### 3.4.4 Possible side effects

- After-image
- Headache
- Vertigo
- Tearing eyes

### 3.4.5 Intended users

The Pentacam® und Pentacam® HR is intended exclusively for use in:

- ophthalmology practices
- clinics
- opticians or optometrists

The Pentacam® / Pentacam® HR is intended for use by trained personnel:

- Who can ensure proper handling based on their knowledge, training and practical experience.
- Who have been instructed by OCULUS personnel or an authorized dealer prior to commissioning.

### 3.4.6 Patient group

- Children from 3 years up to not limited.
- No restrictions on weight.
- Patient is awake.
- Patient is able to understand and to look into a fixation target.

## 4 Setup and Installation

- Only OCULUS or an authorized dealer is allowed to set up and to connect the device.
- Do not use or store the device in rooms that are humid, see chap. "10 Transport and Storage" on page 58.
- Keep the device away from water that may drip, splash or spray on it, and make sure that no liquids can get into the device. Do not place any containers holding liquids in the vicinity of the Pentacam® und Pentacm® HR.
- Germany: Only operate the device in rooms used for medical purposes after they have been set up according to the VDE Regulation 0100-710.
- Do not operate the devices included in the delivery in areas where explosions may occur, or in proximity to flammable anesthetics or volatile substances such as alcohol, benzine or similar products.
- Set up the device so that the power plug is easy to access. That way, you can easily disconnect it from the power supply for any repairs or maintenance work.
- ➔ Do not force any plug connections.  
If you are unable to make a plug connection, check whether the plug fits the socket.  
If you detect damage to the connection, you should let authorized dealer repair the damage.
- ➔ Only use a device which is mounted at the lifting table properly.

### 4.1 Installation and Operating Conditions

Temperature	+10 – +35 °C
Air humidity	30 – 90 %
Air pressure	800 – 1060 hPa

- Before installation, compare the transport and storage temperature with the temperature in the intended installation room.
- The difference between the transport and storage temperature and the installation room should not be more than 10°C to prevent the internal optics from misting up.

### 4.2 Patient environment information

Patient environment is the area where patients can come into contact with any part of a medical electrical equipment (ME equipment) or with another person being in contact with the ME equipment.

- ➔ Only use devices that are compliant with IEC 60601-1 in the patient environment.




Fig. 4-1: Patient environment

### 4.3 Setup



#### Attention

Incorrect installation may result in incorrect measurements or even damage to the device.

- Have the device set up and connected by our service department or by a specialist authorized by OCULUS.
- Set up the device so that it cannot fall over. Mount the device on an examination table.
- Set up the device so that it is protected from dripping, splashing or spray water.



#### Attention

Damage to the device due to incorrect handling

- Do not expose the device to any vibrations, shocks, contaminants, moisture, or high temperatures.
- Handle the device with care.

- Set up the Pentacam® und Pentacm® HR so that the power plug is easy to access. That way, you can easily disconnect it from the power supply for any repairs or maintenance work.
- Place the device so that direct light cannot affect the measurement.
- Make sure the examination is free from light reflections. To achieve this, darken the examination room.

## 4.4 Electrical Connection

### Caution

If there are temperature differences between the transport and storage temperature and the installation room, especially if the temperature exceeds 10°C, the optics may fog up and/or condensation may form.

- Leave the device unused in the installation room for at least 6 hours before connecting it so that the device can adjust to the ambient temperature.

### Caution

Electrical safety hazard

- Do not use the device directly next to other devices.
- Do not stack the device with other devices.
- Only use the power supply unit supplied or one identical to the one specified in the Tech. Data (see chap. "12 Technical Data" on page 61).
- Do not place any heavy objects or the device itself on the power cable.
- If you use a multiple socket to connect the device, the multiple socket must comply with the requirements of IEC 60601-1.
- Do not place the power strip on the floor.
- Do not expose the mains cable or the multiple socket to high temperatures. Do not place on heaters!
- Use a maximum of one multiple socket.
- Only connect the device and, if applicable, the associated computer to this multiple socket.
- Use a socket outlet that has a faultless protective conductor connection.



No.	Description
1	On / Off switch
2	Plug of the Y-cable

Fig. 4-1: Electrical Connection



### Attention

If the device is not connected correctly and voltage is applied, the device may be damaged after a short time.

- Do not use excessive force to connect electrical connections.
- Observe the specifications on the name plate.
- If the plug is defective, contact OCULUS Service or an authorized dealer to repair the damage.

1. Plug the connector of the Y-cable into the jack and tighten the connection. Make sure that the plug is inserted in the correct position.
2. Firmly tighten the connection.
3. Connect the Y-cable to the PC/laptop and the power adapter.

## 4.5 Switching On

### Caution

Risk of incorrect measurements due to improper setup

- Before taking measurements, the device has to be switched on, at least for one hour.

1. The first step is to switch on the PC or laptop.
2. Then turn on the device with the on/off switch (position ON). The LED on the switch lights up green, [fig. 4-1, page 23](#).

## 4.6 Switching Off

1. Close the Pentacam program and close the Patient Data Management.
2. Shut down the Windows operating system.
3. Turn the device off with the on/off switch (OFF position)

## 4.7 Update Software or Install on separate PCs

The Pentacam-Software is network compatible. This makes it possible to install the Pentacam® Software on several PCs, connected in a local network.

Make sure, that all PC in a network do have the same Pentacam® software version installed.

A Floating License Key is part of every Pentacam® shipment. Please ensure a proper installation based on the Manual for the Floating License Key.



The software or software updates may only be installed by OCULUS Support or an authorized dealer.

The latest version of the software is already installed upon delivery.

This allows an interactive parallel evaluation of the Pentacam® exams based on the released optional software packages and modules.

The demo examinations included can be viewed on every computer having the Pentacam® software installed.

For more detailed information contact your authorized dealer or our service department.

## 5 Pentacam® Program

The device is safe to use if no error message is displayed after starting the software when the device is connected and switched on (e.g. component failure, camera not recognized, missing reference data, etc.).

After starting, the patient data management opens (see chap. "7 Patient Data Management" on page 45).



The Pentacam® program is not to be used to dictate therapy without professional correlation with other clinical findings and diagnostic tests.

### 5.1 Start screen

To start the Pentacam® program from the patient data management click on the Pentacam® button.




Fig. 5-1: Start screen: Overview

No.	Description
1	open 'Scan' screen
2	Examination date and Patient data

5.2 Scan Screen

→ menu bar [Examination] > [Scan]




Fig. 5-1: "Scan" screen (HR model)

No.	Description	
1	Front image / Pupil image	shows a live image of the patients eye and the vertical / horizontal positioning of the Pentacam The object is to position the yellow dot marking the apex of the cornea in the centre of the crosshairs. The blue dot marks the centre of the pupil, while the blue ring frames the pupil.
2	Eye	The eye currently being examined is usually detected automatically. If not, you can select it manually here.
3	Scheimpflug Image	If this option is enabled, only a single Scheimpflug image will be captured. You can freely select the desired camera position by clicking on the white rings in the "Orientation" field.
4	Enhanced Dynamic Scheimpflug Image	Use this option to make the camera record either 5, 10 or 15 Scheimpflug images, with the camera remaining in the same position. Image averaging is carried out on the recorded images to minimise background noise. Only one Scheimpflug image is displayed. You can freely select the camera position you require by clicking the white rings in the "Orientation" field. This type of image is suitable for a purely densitometric assessment of the lens.

No.	Description	Description
5	Cornea Fine	(only Pentacam HR) 3D scan with 50 or 100 cross-sectional images (instead of the standard 25 cross-sectional images) This option allows a more detailed image of the cornea. The deeper layers of the anterior segment are not captured in this process. You can choose between 50 Scheimpflug images in 1 second or 100 Scheimpflug images in 2 seconds.
6	3D Scan	Select the number of images to be captured per scan. The difference lies in the duration of the examination and the number of measurement points analyzed. A scan with 50 images therefore takes longer, but achieves the highest accuracy when the patient is well-positioned. This type of examination is selected for evaluating the cornea and the anterior chamber.
7	[Slit Light]	Activate or deactivate the illumination of the eye with blue light.
8	[Automatic Release]	Activate or deactivate automatic measurement.
9	[Start Scan]	Activates manual measurement, when [Automatic Release] is deactivated. You can also use the Return key.
10	Adjustment window	The arrows show the direction in which you must move the device to activate automatic measurement (Automatic Release).
11	Fixation target	(only Pentacam® HR) appears as a red LED in the center of the blue slit light The "fixation target" helps the patient fixate more effectively. To do this, the active "fixation target" can be adjusted in 0.5 dpt increments. The goal is to compensate for the patient's refractive error and ensure easier fixation.
12	[Quit] button	cancel current measurement.
13	Device Notifications	appear as needed, e.g., when a service is due
14	Scheimpflug Image	shows the distance between the device and the patient. The objective of this setting is to move the red dot on the front surface of the cornea so as to coincide with the red line
12	"Orientation" field	shows the respective position of the camera and the eye, which is currently being examined.


### 5.3 Loading existing examinations

1. Select the menu item [Examination] and click [Load].  
The dialog box "Load Examination" appears.
2. Make a selection by clicking the required examination.
3. Confirm your selection by clicking [OK], or by double clicking.  
The Pentacam® program will load the examination you have selected.

## 5.4 Online help

Direct help can be accessed by clicking on a small yellow mark in the corners of buttons and fields.

Examples:



## 5.5 Recording Scheimpflug images for specific diagnostic purposes

Type of examination	Examination mode	Images	Automatic measurement	Notes
Topography	3D scan	25-50	Yes	
Pachymetry	3D scan	25-50	Yes	
Analysis of the anterior chamber	3D scan	25-50	Yes	Do not apply mydriatic drops.
Artificial lenses (general)	Enhanced dynamic (for HR in the "3D scan" examination mode)	15	Yes	If the pupil is insufficiently dilated, apply mydriatic drops. Use 3D scan for measurements.
Measuring functions	3D scan	25-50	Yes	If the pupil is insufficiently dilated, apply mydriatic drops.
Densitometry	3D scan Enhanced dynamic	25-50 5-15	No	Use the same number of images to enable a progress check and apply mydriatic drops.

### Specific instructions for the Pentacam® HR

Type of examination	Examination mode	Images	Automatic measurement	Notes
IOLs, ICLs, PIOLs	3D scan for PIOLs, possibly longer exposure	25-50	Yes	If the pupil is insufficiently dilated, apply mydriatic drops

## 6 Measuring Procedure



### Note

Risk of incorrect measurement due to incorrect use

- Before first use: Let OCULUS or an authorized dealer train you in the operation of the Pentacam® und Pentacam® HR.



### Note

Risk of incorrect measurements due to improper setup

- Before taking measurements, the Pentacam® und Pentacam® HR has to be switched on, at least for one hour.



### Note

Due to small movements of the patient or caused by a wheelchair, the patient is no longer positioned appropriately to the Pentacam® / Pentacam® HR, which can lead to incorrect measurements.

- Perform a Pentacam® und Pentacam® HR scan only if the patient sits in a stationary chair.
- In cases of wheel chairs lock the brakes.



The instructions for use focus on the operating concept of the Pentacam® / Pentacam® HR. The functional description of the Pentacam program is limited to initiating a measurement and loading existing examinations.

For detailed information about the features of the Pentacam program, see the user manual.

- Before first use: Let OCULUS or an authorized dealer train you in the operation of the device.
- Never operate a damaged device.
- Only operate the device with the original accessories supplied by us and only when the unit is in technically perfect condition.
- Do not cover the ventilation openings.
- Do not touch the patient and the device at the same time.
- Make sure that the device cannot tip over by leaning against it or sitting on it.
- Do not put the device, including rechargeable battery or cable, down onto devices that produce heat, heaters (for example radiators), microwaves or similar.
- Only operate the device if you have understood the operating instructions.



### ISO 15004-2:2007 Group 2 instrument

A graph showing the relative spectral output of the instrument between 305 nm and 1100 nm when the instrument is operating at maximum light intensity and maximum aperture will be provided by OCULUS on request. The spectral output will be shown for the beam after it exits the instrument.



### Caution

The light emitted from this instrument is potentially hazardous. The longer the duration of exposure and the greater the number of pulses, the greater the risk of ocular damage. Exposure to light from this instrument when operated at maximum output will exceed the safety guideline after 225 treatments.

## 6.1 Preparations

1. Start the Scan menu in the Pentacam® Program:  
Menu [Examination] > [Scan]  
The blue slit light is activated and the Scan menu opens.
  2. If necessary make changes to the image options for the particular part of the front of the eye that is to be examined.  
The default settings in the "3D Scan" options are "25 images/1 second".
  3. Adjust the table height.
  4. Check that
    - fresh paper has been put onto the chin rest or
    - the chin rest has been cleaned and disinfected.
  5. Ask the patient to place his or her head on the chin and forehead rest.
- Do not touch the patient and the Pentacam® und Pentacm® HR simultaneously.

## 6.2 Darken the room

6. Dim the lights in the room.
- If the lighting in the examination room cannot not be dimmed or switched off, use the dark sheet supplied to cover the patient and the device.



Fig. 6-1: Patient and device covered with dark sheet

### 6.3 Rough adjustment

7. Position the patient correctly.  
Adjust the table height if necessary.  
The patient is seated correctly when
  - the chin rests in the indentation of the chin rest
  - the forehead rests against the forehead rest
  - the eyes are level with the mark.





Fig. 6-1: Rough adjustment of the chin and forehead rest

No.	Description
1	Black ring for eye level orientation
2	Twist grip for adjusting the height of the chin rest

6.4 Align device



No.	Description
1	Cross slide
2	Left marking
3	Right marking
4	Marking on the cross slide
5	Joystick
6	Sliding plate

Fig. 6-1: Markings on the cross slide

8. Move the cross slide (example for the right eye) until the marking at the end of the cross slide roughly coincides with the circle R on the sliding plate.
9. Look at the patient's eye you are examining from one side and make sure that the blue slit light illuminates the cornea.  
If necessary, adjust the position of the cross slide to the left or right.




Fig. 6-2: Slit light on the cornea



If the blue slit light is not visible, ensure that you have activated the [Slit Light] checkbox on the "Scan" screen.

## 6.5 Fine Adjustment and start the measurement

- The settings for the desired measurement have been configured on the scan screen → [chap. 5.2 Scan Screen, page 27](#).
- 10.** Move the cross slide towards the patient until the Scheimpflug image shows the cornea of the eye that you are examining.





Fig. 6-1: Scheimpflug image (left) and pupil image (right)

- 11.** Focus the pupil image by moving the joystick towards the Pentacam® or away from it. The image is sharpest when the red dot coincides with the red line in the Scheimpflug image.
- 12.** Adjust the left/right position of the Pentacam® and its height setting. Move the joystick to the left or right and rotate the joystick clockwise or anti-clockwise. The tentative final position of the camera is reached when the yellow dot is in the centre of the crosshairs.
- 13.** The patient should blink once more, then open their eyes wide and stop blinking.
- 14.** Make any final adjustments to the device's alignment.

Once the position is sufficiently accurate, a cross surrounded by four bars will appear in the center of the ring.

- The device will automatically trigger the measurement.

To trigger the measurement manually, press the [Scan] button or the Return key.





A measurement taken with manual triggering may not be reproducible.

- 15. After the measurement, the patient can remove their head from the chin and forehead rests.  
After the measurement, the overview screen of the Pentacam® program opens automatically.
- 16. Check the measurement results by referring to the quality specifications ( → chap. 8 Cleaning, Disinfection and Maintenance, page 52).

## 6.6 Checking the quality (QS) of a measurement and identifying measurement errors

The analysis using the “QS” field helps evaluate the quality of the measurement and identify errors in the measurement process



### Note

All examinations are automatically saved, regardless of the quality of the measurement taken.

QS field	Meaning
OK	Measurement is correct and can be reproduced.
red	Repeat measurement.
yellow	Measurement not optimal: Check measurement results! Click on field.



If the “QS” field is highlighted in yellow or red, check the QS values.

Clicking the “QS” field opens the detailed quality specifications.




Fig. 6-1: Pentacam® program with "QS" field (1)

17. Check the measurement results.  
If you have any doubts, repeat the measurement.





Fig. 6-2: Examination Quality Specification

QS Parameter	If an error occurs here, ...
Analysed Area	... the measured area of the cornea is too small. ➔ The patient must open their eye wider. If necessary, assist the patient by holding their upper eyelid.

QS Parameter	If an error occurs here, ...
Valid Data	<p>... no continuous data points could be determined in the Scheimpflug images, which can sometimes be the case with irregular or very cloudy corneas. In such cases, even repeated measurements will yield values highlighted in yellow at most. Or an ambient light source interfered with the recording.</p> <p>→ Completely darken the room..</p>
Lost Segments Lost Seg. Continuous	<p>... the patient blinked or the shadow cast by the nose is too large.</p> <p>→ Before the measurement process begins, the patient should blink again and then, without blinking, focus on the red LED or black ring on the device during the measurement process.</p> <p>→ If the error was caused by the nose covering the camera line, you must slightly turn the patient's head so that the nose is positioned away from the camera..</p>
Alignment (XY) and Alignment (Z)	<p>... the device was moved during the measurement trigger..</p> <p>→ Repeat the measurement.</p>
Eye Movement	<p>... the patient has not focused correctly on the target.</p> <p>→ Before the measurement process begins, the patient should blink once more and then, without blinking, focus on the red LED or black ring on the device during the measurement process.</p>
CSP Fixation	<p>→ If this value exceeds the limit, the measurement must be repeated. If necessary, explain to the patient that they must focus on the black ring.</p>

## 6.7 End measurement

18. Close the window by clicking [Cancel].
19. End the current measurement or prepare to measure a new patient via the menu [Patient] > [New Patient/End].

The Pentacam® program will close. You will return to the patient data management screen, where you can create or select a new patient.

## 6.8 Measuring Tomography

A tomography measurement using the Tomography Sequence menu described below is only possible with a CSP Pro license.

1. To open the tomography measurement in the Pentacam® program:  
Menu [Examination] > [New Tomography Sequence]  
The blue slit light is activated, and the scan screen for the tomography measurement opens.
2. Prepare the measurement, adjust the patient and perform the measurement → chap. 6.1 Preparations, page 31 to → chap. 6.5 Fine Adjustment and start the measurement, page 34




Fig. 6-1: Scan screen "Tomography"

No.	Designation	Description
1	[Automatic Release]	If the checkbox is selected, the measurement will trigger automatically as soon as the position is sufficiently accurate
2	[Start Scan]	
3	[Zurück]	
4	[CSP Pro]	Enable/Disable CSP Pro Measurement → chap. 6.9 CSP Pro measurement, page 40
5	Live Scheimpflug image	
6	Tomography parameter	<ul style="list-style-type: none"> <li>■ <b>Use:</b> Examination can be used for the evaluation. Click the checkbox of the respective measurement to use it for the evaluation.</li> <li>■ Use only one measurement for the full sequence examination.</li> <li>■ <b>QS:</b> Quality specifications, see → Chap. 6.6 (Page 35)</li> <li>■ <b>K1:</b> Flat radius of the corneal curvature</li> <li>■ <b>K2:</b> Steep radius of the corneal curvature</li> <li>■ <b>CCT:</b> Central corneal thickness</li> <li>■ <b>ACD:</b> Anterior chamber depth</li> </ul>
7	Current examination mode	= Tomography

- Check the measurement results by referring to the quality specifications ( → chap. 6.6 Checking the quality (QS) of a measurement and identifying measurement errors, page 35).

After the scan, the overview image opens.




Fig. 6-2: Result display of tomopgraphy scan

No.	Description
1	"QS" value for quality specification
2	[+New Scan] button
3	[Use] checkbox



**Note**

All examinations are automatically saved, regardless of the quality of the measurement taken.

QS field	Meaning
OK	Measurement is correct and can be reproduced.
red	Repeat measurement.
yellow	Measurement not optimal: Check measurement results! Click on field.



If the "QS" display is highlighted in yellow or red, check the QS values → chap. 6.6 Checking the quality (QS) of a measurement and identifying measurement errors, page 35

## 6.9 CSP Pro measurement

In a CSP Pro measurement, not only the cornea, but also parts of the sclera are measured. This means larger diameter contact lenses such as scleral lenses can be fitted.

By default, always a tomography measurement is carried out ( → chap. 6.8 Measuring Tomography, page 38). To carry out a CSP Pro measurement, proceed as follows:

1. .
2. Make sure that the [Automatic Release] checkbox (1) is selected.
3. Prepare the measurement, adjust the patient and perform the measurement → chap. 6.1 Preparations, page 31 to → chap. 6.5 Fine Adjustment and start the measurement, page 34.



Fig. 6-1: "CSP Pro" examination screen

No.	Designation	Description
1	[CSP Pro]	swipe right to enable CSP Pro measurement
2	Parameter CSP Pro	<ul style="list-style-type: none"> <li>■ <b>Use:</b> Examination can be used for the evaluation. Click the checkbox of the respective measurement to use it for the evaluation.</li> <li>■ Use only one measurement for the full sequence examination.</li> <li>■ <b>QS:</b> Quality specifications, see → Chap. 6.6 (Page 35)</li> <li>■ <b>K1:</b> Flat radius of the corneal curvature</li> <li>■ <b>K2:</b> Steep radius of the corneal curvature</li> <li>■ <b>CCT:</b> Central corneal thickness</li> <li>■ <b>ACD:</b> Anterior chamber depth</li> </ul>
3	Current examination mode	= CSP Pro

4. Check the measurement results by referring to the quality specifications ( → chap. 6.6 Checking the quality (QS) of a measurement and identifying measurement errors, page 35).

After the scan, the overview image opens.




Fig. 6-2: Results display for a CSP Pro measurement

No.	Description
1	[+New Scan] button
2	"QS" value
3	Value for the overall coverage
4	Display sum of all individual measurements
5	Individual measurements
6	[Add measurement] button
7	[Repeat measurement] button



**Note**

All examinations are automatically saved, regardless of the quality of the measurement taken.

QS field	Meaning
OK	Measurement is correct and can be reproduced.
red	Repeat measurement.
yellow	Measurement not optimal: Check measurement results! Click on field.



If the "QS" display is highlighted in yellow or red, check the QS values → [chap. 6.6 Checking the quality \(QS\) of a measurement and identifying measurement errors, page 35](#)




Fig. 6-3: Error message after a CSP Pro measurement

No.	Description
1	Error message
2	Button [X] = Delete measurement
3	Button [+] = Add measurement

If the QS value of the first measurement is not sufficient to use this for evaluation, repeat the measurement:

1. Click on the [Repeat measurement] button.  
The measurement is repeated and the previous measured values automatically get deleted.



The eye needs to be held open wide enough for the desired measurement area not to be covered by the eyelids to achieve good coverage of the cornea and sclera. We recommend holding up the upper eyelid with the LidStick® or alternatively a long cotton swab. The patient can carefully pull down his or her lower eyelid with a finger.

2. Click the [Add measurement] button to the right of the last measurement taken to add a measurement.  
The coverage area of each measurement is shown in a different color in the corresponding display.
3. Carry out additional examinations until the required measurement area is achieved i.e. a complete corneal scleral profile.

4. If necessary, delete measurements with yellow or red QS. This is also necessary if more than 4 individual measurements have to be carried out in order to obtain coverage > 95%.




Fig. 6-4: Coverages of individual measurements and total coverage

No.	Description
1	Coverage areas of the single measurements Scan 1 and Scan 2 superimposed
2	Single measurement scan 1
3	Single measurement scan 2
4	[Finish] button

5. Click the [Finish] button to complete the CSP Pro measurement.




Fig. 6-5: Completed CSP Pro measurement

No.	Description
1	Parameters of the CSP Pro measurement
2	Scheimpflug images
3	Elevation map

## 7 Patient Data Management

Once you have completed an examination, you can do the following with the patient data:

- Rename it  
(chap. "7.5 Rename Patient Data" on page 47)
- Export it  
(chap. "7.6 Exporting Patient Data" on page 47)
- Import it  
(chap. "7.7 Importing Patient Data" on page 49)
- Backup  
(chap. "7.8 Data Backup" on page 50)



For more information on Patient Data Management, refer to the [User Guide](#).

### 7.1 Starting Patient Data Management

After you have switched on the PC, it first loads the operating system.

1. Press on the Pentacam® icon at the desktop.

The user interface for the Patient Data Management appears:





Fig. 7-1: Patient Data Management user interface

Nr.	Designation	Description
1	"Functions" area	Auflistung aller zur Verfügung stehenden Geräte <ul style="list-style-type: none"> <li>■ Einstellungen</li> <li>■ Backup: chap. "7.8 Data Backup" on page 50</li> <li>■ Beenden</li> </ul>
2	Previous examinations	all existing test results for the selected patient
3	[Delete exam.] button	
4	[Delete Patient] button	
5	[Export] button	chap. "7.6 Exporting Patient Data" on page 47
6	[Import] button	chap. "7.7 Importing Patient Data" on page 49
7	Patient list	
8	„Patient“ area	Patient data <ul style="list-style-type: none"> <li>■ New: chap. "7.2 Entering a New Patient" on page 46</li> <li>■ Change</li> <li>■ Search</li> </ul>

## 7.2 Entering a New Patient

1. Press the [New] button to enter a new patient in the Patient Data Management system.
2. Enter the patient's last name, first name and date of birth in the patient window.



The screenshot shows a window titled "Patient:" with the following fields and buttons:

- Last name:
- First name:
- Dat. o. B.:
- ID Number:
- Comment 1:
- Comment 2:
- Save button
- Cancel button

Fig. 7-1: Entering patients

- Optionally you can enter an ID number for the patient.
3. To save your entries use the [Save] button.  
The patient you have just entered now appears in the patient list and is automatically selected..

## 7.3 Selecting an Existing Patient

The patient data list on the left-hand side of the screen displays all previously examined patients in alphabetical order.

Patient list:					
Lastname	Firstname	Date of birth	ID	Comment 1	
Mustermann	Annie	23.11.1986			
Demo	Normal	25.01.1966			
Demo	Normal	17.11.1976			
Demo	Patient	31.12.1960			
Demo	Patient	20.05.1926			
Demo	Patient	13.01.1950			
Demo	Patient	09.12.1964			
▶ Demo	Progression 18	05.09.1939			

Fig. 7-1: Patient list

1. Select the patient from the patient list.  
or  
Click the [Search] button and enter the last name, first name, ID number, or date of birth of the patient you are looking for.
2. Click the patient's name to transfer the data to the "Patient" window.  
The window on the right will display a list of all existing examinations for this patient.

## 7.4 Extended Patient Search: [Extended] Checkbox

➔ Click on the [Extended] checkbox.

The screen displays additional search parameters which reference previous examinations. Proceed as for the input of a patient name.




Fig. 7-1: Advanced search

## 7.5 Rename Patient Data

After creating of the patient data, you can edit it.

1. Press the [Change] button.  
The input boxes for patient data are now enabled, and the cursor jumps to the "Last name" field.
2. Change the entries in the individual boxes.
3. Press the [Save] button.

## 7.6 Exporting Patient Data

Patient and examination data can be exported, for example, for forwarding to another clinic.

1. Select the patient and also one of the examinations in the respective list as required.
2. Click [Export] button below the patient list. The following dialog appears:




Fig. 7-1: "Export patient data" dialog

No.	Description
1	Saving destination selection
2	[...] button
3	[Cancel] and [Export] buttons
4	Options for data export

The default options for import and export of data are configured in the "Settings" field, see also the [User Guide](#).

Depending on the settings, you may not have to perform all of the following steps (e.g. selection of the directory).

3. Select the "Target" (1) where you would like to export the data.

We recommend exporting the patient data using the "Single file (U12)" option.

4. Press the [...] button. (2).
5. In the dialog that appears, select the folder or the file to which the patient data should be exported.
6. Specify the name and destination of the file you are saving.
7. Make sure you have selected [Including Bitmaps].
8. Click [Export].  
The patient and examination data have now been saved at the destination specified.

You can send data stored on the hard drive as an e-mail attachment.

Requirements for transferring data to another PC:

- The Pentacam® program must be installed on the other PC. If the program is updated on the Pentacam® PC (sender), the program on the other PC (recipient) must also be updated.
- Make sure, that the PC is connected to a local network controlled by the Floating License Key or that a single license key is connected to the PC in order to interactively evaluate the exams.

## 7.7 Importing Patient Data

If you receive patient data, for example, on a USB stick, you can import these data.



Risk of loss of data due to computer viruses

Computer viruses can cause loss of data.

➔ Run a virus check before importing data from the USB flash drive.

1. Press the [Import] button. The following dialog appears:




Fig. 7-1: "Import" dialog

No.	Beschreibung
1	Select data source
2	[...] button
3	Previous examinations
4	[Import] button
5	Patient list

The default options for import and export of data are configured in the "Settings" field, see also the [User Guide](#).

➔ Depending on the settings, you may not have to perform all of the following steps (e.g. selection of the directory).

2. Select the option (1) where the source data are contained ("Folder" or "U12" (single file)).

Recommendation: Export the patient data using the "Single file (U12)" option.

3. Press the [...] button. (2).

4. In the dialog box, select the directory or the file where the patient data are located.
5. Confirm your selection with [OK] or [Open].
6. The patients and the associated examinations that are found are displayed in the lower part of the dialog.
7. To import the data, press the [Import] button (4).  
The data will then be available in the Patient Data Management system.

## 7.8 Data Backup

You should make a backup copy of patient and examination data at regular intervals. In case of loss of data, you can reconstruct the data from a previously created backup with the help of this function. Since data backup takes several minutes, depending on the scope of the database and the data to be backed up, a backup should be carried out when the PC and the device will not be needed.



### Note

Risk of loss of data due to computer viruses

Computer viruses can cause loss of data.

→ Run a virus check before making a backup to an external a USB flash drive.

The general rules for security backups apply to the creation of backup copies created with the help of the Patient Data Management system. Storage of backup files should always be done on a separate system (for example an external hard drive or a USB flash drive with adequate capacity).

### 7.8.1 Backup Data

1. Press the [Backup] button at the top right of the Patient Data Management system. The following dialog is displayed:




Fig. 7-1: "Backup" dialog

No.	Description
1	Backup data selection
2	[Restore] button
3	Display free storage space
4	Backup directory and button [...]
5	[Save] button

2. Select whether all of the data or only changed data should be backed up.

---

The Patient Data Management system internally tags all saved data records. If you select the option "Changed and new data only", only the data records that were not saved during a previous backup will be backed up.

---

3. Press the [...] button to the right of the "Backup directory" box (4).
4. In the dialog that appears, select the folder to which the data should be backed up.
5. Confirm your selection with [OK].
6. To back up the data, press the [Save] button (5). The previously selected data will then be backed up to the corresponding folder.

### 7.8.2 Reconstructing Data

If a loss of data occurs, the data from a previous backup can be re-imported into the Patient Data Management system.

1. Press the [...] button.
2. In the dialog that appears, select the folder which contains the backup data.
3. Confirm your selection with [OK].
4. To import the data, press the [Restore] button (2). All data in the appropriate directory are copied to the Patient Data Management system.

### 7.8.3 Automatic Backup

In addition to the manually performed backup, it is also possible to automatically run a backup when exiting the Patient Data Management system. The settings required for this can be made in the "Settings" area, see [User Guide](#).

In the patient data management you record and manage patient data.

## 8 Cleaning, Disinfection and Maintenance

Cleaning and disinfection of the Pentacam® and Pentacam® HR is described in this chapter.

Sterilization is not required.

- Heed the product descriptions and Instructions for Use of products and equipment you use to care for, clean, and disinfect the unit and/or its accessories.
- Do not clean the Pentacam® and Pentacam® HR with aggressive, chlorine containing, abrasive or sharp cleaning agents.

### 8.1 Intervals for Cleaning, Disinfection and Maintenance

Interval	Activity
before each use	Place fresh paper on the chin rest or disinfect the chin rest if it is used without paper
before each use	Disinfect forehead rest
monthly	Clean device (housing, illuminated slit)
monthly	Perform test measurement
every 2 years or after 25 000 measurements	Maintenance by OCULUS Service or an authorized dealer

### 8.2 Consumables

Chin rest paper	400 pcs., Order No. 65313
LidStick®	2 rolls with 100 pieces each, Order No. 77502
Desinfection wipes	mikrozyd® sensitive wipes premium Comp. Schülke & Mayr GmbH various pack sizes, e.g.: 2x 50 pcs. Softpack, Order No. 59882

### 8.3 Cleaning



#### Caution

Risk of electric shock if the Pentacam® and Pentacam® HR is not completely disconnected from the mains for these jobs.

- Turn the Pentacam® and Pentacam® HR off, → [chap. 4.6 Switching Off, page 24](#).
- Pull the power plug before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.

Required materials:

- Antistatic cleaning agent

- Cleaner for painted surfaces: Mixture of equal parts of alcohol and distilled water, possibly with a few drops of household detergent
- Soft, lint-free cloth
- Purified compressed air

### 8.3.1 Cleaning the Housing

- ➔ It is best to clean the housing surfaces with a soft cloth and an anti-static cleaning agent.
- ➔ Wipe off any residue from painted surfaces with the mixture for painted surfaces.

### 8.3.2 Clean chin rest and forehead rest

- ➔ Make sure that no liquid gets into any of the openings of the Pentacam® und Pentacm® HR.  
Do not use liquid from aerosol cans.
- ➔ Clean the chin-forehead rest with a soap solution (or with alcohol, if very dirty).
- ➔ Use a lint-free, damp cloth.

### 8.3.3 Cleaning the illuminated slit

The optical components illuminating the slit, and the lens in front of the camera are precision parts and sensitive to pressure. The surfaces of these components are susceptible to scratching.



Damage to the optics due to improper cleaning

- ➔ When cleaning the illuminated slit, do not use a cloth or cleaning agents.

- ➔ Carefully clean the centre of the illuminated slit with purified compressed air.
- ➔ Clean the lens in front of the camera using a dry, lint-free cloth.

## 8.4 Disinfection



Equipment damage caused by disinfectant solution

The disinfectant solution may damage the surface of the device if it is sprayed directly on it.

- ➔ Spray the disinfectant solution onto a cleaning cloth, do not spray it directly on the device.

- ➔ Disinfect the forehead rest after each examination.
- ➔ If you do not use paper for the chinrest, disinfect the chinrest after each examination.

## 8.5 Test measurements

The device is tested and calibrated at OCULUS prior to shipment.

In addition, OCULUS recommends performing regular test measurements.

- Begin the test by taking a measurement on a human eye.  
Perform at least 5 consecutive measurements per eye.  
Calculate the arithmetic mean and record the values.
- These measurements (as described above) should be performed once a month on the same eye.
- Compare the arithmetic mean of the initial measurement with the current measurement.

The following table describes the tolerance range between the result of the initial measurement and the result of the current measurement

		Tolerance Range
Tomography	Basic	+/- 0.2 dpt
	HR	+/- 0.1 dpt
Pachymetry		+/- 10 µm



If the difference between the initial and the current measurement is outside the tolerance range, please inform us or contact your authorized Oculus distributor to get further support.

## 8.6 Maintenance

In order to retain the high measurement accuracy of the Pentacam<sup>®</sup> und Pentacam<sup>®</sup> HR OCULUS Optikgeräte GmbH recommends to perform a maintenance service every two years or after 25000 scans. A corresponding message appears.

- By daily pop up window
- In the settings, [see User Guide](#)  
Date of next service and number of performed examinations
- In the scan menu ( → [chap. , page 29](#) ) as preliminary information (3 month before)

bzw.  
Information when service is due.



- During an examination (which is saved), an indicator appears next to the QA field.



- ➔ Let the Pentacam® und Pentacam® HR checked by our service department or an authorized dealer.

## 8.7 Attaching Paper to the Chin Rest

If you want to attach new chin rest paper, follow these instructions:




Fig. 8-1: Fasten chin rest paper

Nr.	Beschreibung
1	Pins
2	Paper for chin support
3	Chin rest

1. Pull the two pins out of the chin rest.
2. Place the chin rest paper in such a way that the holes of the paper and those in the chin rest are aligned.
3. Insert the two pins in the chin rest.

## 9 Troubleshooting



### Caution

Incorrect troubleshooting may result in personal injury or damage to the device.

- If a fault cannot be rectified, label the device as “out of order” and contact the service department or an authorized dealer.

Contact details:

- Telephone (urgent cases): +49 641 2005-800  
Have TeamViewer ready and provide the following information:
  - Customer number
  - Serial number
  - Software version
  - Description of the fault
  - Steps already taken
- E-Mail: [service@oculus.de](mailto:service@oculus.de)  
Provide the same information.
  - If necessary, include: U12 files, images
  - Large files can be sent via WeTransfer

### 9.1 Basic Troubleshooting

Action	Description
Restart	<ul style="list-style-type: none"> <li>→ Switch off the device using the on/off switch.</li> <li>→ Wait 15 seconds.</li> <li>→ Switch on the device using the on/off switch.</li> </ul>
Check connections	<ul style="list-style-type: none"> <li>→ Check that all cables are properly connected.</li> <li>→ Check the connection between the Y-cable and the power supply.</li> <li>→ Check the connection between the Y-cable and the PC.</li> <li>→ Check whether a USB extension cable is being used. This must be a repeater cable.</li> <li>→ Check whether an active USB hub (with its own power supply) is being used.</li> </ul>
Check energy options	<ul style="list-style-type: none"> <li>→ Go to Power System Settings.</li> <li>→ Deselect the [Turn on fast startup (recommended)] option.</li> <li>→ Go to Device Manager.</li> <li>→ Click the [Power Management] tab.</li> <li>→ Deselect the [Allow the computer to turn off this device to save power] option.</li> </ul>
Check USB ports	<ul style="list-style-type: none"> <li>→ Replace the USB ports on the PC.</li> </ul>
Check XY base	<ul style="list-style-type: none"> <li>→ Check whether the device is positioned parallel to the XY axis.</li> </ul>

## 9.2 Device-Specific Troubleshooting

Erros	Possible cause	Solution
After you have started the Pentacam® program, the following dialog box appears: "No communication with Pentacam!".	Power supply without voltage	<ul style="list-style-type: none"> <li>■ Check whether the indicator light on the power supply is on. If not, connect the power supply to the mains.</li> </ul>
	Power cord of the Pentacam® und Pentacam® HR not plugged in properly	Check whether <ul style="list-style-type: none"> <li>■ the power supply cable is correctly attached to the Pentacam® und Pentacam® HR.</li> <li>■ the blue slit light is visible in the Scan menu ( → chap. , page 29).</li> <li>■ the USB connector is properly inserted.</li> </ul>
	Software/hardware problems	<ul style="list-style-type: none"> <li>■ Switch the Pentacam® off and restart the PC. Switch the Pentacam® on as soon as Patient Data Management becomes active. When you start the Pentacam® program, the message, "Load Bootloader" must appear.</li> <li>■ Contact the service department or your authorised dealer</li> </ul>

## 10 Transport and Storage

The Pentacam® und Pentacam® HR, must be properly dismantled and packed before being transported or stored.



### Caution

Risk of equipment damage due to incorrect shipment or from improper storage

- Avoid shocks, vibrations, and contamination.
- Avoid high temperatures and humidity.

- Transport the Pentacam® und Pentacam® HR carefully.
- Do not hold the device by the joystick to carry it.
- Store the Pentacam® und Pentacam® HR in compliance with the storage conditions.
- Avoid placing near heaters and moisture.

### 10.1 Storage conditions

Ambient temperature range	-10°C to +55°C
Relative humidity, including condensation	10% to 95%
Air pressure range	700 hPa to 1060 hPa

### 10.2 Transport conditions

Ambient temperature range	-40°C to +70°C
Relative humidity, including condensation	10% to 95%
Air pressure range	500 hPa to 1060 hPa

### 10.3 Disassembly

1. End the current session.
2. Switch off the device.

3. Disconnect the cable from the computer/laptop and the power adapter.
4. Loosen the screw connection of the Y cable and pull it out.
  - ⚠ **When disconnecting electrical connections, pull on the respective plug and not on the cable itself.**



Fig. 10-1: Disassembly

## 10.4 Disposal



In accordance with Directive 2012/19/EC of the European Parliament and the Council and, in accordance with German law governing the marketing, return and environmentally compatible disposal of used electrical and electronic devices, such appliances must be recycled and may not be discarded as household waste.

- ➔ Dispose the Pentacam<sup>®</sup> und Pentacm<sup>®</sup> HR in a compliant manner. Dispose of the device according to legal regulations.

## 11 Terms of Warranty and Servicing

Please note our General Terms and Conditions (GTC) on our website [www.oculus.com](http://www.oculus.com).

## 12 Technical Data

### Measuring equipment

	Pentacam®	Pentacam® HR
Camera	digital CMOS camera	digital CMOS camera
Light source	Blue LED (475 nm, UV-free)	Blue LED (475 nm, UV-free)
Processor	DSP with 400 Mio. operationen/s	DSP with 400 Mio. operationen/s
Speed	50 images in 2 seconds <sup>a</sup>	100 images in 2 seconds <sup>b</sup>
Number of evaluated measuring points	max. 25.000	max. 138.000
Dimensions (W x D x H)	275 x 320 – 400 x 500 – 530 mm (10.8 x 12.6 – 15.7 x 19.7 – 20.9 in)	275 x 320 – 400 x 500 – 530 mm (10.8 x 12.6 – 15.7 x 19.7 – 20.9 in)
Weight	7,2kg <sup>c</sup> (15.9 lbs)	7,8 kg <sup>c</sup> (17.2 lbs)

a. Scheimpflug image of the entire anterior segment

b. Cornea fine scan

c. Weight without base

### Measuring range

	Pentacam®	Pentacam® HR
Curvature	3 – 38 mm 9 – 99 D	3 – 38 mm 9 – 99 D
Accuracy	± 0.2 D	± 0.1 D
Reproducibility	± 0.2 D	± 0.1 D
Working distance	80 mm (3.1 in)	80 mm (3.1 in)

### Power adapter

Power adapter	HEMG49-S240210-7 (05150150)
Mains connection	100 – 240 V AC
Frequency	50/60 Hz
Power input, max.	120-150 VA
Output voltage	24 V DC
Fuses	Integrated overcurrent shut-off

### Power supply

Output voltage		24 V DC
Max. electric consumption	Pentacam®	35 W
	Pentacam® HR	35 W

### Life expectancy

Lifecycle expectancy	up to 10 years
----------------------	----------------

### Classification

according to IEC 60601 - 1	
Type of protection against electrical shock	Protection class 2
Insulation of applied parts	Type B

**IT requirements** The IT equipment (computer, monitor etc.) must comply with the requirements of IEC 62368-1.

Recommended PC specifications	Intel® Core™ i5, 500 GB SSD, 8 GB RAM, Windows® 11, Intel® HD Graphics
-------------------------------	--

**Softwareversion**

Pentacam® Software	ab Version 1.34
--------------------	-----------------

**CE-Marking**



The unit is a class IIa product.

Conformity assessment procedure in accordance with Regulation (EU) 2017/745 on Medical Devices (MDR), Annex IX, chapter I and III.

## Annex

### Annex A Electromagnetic Compatibility

Medical electrical equipment is subject to special precautionary requirements with respect to EMC, and must be installed and operated according to the EMC-Instructions contained in the accompanying paperwork.

OCULUS devices and systems are suitable for use in professional healthcare establishments, e.g. doctor's practices or clinics; however, they may not be used in proximity to HF surgical units or inside of the HF shielded room of an ME system for magnetic resonance imaging. Portable and mobile HF communications appliances can affect medical, electric equipment.

No special measures need be observed in respect of OCULUS devices and systems.



#### Attention

**Portable and mobile RF-communications devices can interfere with electrically operated medical devices and affect the performance of the device.**

The device is intended for use in an electromagnetic environment in which the radiated RF interference is uncontrolled. The user can help prevent electromagnetic interference by maintaining the following minimum distance, based on the maximum output of the communication equipment, between portable and mobile RF communication devices (transmitters) and the device:

- Portable RF communication devices (including peripheral devices such as antenna cables and external antennas) should not be any closer than 30 cm (12 inches) to any part of the device.

#### Minimal performance quality and essential performance criteria

- A slightly disturbance of the analog camera of the Pentacam<sup>®</sup> und Pentacm<sup>®</sup> HR (slightly image noise on screen) during the examination is permissible because it will not affect the diagnosis, treatment and observation.
- A short flicker of the illumination of the Pentacam<sup>®</sup> und Pentacm<sup>®</sup> HR during the examination is permissible because it will not affect the diagnosis, treatment and observation.
- A short interruption of the USB connection during the examination is permissible because it will not affect the diagnosis, treatment and observation.



#### Caution

The use of accessories, transducers and cables specified by OCULUS with devices other then the Pentacam<sup>®</sup> und Pentacm<sup>®</sup> HR may result in increased emissions or decreased immunity of the other device.

- Do not use the accessories, transducers and cables specified by OCULUS with devices other then the Pentacam<sup>®</sup> und Pentacm<sup>®</sup> HR.
- Use only the original accessories, transducers and cables specified by OCULUS.

To be in compliance with the requirements of the IEC 60601-1-2 the following types of equipment, accessories, power adapters and cables must be used.

Model number	Article number	Description	
70700	10037829	Pentacam®	
70900	10010403	Pentacam® HR	
	05200320	Cable with connector plug, EU standard	2.5 m (98.4 in)
	05200210 (110 Volt)	Cable with connector plug, US standard	2.5 m (98.4 in)
	05150150	Power adapter HEMG49-S240210-7	24 V, 2,1A
70002	10040099	Med. secure isolator	2 m

## Annex B Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity


### Guidance and manufacturer's declaration electromagnetic emissions IEC 60601-1-2

The OCULUS Pentacam® und Pentacm® HR is intended for operation in the electromagnetic environment specified below. The user of the Pentacam® und Pentacm® HR should ensure that it is being used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Pentacam® und Pentacm® HR uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
HF-emissions CISPR 11	Class B	
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	complies	

Electromagnetic immunity, IEC 60601-1-2			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical Fast transient/bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines	± 2 kV ----- ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% $U_T$ ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree  0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single-phase: at 0 degree  0% $U_T$ ; 250/300 periods	0% $U_T$ ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree  0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single-phase: at 0 degree  0% $U_T$ ; 250/300 periods	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Pentacam® und Pentacam® HR requires continued operation during power mains interruptions, it is recommended that the Pentacam® und Pentacam® HR be powered from an uninterruptible power supply or battery.
Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			

**Electromagnetic immunity, IEC 60601-1-2**

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V<sub>eff</sub> 150 KHz to 80 Mhz 6 V in ISM- and amateur radio frequency bands between 150 kHz and 80 MHz 80% AM to 1 kHz</p> <p>3 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz</p>	<p>V<sub>eff</sub> = 3 V</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of Pentacam® and Pentacam® HR, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3,5}{(V_1)} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{(E_1)} \right] \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{(E_1)} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interferences may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 Hz and 800 MHz, the higher frequency range applies.  
 Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Pentacam® and Pentacam® HR is used exceeds the applicable RF compliance level above, the Pentacam® and Pentacam® HR should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Pentacam® and Pentacam® HR.  
 b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the Pentacam® und Pentacam® HR, IEC 60601-1-2**

The Pentacam® und Pentacam® HR in intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Pentacam® und Pentacam® HR can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pentacam® und Pentacam® HR as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 Mhz $d= 1.2 \sqrt{P}$	80 MHz to 800 MHz $d= 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d= 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.80	3.80	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Annex C Description of the Connection



## Annex D Data sheet HEMG 49-S240210-7 (05150150)

# HiTRON

## UNIVERSAL INPUT AC-DC MEDICAL & ITE APPLICATION EXTERNAL DESKTOP SWITCHING ADAPTER 48-50 WATTS GREEN POWER SINGLE OUTPUT HEMG49 SERIES



### FEATURES:

- ACCOMMODATE UNIVERSAL AC INPUT
- MEET MEDICAL STANDARDS IEC60601-1 & ITE STANDARDS IEC60950-1
- EMI MEET EN 55011 & EN55022 / FCC CLASS B
- MEET ENERGY STAR LEVEL V & CEC LEVEL IV
- CE MARKING COMPLIANCE

### SPECIFICATION

#### INPUT SPECIFICATION

**Input Voltage:** Typical 90-264Vac.  
**Input Connector:** 3 pole AC inlet IEC320-C14(DT7) / 2 pole AC inlet IEC320-C8(DT8).  
**Input Frequency:** 47-63Hz.  
**Inrush Current:** 12Arms (52Apk) at 230Vac.  
**Input Current:** Typical 0.91A at 115Vac/ 0.57A at 230Vac.  
**Dielectric Withstand:** Meet IEC60601-1 & IEC60950-1.  
**EMI:** Meet EN55011 & EN55022 / FCC Class B.  
**Hold-up Time:** Typical 12mS at 115Vac. Typical 70mS at 230Vac.  
**Over Temp. Protection:** Optional (NTC circuit).  
**Earth Leakage Current (Class I) :** Less than 0.3 mA.  
**Touch Leakage Current (Class I & II) :** Less than 0.1mA.  
**No Load Power:** Less than 0.3W at 230Vac

#### OUTPUT SPECIFICATION

**Output Voltage:** See Ratings Chart.  
**Output Current:** See Ratings Chart.  
**Output Wattage:** Typical 48-50Watts.  
**Output Connector & Cord:** Optional.  
**Line Regulation:** Typical 0.1%.  
**Load Regulation:** Typical  $\pm 1.5-3.0\%$ .  
**Noise & Ripple:** 1.0% peak to peak.  
**OVP:** Built-in by latch circuit.  
**Adjustability:** Factory set.  
**Over Current Protection (OCP):**  
 Fully protected against output overload and short circuit.  
 The PSU will shut down after OCP is activated.  
 Consult the factory for OCP setting.

#### GENERAL SPECIFICATION

**Efficiency:** Typical 87%-88% (various with the output voltage)  
**Switching Frequency:** Typical 65KHz.  
**Circuit Topology:** Fixed Frequency Flyback circuit.  
**Transient Response:** Output voltage returns in less than 5.5mS following a 50% load change.  
**Safety Standard:** Meet Medical IEC60601-1 & ITE IEC60950-1, Class I for DT7(C14) or Class II for DT8(C8)

**Operating Temperature:** 0°C to +40°C.  
**Storage Temperature:** -20 to +85°C.  
**Cooling:** Free air convection.  
**Construction:** Impact resistant thermo-plastic enclosure case.  
**Power Density:** 3.14-3.27Watts. / Cubic inch.  
**Desktop Format.**

NOTE: (1) All measurements are at nominal input, full load, and +25°C unless otherwise specified.  
 (2) Load regulation is measured at 115Vac or 230Vac in percentage to indicate the change in output voltage as the load varied from half load to full load ( $\pm\%$ ).  
 (3) The exact obtainable load regulation depends upon the output cord selected and load current.  
 (4) Due to requests in market and advances in technology, specifications subject to change without notice.



For the details of safety approval, please consult the factory.


**OUTPUT VOLTAGE / CURRENT RATINGS CHART**

**SINGLE OUTPUT**

MODEL NO.	AC INLET	O/P VOLTAGE	O/P CURRENT
HEMG49-S120400-7	IEC320-C14(DT7)	12.0Vdc	4.0A
HEMG49-S120400-8	IEC320-C8(DT8)	12.0Vdc	4.0A
HEMG49-S150330-7	IEC320-C14(DT7)	15.0Vdc	3.3A
HEMG49-S150330-8	IEC320-C8(DT8)	15.0Vdc	3.3A
HEMG49-S240210-7	IEC320-C14(DT7)	24.0Vdc	2.1A
HEMG49-S240210-8	IEC320-C8(DT8)	24.0Vdc	2.1A

**MECHANICAL DIMENSIONS: MM [INCHES]**

**WEIGHT: 373.0g (13.2 Oz.)**



## Annex E Instructions for integration into an IT-Network

The device, in conjunction with the connected computer and the device software running on it, forms a programmable electrical medical system (PEMS) according to IEC 60601-1.

It is essential that you follow the chap. "2.3 Instructions on Cyber Security" on page 14 of "Safety Instructions for Use" (→ page 12) in the device Instructions for Use.

Please note the following information for implementing an integration of the PEMS into an IT-Network:

### **The purpose of integrating the PEMS into an IT-Network can be:**

- Licensing by local license server
- Storage and retrieval of the examination data on a local network drive
- Printing
- Data export
- DICOM workflow

### **Required characteristics of the IT-Network into which the PEMS is to be integrated:**

- Prefer wired LAN connection
- IPv4 network
- Fast-Ethernet (100 Mbit/s minimum)

### **Required configuration of the IT-Network into which the PEMS is to be integrated:**

- Licensing: required opened ports: 3968 TCP; 51371 - 51372 UDP
- Storage, printing, data export: File and Printer Sharing for Microsoft Networks (SMB 3.0 or higher - required opened port: 445]
- DICOM Storage Service Class = PACS
- DICOM Worklist Management Service Class (Modality Worklist Server)

### **Technical specifications of the network connection with the PEMS, including the specifications of data security:**

- Refer to the „Instructions on Cyber Security" section (chap. "2.3 Instructions on Cyber Security" on page 14).
- Refer to the "Floating License Key - License management for software options" Instructions for Use (if applicable)
- Refer to device specific DICOM interface description (if applicable)

### **The intended flow of information between PEMS, the IT-Network and other devices in the IT-Network and the intended routing through the IT-Network**

- License handling from local license server to PEMS and vice versa
- Storage and data export to local network storage and loading from local network storage
- Printout to local printer

### **List of hazard situations that result from the IT-Network being unable to provide the features that are required to fulfill the purpose of integrating the PEMS into to meet the IT-Network:**

- Loss of data
- unsuitable data exchange
- Data corruption
- unsuitable temporal data allocation
- unexpected data reception
- unauthorized access to data



Connection of the PEMS to an IT-Network that includes other equipment could result in previously unidentified risks to patients, operators or third parties.

The responsible organization should identify, analyze, evaluate and control these risks.

Subsequent changes to the IT-Network could introduce new risks and require additional analysis.

**Changes to the IT-Network include:**

- changes in IT-Network configuration
- connection of additional items to the IT-Network
- disconnecting items from the IT-Network
- update of equipment connected to the IT-Network



WWW.OCULUS.DE

**OCULUS Optikgeräte GmbH**  
Münchholzhäuser Straße 29 • 35582 Wetzlar • GERMANY  
Phone +49 641 2005-0 • Fax +49 641 2005-255  
E-mail: sales@oculus.de • www.oculus.de

10028075 / EN / Rev05  
Lot:

