

OCULUS Corvis® ST



INSTRUCTION MANUAL
Tonometer - Pachymeter

Preface

Thank you for your purchase and the trust you have placed in this OC-ULUS product. The Corvis® ST has been manufactured and tested according to strict quality criteria. You have selected a modern and well-engineered product.

To ensure safe operation, it is essential that you use the device correctly. For this reason you should familiarise yourself thoroughly with the contents of this instruction manual before operating the device. In particular, pay attention to the safety instructions.

- This instruction manual describes the measuring procedure of the Corvis® ST.

Due to ongoing development, the diagrams shown may depict minor changes to the actual device delivered.

If you have any questions or wish to obtain further information on your device, please contact us by phone, email or fax. Our service team is pleased to assist you.

OCVLUS Optikgeräte GmbH

Item Number: G/72100/XXXX/EN

Revision: 06

Release: 01.05.2024

Table of contents

1	Scope of Delivery	9
2	Symbols	10
2.1	On the Device / Type plate	10
2.2	On the Packaging	11
3	Structure of the Documentation	12
4	Safety Instructions	13
4.1	Graphic symbols used.....	13
4.2	Safety Instructions for Use	14
4.2.1	Instructions for Operating Personnel.....	14
4.2.2	Transport and Storage Instructions	14
4.2.3	Instructions for Setup and Connection.....	14
4.2.4	Patient environment information	15
4.2.5	Information about the operation of an ME system.....	16
4.2.6	Instructions for Operation	16
4.2.7	Instructions for Maintenance.....	16
4.2.8	Instructions for Disassembly and Disposal.....	16
4.2.9	Instructions on Electrical Safety.....	17
4.3	Cyber Security Instructions	18
4.3.1	Data responsibility:.....	18
4.3.2	Device Security.....	18
4.3.3	User Responsibility	18
4.3.4	Reporting Device Security or Privacy Breaches	18
4.3.5	Recovering from Compromised Accounts or Devices	19
4.3.6	Unavailable Service.....	19
4.3.7	Precautions.....	19
4.3.8	Precautions for access control of the computer	19
4.3.9	Precautions if the computer is connected to a LAN or internet network.....	20
5	Intended Use.....	21
5.1	Intended purpose	21
5.2	Intended medical indication.....	21
5.3	Contraindication.....	21
5.4	Possible side effects.....	21
5.5	Intended users	21
5.6	Patient group.....	21
5.7	Disclosure of residual risks.....	22
6	Transport to Installation Location.....	24

7	Device Description	25
7.1	Overview of Device Components	25
7.2	Mode of Operation of the Corvis® ST	27
7.3	Applied parts	28
8	Start-up	29
8.1	Initial Start-up	29
8.2	Set-up Jobs for Initial Start-Up.....	29
8.3	Adjustments after an in-house transport	30
8.3.1	Device set-up.....	30
8.3.2	Unlock transport safety device	30
8.4	Electrical Connection.....	31
8.5	Daily operation	32
8.5.1	Switching on the Corvis® ST.....	32
8.5.2	Setting the Safety Stop	32
8.5.3	Adjust the break	33
8.5.4	Switching off the Corvis® ST	33
9	Preparing a measurement.....	34
9.1	Using the control unit.....	34
9.2	Display with Touch Screen	36
9.3	Start Patient Data Management	37
9.4	Importing patient data	38
9.4.1	Enter a new patient	39
9.4.2	Select an existing patient.....	40
10	Perform a Measurement	41
10.1	Start the Corvis® ST Program at the Computer.....	41
10.2	Adjusting the Corvis® ST	42
10.2.1	Rough Adjustment.....	43
10.2.2	Fine Adjustment.....	45
10.3	Triggering the measurement	47
10.4	Perform a second measurement.....	47
10.5	Saving data.....	48
10.6	Complete measurement.....	48
10.7	Using Patient Data Management	49
10.7.1	Rename Patient Data	49
10.7.2	Exporting Patient Data.....	50
10.7.3	Data backup	51
10.7.4	Backup data.....	51
10.7.5	Reconstruct data	52
10.7.6	Automatic backup.....	52
11	Working with the Corvis® ST program.....	53
11.1	View patient and examination data (1)	53
11.2	Working with the Menu Bar (7).....	53

- 11.2.1 Menu item "Patient" 54
- 11.2.2 Menu item "Examination" 54
- 11.2.3 Menu item "Display" 54
- 11.2.4 Menu item "Export" 55
- 11.2.5 Menu item "Settings" 55
- 11.3 Use video function (4)..... 55
- 11.4 View specific images of the cornea (4)..... 55
- 12 Perform a measurement with the Corvis® ST unit 56**
- 12.1 Change Settings 56
- 12.2 Entering Patient Data 57
 - 12.2.1 Entering new Patient..... 58
 - 12.2.2 Select an existing patient..... 59
- 12.3 Perform a Measurement 60
 - 12.3.1 Selecting a Measurement Mode 60
- 12.4 Complete measurement..... 60
- 12.5 Re-use examination data..... 61
 - 12.5.1 Delete a Patient or an Examination..... 61
 - 12.5.2 Re-using data with the USB flash drive (optional) 61
- 13 Cleaning, Disinfection and Maintenance 63**
- 13.1 Cleaning..... 63
 - 13.1.1 Clean the front panel..... 64
 - 13.1.2 Clean head (1) and chin (4) rest..... 64
 - 13.1.3 Clean lens protection glass (3)..... 65
 - 13.1.4 Clean the air nozzle (2) 65
 - 13.1.5 Clean the housing..... 65
 - 13.1.6 Clean the display 66
- 13.2 Disinfection 66
- 13.3 Maintenance 67
- 13.4 Attach paper to the chin rest 67
- 13.5 Insert new printing paper roll 68
- 14 Dismantling, Transport and Storage..... 70**
- 14.1 Parking position 70
- 14.2 Lock the transport safety device..... 71
- 14.3 Joystick lock 72
- 14.4 Transport and Storage Information..... 73
- 14.5 Transport and Storage..... 73
- 15 Disposal of Used Devices..... 74**
- 16 Troubleshooting..... 75**
- 17 Terms of Warranty and Servicing..... 76**
- 17.1 Liability for proper function or damages..... 76

18 Technical Data	77
19 Appendix.....	80
19.1 Correction Calculation of the Tonometrically calculated IOP.....	80
19.2 Electromagnetic Compatibility	80
19.3 Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity for the Corvis® ST.....	82
19.4 Description of the Connection	86
19.5 Data sheet Power Supply Adapter (05150285)	87
19.6 Instructions for integration into an IT-Network.....	90

1 Scope of Delivery

Component	Order number
■ Corvis® ST	72100
■ Corvis® ST package with laptop en	72200
■ Dust protection cover	026010005001
■ Paper for chin rest	65313
■ Paper roll (3 rolls)	65311
■ USB cable	05200601
■ USB FS MED Isolator	015692000010
■ Power supply GSM90B15-P1M	05150285
■ Instruction Manual	G/72100/XXXX/EN
■ User Guide	BH/72100/EN...
■ Software Installation	SI/50000...en



We reserve the right to change the scope of delivery in line with ongoing technical development.

- ➔ When checking the delivery, if you discover transport damages, immediately make your claim with the transport company.
- ➔ Have the damage confirmed on the bill of lading, so that a proper claim settlement is possible.
- ➔ Keep the packing material.

2 Symbols

2.1 On the Device / Type plate




Fig. 2-1: Type plate Corvis ST (example)




Fig. 2-2: Additional type plate Corvis ST (example)

No.	Symbol	Description	No.	Symbol	Description
1		Logo of manufacturer	9		Device type
2		Name and address of manufacturer	10		Attention
3		Date of manufacture	11		Applied part Type B
4		UDI number (example), consisting of: UDI-DI (Device-Identification) UDI-PI (Product Identifier) and machine-readable matrix code	12		Disposal in household trash is prohibited
5		Conformité européenne	13		Protection class












No.	Symbol	Description	No.	Symbol	Description
6		Medical device	14	IPxx	Type of protection
7		Serial number	15		Follow instruction for use
8		Article number			



Fig. 2-3: Type plate measure head

There are no temporarily fixed markings on the device.

2.2 On the Packaging

Symbols packaging	
	Keep dry
	This way up
	Fragile
	Limit of temperature for transport
	Limit of temperature for storage
	Limit of humidity
	Limit of air pressure

3 Structure of the Documentation

A folder containing documentation is supplied with your Corvis® ST:

- **Instruction Manual:** The design of the unit is described in detail in this document. The instruction manual also gives you general information about working with the Patient Data Management system and all safety-related instructions for use of the Corvis® ST.
- **User Guide:** All features of the examination and analysis software are described in the User Guide, along with detailed information about the Patient Data Management system.
- **Software Installation:** The introduction to the Software Installation describes how to install the Corvis® ST software and the associated drivers.

4 Safety Instructions



Attention

All safety-related instructions for use of the Corvis[®] ST are given in the Instruction Manual for the unit. It is imperative that you read and understand the whole Instruction Manual before you use the Corvis[®] ST.

- ➔ Read the instruction manuals carefully.
- ➔ Carefully store the instruction manuals near the device.
- ➔ Observe the legal regulations with regard to accident prevention. If standards are named without an issue date, the current version always applies.

4.1 Graphic symbols used



Warning

Identifies a potentially dangerous situation which may cause serious bodily injuries.



Attention

Identifies a potentially dangerous situation which may cause minor injury or damage to property.



Note

Identifies Instructions for use, and useful or important informations.



Identifies important information about the product or its operation, which require special attention.

- > This symbol denotes menu paths and screen shots. Example for starting a new examination:
Corvis[®] ST > Examination > Scan
which means:
 - ➔ Select the "Examination" menu from the menu bar.
 - ➔ Select the menu item "Scan".

4.2 Safety Instructions for Use



Attention

Personal injury or property damage due to improper operation

- Observe the following safety instructions.



Attention

Personal injury or property damage due to equipment modifications that could jeopardize safety.

- No modifications may be made to this device without the permission of the manufacturer. Only the OCULUS service are allowed
 - to modify the device or the associated lifting table
 - to install software and software updates.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer (vigilance@oculus.de) and the competent authority of the region in which the user and/or patient is established.

4.2.1 Instructions for Operating Personnel

- Refer to the notes in chapter '5.5 Intended users' on page 21.

4.2.2 Transport and Storage Instructions

- Refer to the notes in chapter '14 Dismantling, Transport and Storage' on page 70.

4.2.3 Instructions for Setup and Connection

- Only OCULUS or an authorized dealer is allowed to set up and to connect the Corvis® ST.
- Do not use or store the Corvis® ST in rooms that are humid, see the transport and storage conditions in chapter '14 Dismantling, Transport and Storage' on page 70.
- Keep the Corvis® ST away from water that may drip, splash or spray on it, and make sure that no liquids can get into the Corvis® ST. Do not place any containers holding liquids in the vicinity of the Corvis® ST.
- Germany: Only operate the Corvis® ST in rooms used for medical purposes after they have been set up according to the VDE Regulation 0100-710.
- Do not operate the devices included in the delivery in areas where explosions may occur, or in proximity to flammable anesthetics or volatile substances such as alcohol, benzine or similar products.
- Set up the Corvis® ST so that the power plug is easy to access. That way, you can easily disconnect it from the power supply for any repairs or maintenance work.

- ➔ Do not force any plug connections.
If you are unable to make a plug connection, check whether the plug fits the socket.
If you detect damage to the connection, you should let authorized dealer repair the damage.
- ➔ Only use a device which is mounted at the lifting table properly.

4.2.4 Patient environment information

Patient environment is the area where patients can come into contact with any part of a medical electrical equipment (ME equipment) or with another person being in contact with the ME equipment.



Attention

In the patient environment, use devices that conform to IEC 60601-1. If a multiple power socket is to be used, or if a device is to be used that does not meet the IEC 60601-1 standard, use an isolating transformer.




Fig. 4-1: Patient environment

4.2.5 Information about the operation of an ME system

The Corvis® ST and a connected computer form a medical electrical system (ME system) according to IEC 60601-1. If you connect additional devices, such as, for example a printer, those devices become part of the ME system.

- Make sure that all devices of the ME system meet the requirements of IEC 60601-1 or IEC 62368-1.

4.2.6 Instructions for Operation

- Before first use: Let OCULUS or an authorized dealer train you in the operation of the Corvis® ST.
- Never operate a damaged Corvis® ST.
- Only operate the Corvis® ST with the original accessories supplied by us and only when the unit is in technically perfect condition. Use only the power supply unit specified in the scope of delivery.
- Do not cover the ventilation openings.
- Do not touch the patient and the device at the same time.
- Make sure that the device cannot tip over by leaning against it or sitting on it.
- Do not put the Corvis® ST, including rechargeable battery or cable, down onto devices that produce heat, heaters (for example radiators), microwaves or similar.
- Only operate the device if you have understood the operating instructions.

4.2.7 Instructions for Maintenance

- Germany: The operating company must ensure that the device undergoes technical measurement testing every 2 years according to MPBtreibV, Appendix 2 Tonometer.
- Only authorized persons are allowed to perform the software updates.

To ensure that it functions correctly and safely we recommend the following:

- Have the Corvis® ST checked every two years by our service department or an authorized dealer.

If an error occurs which you cannot correct, label the Corvis® ST as being "out-of-order" and contact our service department.

4.2.8 Instructions for Disassembly and Disposal

- When disconnecting electrical connections, pull on the respective plug and not on the cable itself.
- Dispose of the device according to legal regulations.

4.2.9 Instructions on Electrical Safety


Attention

Risk of personal injury or damage to property due to an incorrect level of safety

Connecting the Corvis® ST with its non-medical electrical equipment (e.g. data processing equipment) to a medical electrical system must not result in a patient safety level below that prescribed by IEC 60601-1. If making this connection leads to the leakage current threshold being exceeded, protective measures including a circuit breaker must be in place.

- Ensure that connections with non-medical devices are made correctly.
- Only use the power adapter listed in the packing list.
- Use only a computer that meets the specifications given in this instruction manual, *sec. 18, page 77*.


Attention
Use of a multiple socket extension cord

Risk of personal injury or material damage caused by unsafe multiple socket extension cord

If you use a multiple socket extension cord to connect the Corvis® ST to the power supply, you must heed the following information:

- Use an extension cord that complies with the requirements of IEC 60601-1, section 16.
- Do not place the multiple socket extension cord on the floor.
- Do not use more than one multiple socket extension cord.
- Plug only the Corvis® ST and the computer that is being used with the unit (if applicable) into the multiple socket extension cord.

If you are using a multiple socket extension chord it has to be supplied with a isolation transformer.

If you are using a new computer for the Corvis® ST, you must have the electrical safety checked. Call OCULUS Service for this purpose.


Attention
Electromagnetic Compatibility (EMC) / Cables

Risk of personal injury or damage to property due to electromagnetic interference

Portable and mobile RF communications equipment can affect medical electrical equipment *sec. 19.2, page 80*.

- Make sure that portable and mobile RF communications equipment do not cause interference.
- Recommendation: Maintain a minimum distance of 4 m. If the distance is shorter, you must ensure that the Corvis® ST functions correctly.

4.3 Cyber Security Instructions



The device itself was not designed to connect via a coupled computer to the internet or other network, or to portable media, as the device does not require a network or internet connection to function.

User who elect to connect the device coupled computer to the internet or other network for other purposes are responsible for ensuring that it is done so in a controlled manner.

4.3.1 Data responsibility:

The device itself is not designed to connect with the internet, but only with a computer. It does not require Internet access for regular use.

Do not connect with the internet while using the device. It is considered misuse.

If you elect to connect the computer to the internet for other purposes you are responsible for ensuring data security.

4.3.2 Device Security

It is the responsibility of the authorized user to ensure that the Corvis® ST device is not left unlocked, or otherwise unsecured when not in use, to ensure that non-authorized medical, professional, or otherwise unapproved personnel are not exposed to, or gain access to, ePHI.

4.3.3 User Responsibility

User names or passwords must not be shared with colleagues or others, even if they are permitted by law and provider policy to view the same type of information (e.g. two operators reviewing the same patient samples).

Operators have access to patients' ePHI and must not take snap-shots, screenshots or pictures (e.g. using another device) of any information viewed through the device.

Operators should not enter identifying data into the device. All the data on the device should be de-identified and relate to the sample id and not to the patient.

4.3.4 Reporting Device Security or Privacy Breaches

Operators must contact their local IT department and disclose any suspected or confirmed compromised user accounts, and any other privacy or security breaches.

4.3.5 Recovering from Compromised Accounts or Devices

When accounts are considered compromised, devices are lost, or unauthorized access is discovered or suspected, the Healthcare Organization's IT network administrators suspend and modify the user login criteria and issue new login credentials for the user to access their account securely.

4.3.6 Unavailable Service

Users should report unavailable service or prohibited access to information to their local Healthcare Organization's IT department.

4.3.7 Precautions

- To ensure cyber security in order to the usage of the device, the following security measures should be considered. Contact your computer administrator:

4.3.8 Precautions for access control of the computer

- Secure the computer with a password (for example at Windows start up).
- Choose a complex password: A good password should be at least eight characters long and are not in the dictionary. In addition to letters, it should also include numbers and special characters.
- Do not choose a name or device name for a password (for example "Corvis").
- Change the password regularly.
- Do not note the password in an accessible location.
- Use different passwords for different users.
- Enable the screen saver and use the option for the necessity of reentering the password when exit the screen saver.
- Choose an adequate time setting for starting the screen saver if software session is inactive (e.g. 10 minutes). Adequate time setting should consider duration of examination, number of patients, time between examinations, use of other devices in the examination room, several user, etc.
- Lock the computer if you are leaving the workstation (shortcut: 'windows logo key' + 'L').

4.3.9 Precautions if the computer is connected to a LAN or internet network

- If you elect to connect the computer to the LAN or internet, you are responsible for ensuring data security.
- Prefer wired connections of the computer to the network.
- If you are using Wi-Fi connections nevertheless, please ensure the usage of adequate security methods (for example WPA2/AES – Wi-Fi Protected Access / Advanced Encryption Standard – with a strong network key).
- The usage of a firewall (software or hardware) is recommended.
- Heed the instructions for integration into an IT-Network [sec. 19.6, page 90](#)



Do not use the Corvis® ST with wireless technology, for example with wireless USB (connection between device and computer)

5 Intended Use

5.1 Intended purpose

The product is a non-contact tonometer and pachymeter that provides additional information about the biomechanical response of the cornea. The device is a non-contact tonometer which measures intraocular pressure without contact with the eye by applying an air puff to the eye. During the air puff the eye gets illuminated by a 9 mm slit through the apex.

The Corvis® ST may only be used for the purpose described in this instruction manual.

➔ Heed the safety instructions listed above.

5.2 Intended medical indication

The biomechanical data provided is intended as an additional factor that can be used by ophthalmologists to assess eye diseases such as keratoconus and glaucoma.

5.3 Contraindication

- Ruptured globe
- Unhealed corneal abrasions
- Children under 5 years

5.4 Possible side effects

none known

5.5 Intended users

The Corvis® ST is exclusively intended for use in:

- Ophthalmological practices
- Clinics
- with ophthalmologists, opticians or optometrists

The Corvis® ST is intended for use by trained personnel:

- who can guarantee proper handling due to their knowledge, training and practical experience.
- who have been instructed by OCVLUS staff or an authorized dealer before the initial operation.

5.6 Patient group

Children from 5 years up to geriatric patients. No restrictions on weight, health and condition. Patient is awake and able to understand and to look into a fixation target.

5.7 Disclosure of residual risks

The Corvis® ST is designed to safely fulfill its intended use. Possible risks were identified as far as possible and the probability of occurrence is reduced as far as possible. At least up to the state of the art.

For example Risk control measures:

- Safety by design
- Additional risk minimisation
- Safety instructions and warnings
- Notes and descriptions in the instructions for use.

Nevertheless, these risks, which are inevitable due to the device technology and the use on patients, can only be reduced in their probability of occurrence, not completely ruled out.

The following residual risks are possible when using the device (e.g. if the information in the instructions for use is not followed or the device is not used in accordance with its intended purpose):

- Electric shock
- Fire
- Light hazard
- No, incomplete or invalid examination results
- Pinching of fingers
- Contusion of body parts (hand, arm, leg, foot, etc.)
- Infections
- Skin injuries, irritations, reactions
- Corneal damage for patient
- Discomfort for patients
- Teary eye (short-term from air puff)
- Irritated eye (short-term from air puff)
- Data loss
- Lack in data security (privacy of data)
- Device damage

These residual risks are also addressed in several places in these instructions for use with warnings, notes and descriptions. If the information in the instructions for use is observed and the device is used in accordance with its intended purpose, the clinical benefit clearly outweighs these residual risks.

Clinical Benefit and claims

Apart from the product statements in the intended use and the indications, the manufacturer provided a brochure with the following clinically relevant claims:

- precise measurement of IOP
- precise measurement of corneal thickness
- detailed assessment of corneal biomechanical properties
- calculate a biomechanically corrected IOP (bIOP)
- allows ectatic diseases such as keratoconus to be detected at a very early stage
- biomechanical properties also play an important role in the development and progression of glaucoma

These are defined through the following outcome parameters which can be measured by Corvis ST:

- Assessment of the biomechanically corrected IOP (bIOP)
- Assessment of the Corvis Biomechanical Index (CBI), such as for keratoconus detection
- Assessment of the tomographic/topographic Biomechanical Index (TBI)
- combination of Corvis ST with the Pentacam for enhanced ectasia detection
- CBI-LVC measures biomechanical stability after laser vision correction.
- Assessment of the Biomechanical Glaucoma Factor (BGF)
- Assessment of ectatic diseases at a very early stage
- Assessment of the corneal thickness

From these, a general benefit is defined as:

In accordance with the state of the art the Corvis ST supports eye doctors in evaluating conditions of the eye, such as glaucoma and keratoconus.

6 Transport to Installation Location

The transport and storage conditions see [sec. 14, page 70](#).

- Wait approx. 3-4 hours after transport before operating the Corvis® ST. Extreme temperature changes from cold areas to warm rooms can cause condensation on the optical components.



Note

Equipment damage due to incorrect transport and improper storage

- Avoid shocks and vibration.
- Avoid contamination, high temperatures and humidity.

-
- Transport the Corvis® ST professionally.
 - Store the Corvis® ST according to the storage conditions.
 - Avoid placing near radiators and moisture.



Note

- Keep the packing material. You can then ship or transport the unit in the proper manner for any servicing or repairs that may arise. You can thus avoid unnecessary damage and costs.
-

7 Device Description

7.1 Overview of Device Components

Side - view




Fig. 7-1: Corvis® ST: Side view

1 Measure head

2 Chin rest

3 XY-base

4 Control knob

5 Joystick

6 Function keys

7 Display

8 Printout slot

Front - view



Fig. 7-2: Corvis® ST: Front view and connections

- | | | | |
|---|--------------------------------------|----|------------------|
| 1 | Head rest | 6 | On/Off switch |
| 2 | Air nozzle / window of the slit lamp | 7 | Control LED |
| 3 | Lens protection glass | 8 | Mains connection |
| 4 | LED to illuminate the eye | 9 | USB port |
| 5 | Marking for the eye height | 10 | Safety stop |

7.2 Mode of Operation of the Corvis® ST

The Corvis® ST is a non-contact tonometer equipped with an optical Pachymetry function.

The Corvis® ST measures intraocular pressure without contact with the eye by applying an air puff to the eye. During the air puff the eye gets illuminated by a 9 mm slit through the apex and a built-in high-speed camera records the movement of the eye with more than 4000 images per second.

The high-speed camera uses a sequence of 140 Scheimpflug images of the cornea which are analyzed by a built-in computer.

Intra ocular pressure is determined by the detection of the applanation moments of the cornea.

Based on the Imbert-Fick principle the intraocular pressure is calculated by dividing the amount of air pressure into the area of applanated surface.

The device increases the air pressure puffed onto the cornea in proportion to time. The shape of the cornea changes from the normal convex surface to a concave surface.

This change is optically detected within 140 Scheimpflug images. The device calculates the time required to applanate (plane shape) the cornea with the air puff.

The slit light illuminates a sectional plane from the front surface of the cornea to the back surface during the air puff. The transparent cells of the cornea scatter the slit light such that the sectional plane appears as if it were self-luminous.

This is captured at an angle of 45° through the pupil by a camera, whereby the image plane of the camera is also tilted 45° to the optical axis of the camera lens, in order to sharply focus the light-scattering cornea plane onto the image plane of the camera (Scheimpflug image).

Due to this arrangement, sharp sectional images of the cornea can be obtained.

The pachymetry principle uses also the sectional images of the cornea. Corneal thickness and shape are obtained from images when the Cornea is not influenced by the air puff.

Tonometry and Pachymetry are measured at the same measurement process.

7.3 Applied parts




Fig. 7-3: Applied parts

- 1 Head rest
- 2 Chin rest

8 Start-up

Before you can operate the Corvis® ST for the first time, you must follow the instructions in chapter '8.1 Initial Start-up' on page 29.

If you want to put the Corvis® ST into operation after an in-house transport, follow the instructions in chapter '8.3 Adjustments after an in-house transport' on page 30.

8.1 Initial Start-up

Before you can operate the Corvis® ST for the first time, you must

- set it up and adjust it
- get trained



Attention

Incorrect measurements / equipment damage due to a lack training

- Before first use: Let OCULUS or an authorized dealer train you in the operation of the Corvis® ST.

Incorrect measurements / equipment damage due to incorrect set-up

- Before the first use, make sure the installation and connection of the "Corvis® ST" examination area is completed by our service or by a professional authorized by OCULUS.



Note

- Do not expose the Corvis® ST to any vibrations, shocks, contaminants, moisture, or high temperatures.
- Handle the optical device with care.

8.2 Set-up Jobs for Initial Start-Up

- Wait approx. 3-4 hours after transport before operating the Corvis® ST. If the Corvis® ST was stored in a cold room or vehicle during the cold time of the year, a significant change in temperature may cause condensation to appear on optical parts of the Corvis® ST.
- Check if the transportation safety device is unlocked, [sec. 8.3.2, page 30](#).

8.3 Adjustments after an in-house transport



Note

Equipment damage due to incorrect lifting

If the Corvis® ST is lifted only by the measure head, it can break off.

→ Grab the Corvis® ST from below and the head rest to lift it.

8.3.1 Device set-up

- Place the Corvis® ST on a level surface.
- Place the Corvis® ST so that no direct light can effect the measurement.
- Avoid shocks and vibration.
- Avoid contamination, high temperatures and humidity.

8.3.2 Unlock transport safety device

During transport the Corvis® ST is secured with a transport safety device. This must be unlocked before use.

→ Open the cover with the display.




Fig. 8-1: Open the cover with the display


**Note**

Risk of equipment damage due to incorrect connection

If you do not connect the Corvis® ST properly, and the connection is live, the unit can be damaged within a short period of time.

- Do not use excessive force when connecting the electrical plug.
- Pay attention to the specifications on the nameplate.

If the electrical plug is damaged, contact our service department or an authorized dealer to repair the damage.



- Connect the unit to the power supply with the power cable provided (2).
- If desired, connect the device to your computer/laptop using the USB port (1), with a USB cable with an USB FS MED Isolator.


8.5 Daily operation

8.5.1 Switching on the Corvis® ST



- Make sure that the mains voltage is the same as the voltage specified on the rating plate.
- Switch on the Corvis® ST with the On/Off Switch (1).

8.5.2 Setting the Safety Stop



The safety stop (2) is a lock that prevents the air nozzle from touching the patient's eye.

You can independently determine the position of the stop.

- Press down on the safety stop (2) and hold it.
- Move the Corvis® ST into the desired position.
- Release the safety stop (2).

The safety stop is set. You can only move the Corvis® ST as far as this position. You can move the Corvis® ST towards the rear at any time.


Attention

Risk of contact of the patient's eye with the air nozzle

- ➔ Before starting a measurement, make sure the safety stop is set correctly. This prevents the air nozzle from touching the patient's eye.

8.5.3 Adjust the brake



The brake (1) prevents the Corvis® ST from moving fast and jerkily on the XY-base. This enables you to better control the position of the device.

- ➔ Turn the brake to the desired position.
To the right: Corvis® ST is hard to move
To the left: Corvis® ST is easier to move

8.5.4 Switching off the Corvis® ST

- ➔ End the current session.
- ➔ Switch off the Corvis® ST with the On/Off Switch (1).


Attention

Risk of electric shock if the Corvis® ST is not completely disconnected from the mains for transport, cleaning, maintenance, disinfection and repair.

- ➔ Turn the Corvis® ST off, [sec. 8.5, page 32](#).
- ➔ Pull the power plug before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.

9 Preparing a measurement



In order to transfer data from the device to your computer/laptop you have to install the following software on your computer:

- the Corvis® ST program
- the USB connection must be enabled in the „Settings 2: USB Transfer“ of the Corvis® ST device.
- the Patient Data Management software

For more information see the [User Guide](#).

9.1 Using the control unit


You enter and manage patient data with the control panel. In addition, you start the measurements and can view the results with it.



Fig. 9-1: Functions of the control unit

- | | |
|----------------|----------------------------|
| 1 Display | 4 Joystick button |
| 2 Control knob | 5 Screen-dependent Buttons |
| 3 Joystick | |

Component	Function	Operation
Display (1)	Shows program screens Serves as a touch screen	➔ Lightly press the desired button

Component	Function	Operation
Control knob (2) 	Changes the corresponding Parameters Enables the selected parameters	<ul style="list-style-type: none"> ➔ Turn the knob to the left or right. The selected parameter is highlighted in blue. ➔ Press the control knob downwards. The selected parameter is enabled or disabled.
Joystick (3)	Sets the height, distance and direction to the left and right	➔ Move the joystick forward, back and sideways, turn it, <i>"Fine Adjustment"</i> page 45.
Joystick button (4)	Starts the measurement manually (if the eye-tracking function is turned off)	➔ Press the button.
Buttons (5)	Enables the adjacent button field, depending on the associated screen	➔ Press the desired button.

If you do not work with a PC/laptop/netbook, enter the patient data via the control panel and manage the data there. You can also initiate the measurements and view the measuring results there, see [sec. 12, page 56](#).






9.2 Display with Touch Screen

In addition to the screen-dependent buttons, you can also use the buttons on the touch screen. The buttons change depending on the function of the display.

- ➔ Lightly press the corresponding buttons on the touch screen to enable the function.

Buttons on the touch screen

You can use the following buttons in Patient Data Management.

Button	Function
	Change keyboard
	Delete character
	Escape
	Enter
	Return to upper line




You can start the measurement and edit the patient and examination data via the computer's patient data management. This constellation is described below.

If you work without a computer/laptop, you can perform the measurement at the Corvis® ST unit itself and can save the patient and examination data; see [sec. 12, page 56](#).

9.3 Start Patient Data Management

After you have switched on the computer, it first loads the operating system. Depending on the setting, the Patient Data Management opens automatically.

➔ If necessary, press the Corvis® ST icon: .

The Patient Data Management user interface appears




Fig. 9-2: Patient Data Management User Interface

- 1 "Functions" group box
- 2 Previous examinations
- 3 [Delete exam.] button
- 4 [Delete Patient] button
- 5 [Export] button
- 6 [Import] button
- 7 Patient list
- 8 "Patient" group box

If the Windows desktop appears, you have to start the Patient Data Management from there.

9.4 Importing patient data

You can import patient data from a USB flash drive, see [sec. 12.5.2, page 61](#).



Note

Loss of data due to computer viruses
Computer viruses can cause loss of data.

➔ Run a virus check before importing data from the USB flash drive.

➔ Press the button [Import]. The following dialog appears:




Fig. 9-3: "Import" dialog

- | | |
|---------------------------------|-------------------|
| 1 Select the source of the data | 4 [Import] button |
| 2 [...] button | 5 Patient list |
| 3 Previous examinations | |



The options for import and export of data are set as defaults in the "settings" field, see [User guide](#).

➔ Depending on the settings you may not have to perform all of the following steps (e.g. selection of the directory).

➔ Select the option (1) which contains the source data ("Folder" or "Single file (U12)").



Recommendation: Import the patient data with the option "Folder".

➔ Press the [...] button. (2)

➔ In the dialog box, select the directory or the file where the patient data is located: .DAT and .BMP.

- ➔ Confirm your selection with [OK] or [Save].
The patients that are located and the associated examinations are displayed in the lower part of the dialog.
- ➔ To import the data, press the [Import] button (4).
The data will then be available in the Patient Data Management software.

9.4.1 Enter a new patient



You can only start a measurement from the computer/laptop if you have entered a patient in the patient data management. If you start the measurement from the computer software, the patient data are transmitted (via USB cable) to the Corvis® ST and are saved there. After the measurement with the Corvis® ST has been completed, the examination data for the patient in question are automatically saved to the computer/laptop.

Manually enter a new patient

You can also manually enter a new patient.

- ➔ Press the [New] button to enter a new patient in the Patient Data Management software.
- ➔ Enter the patient's full last name, first name and date of birth in the patient window.



Fig. 9-4: Entering patients

- ➔ Optionally you can enter an ID number for the patient.
- ➔ To save the data you entered, click [Save].
The patient you have just entered now appears in the patient list.

9.4.2 Select an existing patient

The patient list on the left-hand side of the screen displays all previously examined patients in alphabetical order.

Patient list:

	Lastname	Firstname	Date of birth
▶	00_Demo_Corvis	Stable	11.11.1975
	00_Demo_Corvis	Pre/post CXL	24.07.1988
	00_Demo_Corvis	Normal	18.07.1983
	00_Demo_Corvis	Keratoconus	02.05.1996

Fig. 9-5: Patient list

- ➔ Press the [Search] button to quickly find the patient you require in the list.
- ➔ Enter the patient's name or the first letter of the name in the "Last name" field.
Optionally, you can search through the patient's ID number, first name or date of birth, if it was assigned when you first entered the patient.
- ➔ Click the appropriate entry in the list to transfer that patient's name to the patient window. This also brings up a list of any previous examinations for that patient in the examination window (bottom right side).

Extended patient search: [Extended] checkbox

- ➔ Click on the Checkbox [Extended].

The screen displays additional search parameters which, for example, reference previous examinations. Proceed as you did when entering the patient name.

Patient:

Last name:	<input type="text"/>		
First name:	<input type="text"/>		
Det. o. B.:	<input type="text" value="15.08.2019"/>	<input type="text" value="15.08.2019"/>	<input type="button" value="Search by"/>
ID Number:	<input type="text"/>	<input type="text"/>	<input type="button" value="Cancel"/>
Exam date:	<input type="text" value="15.08.2019"/>	<input type="text" value="15.08.2019"/>	<input checked="" type="checkbox"/> Extended
Device:	<input type="text"/>		
Eye:	<input type="text"/>		
Type of exam.:	<input type="text"/>		
Infrated:	<input type="text"/>		
Result:	<input type="text"/>		
	<input type="checkbox"/> Invert result		

Fig. 9-6: Advanced search

10 Perform a Measurement

Perform a measurement for each eye before you save the data. To perform a measurement, you must

1. Start the Corvis® ST program at the computer, [sec. 10.1, page 41](#)
2. Adjust the Corvis® ST, [sec. 10.2, page 42](#)
3. Start the measurement, [sec. 10.3, page 47](#)
4. Perform a second measurement, [sec. 12.3, page 60](#)
5. Save the data, [sec. 10.5, page 48](#)
6. Complete the measurement, [sec. 10.6, page 48](#)

You can use the patient data management, [sec. 10.7, page 49](#)

10.1 Start the Corvis® ST Program at the Computer

- ➔ Select a patient's name.
- ➔ Press the [Corvis® ST] button to start the Corvis® ST program.



Fig. 10-1: Start Corvis® ST program



Condition: Autostart must be activated in the Corvis® ST software settings. See [User Guide](#).

10.2 Adjusting the Corvis® ST

Before starting a measurement, adjust the Corvis® ST.



Attention

Risk of contact of the patient's eye with the air nozzle

→ Before starting a measurement, make sure the safety stop is set correctly, [sec. 8.5.2, page 32](#). This prevents the air nozzle from touching the patient's eye.

Danger of pinching to hands or body parts

→ During a measurement: Make sure that the patient does not put either his hand or any other body part between the measure head and adjusting base.



Note

Incorrect measurements from a dirty air nozzle

→ Before each measurement, check the glass part of the air nozzle from a diagonal angle for dust, dirt, etc.

→ If necessary, clean the air nozzle "[Clean the air nozzle \(2\)](#)" [page 65](#).

10.2.1 Rough Adjustment

- Check that
 - fresh paper has been put onto the chin rest or that the chin rest has been cleaned.
 - that the head rest has been cleaned and disinfect the head rest after each examination, [sec. 13, page 63](#).
- Do not touch the Corvis® ST (including accessories) and the patient simultaneously.
- Ask the patient to place his or her head on the chin and head rest.




Fig. 10-2: Position patient according to the markings

- 1 Head rest
- 2 Mark eye height
- 3 Chinrest

The eye height marking (2) between the chin rest and the headrest should be located roughly at the centerline of the patient's eye.



Attention

Risk of contact of the patient's eye with the air nozzle
Fast and/or uncontrolled movements may cause the air nozzle to touch the patient's eye.

- When operating the Corvis® ST, move it carefully toward the patient's eye.
- If necessary, set the safety stop again, see [sec. 8.5.2, page 32](#).



- If necessary, adjust the height of the chin rest with the buttons. In addition, you can also adjust the height of the measure head by turning the joystick:
Turn it clockwise to move the measure head upwards.
Turn it counter-clockwise to move it downwards, see [sec. 10.2.2, page 45](#).
- To prepare the patient prior to the measurement:
Explain to the patient what will happen next to help him/her relax:
"A little air will be blown into your eye. Don't let it scare you. Please be patient and relax for a moment."




→ Ask the patient to not blink during the measurement, otherwise the measurement results will be falsified.

→ Move the adjusting base until the patient's eye is focused on the display.

→ If necessary: Adjust the height with the buttons.



The auto-tracking function automatically looks for the trigger position.



1 Auto-tracking is activated

Fig. 10-3: Settings 1

→ Proceed as described in [sec. 10.3, page 47](#).

If the auto-tracking function is deactivated, you must make a fine adjustment, [sec. 10.2.2, page 45](#).

10.2.2 Fine Adjustment

- ➔ Use the information on the display and the joystick to make any fine adjustments. Move the joystick in the direction indicated.




Fig. 10-4: Fine Adjustment

- Example
- ➔ Move the joystick to the right.
 - ➔ Turn the joystick clockwise.

Arrow	Camera movement	Joystick movement ^{aa}
➔	right	Move the joystick to the right
➜	left	Move the joystick to the left
⬆	forward	Move the joystick toward the patient
⬇	back	Move the joystick away from the patient
↻	up	Rotate the joystick clockwise
↺	down	Rotate the joystick counter-clockwise

a. If you turn the joystick to its limit, the measuring head and the chin rest move in the opposite direction.

When the position has been reached accurately enough, a cross appears in the center of the ring that is bordered by four bars.





Fig. 10-5: Final position reached

During the measurement process, you will work with the following display:



- 1 Display measuringmode
- 2 Camera image
- 3 Examined eye
- 4 Buttons

Fig. 10-6: Measurement procedure

Element	Function
Displaying mode(1)	Information about the displaying mode, automatically active
Camera image (2)	Camera image of the patient's eye
Examined eye (3)	The examined eye is displayed: left (L) or right (R)
Button (4)	
	Measure mode, automatically enabled
	Back to Patient Data Management
	Adjust height

10.3 Triggering the measurement

A measurement is performed for each eye. The measurement results are then transmitted to the computer/laptop.

Depending on what has been preset (*see User Guide*), the measurement is now either triggered automatically or you must trigger the measurement manually.

Automatic measurement

When the triggering position has been reached (*fig. 10-5, page 45*), the Corvis® ST automatically triggers the measurement.

Manual measurement

- ➔ Trigger the measurement manually by pressing the joystick button (*fig. 9-1, page 34, Item 4*).

10.4 Perform a second measurement

- ➔ Position the patient so that his other eye is in front of the Corvis® ST.

Depending on what has been preset, the measurement is now either triggered automatically or you must trigger the measurement manually.

After the Measurements

- ➔ Press the „Finish“ button.




Fig. 10-7: End examination

The measurements are calculated.

You must wait for this process to end before you can save the measurement results (*sec. 10.5, page 48*) or display the measurement results on the computer (*sec. 11, page 53*). The IOD and pachymetry values are displayed in advance on the device.

10.5 Saving data

Depending on the installation, your data is transferred to a computer or laptop. You can use the patient data management for your data, [sec. 12, page 56](#).



Note

Loss of data due to interrupted save

If you turn off the Corvis® ST while the progress bar is still displayed, data will be lost.

➔ Do not turn off the power until saving is completed, e.g. when the progress bar is complete.



If you have entered the patient data at the Corvis® ST unit itself, the measurement is not automatically saved to a file. You can use a USB stick for this purpose, [sec. 12.5.2, page 61](#).

You can set the data transfer speed from the Corvis® ST to the computer/laptop in the device settings ([see User Guide](#)).

10.6 Complete measurement

With the Corvis® ST program

With this menu item, you can select other patients, or can exit the Corvis® ST program. You then return to the patient data management..




Fig. 10-8: Completing a measurement

➔ Press the [New Patient / Quit] button.

At the Corvis® ST device

- ➔ After each patient remove one of the paper sheets from the chin rest. See also [sec. 13.4, page 67](#).
- ➔ Disinfect the head rest and, if necessary, the chin rest after each patient, [sec. 13.2, page 66](#).



Attention

Risk of infection after examining a sick patient

If you perform a measurement on a sick patient, the air nozzle and the front cover can be contaminated.

- ➔ Clean the air nozzle when you have examined a sick patient, see ["Clean the air nozzle \(2\)" page 65](#).
 - ➔ Disinfect the front cover [sec. 13.2, page 66](#).
-

10.7 Using Patient Data Management

When you have completed an examination, you can continue processing the patient data in Patient Data Management

You can do the following to the patient data

- rename, [sec. 10.7.1, page 49](#)
- export, [sec. 10.7.2, page 50](#)
- save, [sec. 10.7.3, page 51](#)

You can also change the Patient Data Management settings, [User guide](#).

10.7.1 Rename Patient Data


Patient data can be changed retroactively after it has been added.

- ➔ Press the [Change] button.
The input boxes for patient data are now enabled, and the cursor jumps to the "Last name" field.
- ➔ Change the entries in the individual boxes.
- ➔ Press the [Save] button.

10.7.2 Exporting Patient Data

For example, patient and examination data can be exported for forwarding to another clinic.

- ➔ Select the patient and also one of the examinations in the respective list.
- ➔ Press the [Export] button below the patient list. The following dialog appears:



- 1 Select destination
- 2 [...] button
- 3 [Cancel] and [Export] buttons

Fig. 10-9: "Export patient data" dialog



The options for import and export of data are set as defaults in the "settings" field, see [User guide](#).

Depending on the settings you may not have to perform all of the following steps (for example selection of the directory).

- ➔ Select the "Target" (1) where you would like to export the data.



Recommendation: Export the patient data using the "Individual file (U12)" option.

- ➔ Press the [...] button. (2)
- ➔ In the dialog that appears, select the folder or the file to which the patient data should be exported.
- ➔ Confirm your selection with [OK] or [Save].
- ➔ To export the data, press the [Export] button (3).

10.7.3 Data backup

You should carry out a backup of patient and examination data at regular intervals. In case of a loss of data, you can reconstruct the data from a previously created backup with the help of this function. Since data backup takes several minutes depending on the scope of the database and the data to be backed up, a backup should be carried out when the computer and the device will not be needed.



Note

Loss of data due to computer viruses

Computer viruses can cause loss of data.


➔ Run a virus check before importing data from the USB flash drive.



The general rules for security backups apply to backup copies created with the help of the Patient Data Management user interface. Storage of backup files should always be done on a separate system (for example on a USB flash drive with adequate capacity).

10.7.4 Backup data

➔ Press the [Backup] button on the upper right part in the Patient Data Management user interface. The following dialog appears:



- 1 [Restore] button
- 2 Display free storage space.
- 3 Backup folder and button [...]
- 4 [Save] button
- 5 Backup data selection

Fig. 10-10: "Backup" dialog

➔ Select whether all of the data or only changed data should be backed up.



The Patient Data Management function internally tags all saved data records.

If you selecting the option "Changed and new data only", only the data records that were not saved during a previous backup will be backed up.

- ➔ Press the [...] button to the right of the "Backup directory" box (3).
- ➔ In the dialog that appears, select the folder to which the data should be backed up.
- ➔ Confirm your selection with [OK].
- ➔ To back up the data, press the [Save] button (4). The previously selected data will then be backed up to the corresponding folder.

10.7.5 Reconstruct data

If a loss of data occurs, the data from a previous backup can be re-imported into the Patient Data Management user interface.

- ➔ Press the [...] button to the right of the "Backup directory" box (3).
- ➔ In the dialog that appears, select the folder which contains the backup data.
- ➔ Confirm your selection with [OK].
- ➔ To import the data, press the [Restore] button (1). All data in the appropriate directory are copied to the Patient Data Management software.


10.7.6 Automatic backup

In addition to the manually performed backup, it is also possible to automatically run a backup when exiting the Patient Data Management system. The settings required for this can be made in the "Settings" area, see [User Guide](#).

11 Working with the Corvis® ST program

If you started the Corvis® ST program by selecting an examination in the Patient Data Management program, the selected examination is loaded in the Corvis® ST program.

If not, you may need to load an examination first, see [sec. 9.3, page 37](#).



- 1 Patient and examination data
- 2 Curv chart
- 3 Bar chart
- 4 Video function
- 5 Menu bar

Fig. 11-1: Overview Screen; example “Vinciguerra Screening Report”



The most important controls of the Corvis® ST program are described in this section. You will find a detailed description in the [User Guide](#).

11.1 View patient and examination data [1]

The patient and examination data are displayed on every screen of the Corvis® ST program.


11.2 Working with the Menu Bar [7]

You can access the menu bar from any screen of the Corvis® ST.



Fig. 11-2: Menu bar Corvis® ST Program

11.2.1 Menu item "Patient"




End Corvis® ST program, load new patient/examination

The Corvis® ST program ends and switches to Patient Data Management. There you can download a new patient/examination

➔ Click [New Patient / Quit].


11.2.2 Menu item "Examination"



Starting a new examination

If an examination is not automatically initiated, you can initiate a new examination here.

➔ Click [New Examination].




Load an existing examination

You can load an existing examination.

➔ Click [Load Examination].

For more information, refer to the *User Guide*.



Add a comment

You can add or change information and comments.

➔ Click [Add Infotext].

The dialog box "Load Examination" appears.

➔ Enter your text and confirm with [OK].





Fig. 11-3: Enter info text

The new text appears in the Patient Data Management in the list of previous examinations.

11.2.3 Menu item "Display"




Display different presentations

You can display different presentations of an examination.

➔ Press the appropriate button.

For more information about the menu item „Display“, refer to the *User Guide*

11.2.4 Menu item "Export"




You can export different data relating to the examinations, e.g. individual videos, images, or screenshots.

➔ Press the appropriate button.

For more information, refer to the [User Guide](#).

11.2.5 Menu item "Settings"



Change settings


You can change some of the settings, e.g. contrast and brightness, to get a better image quality.

➔ Press the appropriate button.

For more information, refer to the [User Guide](#).

➔ Select the menu item [Miscellaneous Settings] and then click on the desired item.

11.3 Use video function [4]



The video shows the entire deformation of the cornea.

➔ Press this button (1) to start the video.

➔ Press this button (2) to stop the video.

➔ Press this button (3) to return to the start of the video.

➔ Press this button (4) to repeat the video.

11.4 View specific images of the cornea [4]



You can view the images of the cornea at specific times; A1 (Applanation 1), HC (Highest Concavity) and A2 (Applanation 2).

➔ Press the appropriate button.






12 Perform a measurement with the Corvis® ST unit

You can perform an examination with the Corvis® ST unit alone and can directly save the patient and examination data there.

The procedure for this is as follows:

- ➔ Switch on the Corvis® ST, [sec. 8.5.1, page 32](#).
- ➔ Enter the patient data, [sec. 12.2, page 57](#)
- ➔ Prepare and execute the measurement, [sec. 12.2.2, page 59](#)
- ➔ End the measurement, [sec. 12.4, page 60](#)
- ➔ Re-use the patient and examination data, [sec. 12.5, page 61](#)

Buttons and Icons

Button	Function
	This button appears after the measurement to save examination data of the patient
	Switching to the "Dynamic Corneal Response" (DCR) screen
	Switching between eyes/displays
	Printing measurement results ➔ Press this button to start printing
	Measurement is available, displaying the measurement ➔ Press this button to display further measurements
R or L	Display of the examined eye

12.1 Change Settings



You can change the default settings for your own individual measuring mode. For detailed information, refer to the [User Guide](#).

12.2 Entering Patient Data



If the Corvis® ST program is started with the computer, you do not have to enter patient data.

When you turn on the Corvis® ST, you see the Patient Data Management first.




Fig. 12-1: Corvis® ST switch on

Use Patient Data Management to associate examinations with a patient or when you want to save them long-term.



➔ In these cases, it is best to enter the patient's name and date of birth before you conduct the measurement.

12.2.1 Entering new Patient

- ➔ Press the [New Patient] button in the patient data menu to enter a new patient.

The following screen appears:




Fig. 12-2: Touch-screen keyboard, enter patient data

- ➔ Use the touch-screen as described in (sec. 9.2, page 36).
- ➔ Enter the patient's last name and confirm with [Enter].
- ➔ Enter the first name. Confirm by pressing [Enter].

In the "D.o.Birth" field, the touch-screen keyboard changes to a numeric keypad:




Fig. 12-3: Touch screen keyboard, numeric keypad

- ➔ Enter the date of birth and confirm with [Enter].



If you have entered the correction calculation per „Spoerl“ you must enter the patient's date of birth so that the correction is calculated correctly.

You will be asked if you want to save the new patient data.

➔ Select "Yes".

The patient's name appears in the list.




Fig. 12-4: Patient list

➔ Press the [Start] button to switch to measure mode, [sec. 10.1, page 41](#).

12.2.2 Select an existing patient

Choose a patient whose data is already stored, and whom you wish to examine again.

- ➔ In the Patient Data Management ([fig. 12-1, page 57](#)) menu, press the [Patient List] button.
- ➔ Turn the control knob to the desired entry in the list.
- ➔ Press the control knob to select the desired patient.
- ➔ Press the [New Exam/Back] button to switch to measure mode, [sec. 10.1, page 41](#).

12.3 Perform a Measurement

To prepare the measurement, proceed as described in [sec. 9.1, page 34](#) and [sec. 9.2, page 36](#).

Perform the measurement as follows:

- Select the measuring mode, [sec. 12.3.1, page 60](#)
- Adjust the Corvis® ST, [sec. 10.2, page 42](#)
- Start the measurement, [sec. 10.3, page 47](#)

12.3.1 Selecting a Measurement Mode

Depending on the version you can select the measurement mode.



→ Press this button.

Tono + P version

→ Press the [Tono/Pachy] button. Tonometry and pachymetry are measured simultaneously.

Pachy version

→ Press the [Pachy] button.

→ The corneal thickness along the horizontal sectional plane is measured but without the IOP.

12.4 Complete measurement



This button is displayed after a measurement has been conducted.

- Press this button to save the examination data to the patient's record.
- After each patient remove one of the paper sheets from the chin rest. See also [sec. 13.4, page 67](#).
- Disinfect the head rest and, if necessary, the chin rest after each patient, [sec. 13.2, page 66](#).



Attention

Risk of infection after examining a sick patient

If you perform a measurement on a sick patient, the air nozzle and the front cover can be contaminated.

- Clean the air nozzle when you have examined a sick patient, see "[Clean the air nozzle \(2\)](#)" [page 65](#).
- Disinfect the front cover [sec. 13.2, page 66](#).

12.5 Re-use examination data

12.5.1 Delete a Patient or an Examination

If you want to delete a patient or an examination:

- ➔ Select the patient in question.
- ➔ Press the [Pat./Exam Delete] button.




Fig. 12-5: Delete a Patient or an Examination

To delete a patient:

- ➔ Press the [Patient] button.
The patient will be deleted.

To delete an examination:

- ➔ Select the examination that is to be deleted.
The line for the selected examination appears highlighted in blue.
- ➔ Press the [Examination] button.
The examination will be deleted.

12.5.2 Re-using data with the USB flash drive [optional]

If your Corvis® ST is not connected to a computer, your data is stored on a USB flash drive that is connected on the inside of the unit.



- ➔ Press this button. The measurement results are saved with the appropriate patient's data.

You can re-use this data on a computer.

Prerequisite: You must have the Corvis® ST program, and Patient Data Management software installed on your computer/laptop.

There are .DAT and .BMP data records on the USB flash drive. You can import this data with the Patient Data Management program.

- Open the cover with the display.




Fig. 12-6: Open the cover with the display

- Pull out the USB flash drive.
- Plug the USB flash drive into your computer.
Now you can import the data into the Patient Data Management program, [sec. 9.4, page 38](#).

13 Cleaning, Disinfection and Maintenance

This chapter describes how to clean, disinfect, and maintain the Corvis® ST.

Sterilization is not required.

- Heed the product descriptions and instruction manuals of products and equipment you use to care for, clean, and disinfect the unit and/or its accessories.
- Do not clean the Corvis® ST with aggressive, chlorinated, abrasive, or harsh cleansers.



Note

Equipment damage due to moisture

- Make sure that no liquid can get into the Corvis® ST.

13.1 Cleaning



Attention

Risk of electric shock if the Corvis® ST is not completely disconnected from the mains for these jobs.

- Turn the Corvis® ST off, *sec. 8.5.4, page 33*.
- Unplug the power cord before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.


Required materials

- Cleaner for plastic surfaces with anti-static effect
- Cleaner for painted surfaces: mixture of equal parts of alcohol and distilled water, possibly with a few drops of commercial detergent
- Soft, lint-free cloth (for example gauze or microfibre)
- Cleaning alcohol
- Commercial clean agent for the glass
- Cotton swab
- if necessary: bellows

Cleaning intervals

- Clean the components of the Corvis® ST once a month or if necessary.

13.1.1 Clean the front panel



- | | |
|--------------|-------------------------|
| 1 Head rest | 3 Lens protection glass |
| 2 Air nozzle | 4 Chin rest |

Fig. 13-1: Clean lens protection glass and air nozzle

13.1.2 Clean head [1] and chin [4] rest



The Corvis® ST can remain switched on for this cleaning step.



During the measuring process, sweat, cosmetics, etc. from the patient can get on the head and chin rest.

- ➔ Clean these parts before examining the next patient.
Before and after each measurement wipe off the head and chin rest with a clean cloth, for example, with gauze moistened with rubbing alcohol.

Do not wipe more difficult spots repeatedly with a dry cloth. Instead moisten it with rubbing alcohol.

13.1.3 Clean lens protection glass [3]

The openings in the housing for the optics are covered by protective glass covers which must be kept dust and dirt-free.

- Clean the lens protection glass with a lint-free cloth moistened with alcohol.

13.1.4 Clean the air nozzle [2]



Attention

Risk of infection after examining a sick patient

If you perform a measurement on a sick patient, the air nozzle can be contaminated.

- Moisten a cotton swab with rubbing alcohol.
- Clean the nozzle with the cotton swab.
- Disinfect the device as needed, see [sec. 13.2, page 66](#).

- Check the glass part of the air nozzle from a diagonal angle for dust, dirt, etc.
- Blow off dust, foreign particles, etc. with a lens blower.
- Afterwards, carefully wipe the glass with a cotton swab which you have previously moistened with methanol or pure alcohol.



Note

Damage due to improper cleaning

- Carefully wipe the air nozzle and do not rub the surface.
- Make sure there are no foreign particles on the air nozzle when wiping it. Otherwise, the glass could get scratched.

- Check the glass afterwards.

13.1.5 Clean the housing

Wipe the head rest after each examination, and the housing as required.

- Turn the Corvis® ST off, [sec. 8.5.4, page 33](#).
- Unplug the power cord. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.
- When cleaning, use a damp cloth and make sure that no liquid enters the Corvis® ST.
- Clean the lens in front of the camera using a dry, lint-free cloth.
- Clean the plastic surfaces and painted surfaces with the appropriate cleaning agents.

13.1.6 Clean the display

- ➔ Clean the display using a dry, lint-free cloth.

13.2 Disinfection

Required materials:

- Disinfection and Cleaning Kit (included),
Alternatively: Pursept®-A Xpress disinfectant wipes,
Merz+Co company; D-60318 Frankfurt:
Tel: +49 69 1503 1; Fax: +49 (69) 596 21-50
E-mail: merzpr@merz.de
- ➔ Turn the Corvis® ST off, [sec. 8.5.4, page 33](#).



1 Head rest

2 Front cover

Fig. 13-2: Disinfect

- ➔ Disinfect the head rest and, if necessary, the chin rest after each examination.



Attention

Risk of infection after examining a sick patient

If you perform a measurement on a sick patient, the head rest and the front cover can be contaminated.

- ➔ Disinfect the front cover (2) and the housing as required.



Note

Equipment damage due to disinfectant solution

The disinfectant solution may damage the finish if it is sprayed directly on it.


- ➔ Only spray the disinfectant solution onto a cleaning cloth, not directly on the device

13.3 Maintenance

- ➔ Germany: The operating company must ensure that the device undergoes technical measurement testing every 2 years according to MPBtreibV, Appendix 2 Tonometer.
- ➔ To ensure that it functions correctly and safely we recommend the following: Have the Corvis® ST checked every two years by our service department or an authorized dealer.

13.4 Attach paper to the chin rest

If you want to attach new chin rest paper, follow these instructions:



- 1 Pins
- 2 Chin rest paper
- 3 Chin rest

Fig. 13-3: Attach chin rest paper

- ➔ Pull the two pins (1) out of the chin rest.
- ➔ Put the chin rest paper (2) in such a way that the holes of the paper and the chin rest (3) are aligned.
- ➔ Insert the two pins (1) in the chin rest.

13.5 Insert new printing paper roll

➔ Open the cover with the display.



Fig. 13-4: Open the cover with the display

The following screen appears:





Fig. 13-5: Change printer paper

You can advance and reverse the printer paper by pressing the buttons "Printer Feed" and "Feed Back" accordingly.


➔ Press "Feed Back" to reverse, or roll back the printer paper.



- 1 Printer roller
- 2 White slot
- 3 Metal pin
- 4 Holder

Fig. 13-6: Remove and insert printer roller

- ➔ Remove the feed roller (1) from the holder (4).
- ➔ Pull the metal pin (3) out.
- ➔ Push the metal pin into a new printer roller and insert the printer roller into the holder (4).
- ➔ Slide the white paper through the white slot from below.
- ➔ Make sure the paper (1) is correctly aligned.



1 Proper paper guide

2 Wrong paper guide


Fig. 13-7: Insert the paper

- ➔ Press the button "Printer Feed" so that the printer paper is pulled through the opening.
- ➔ Close the cover with the display, see [fig. 13-4, page 68](#).

14 Dismantling, Transport and Storage

Before you transport or store the device you may have to dismantle it properly and lock the transport safety device. To avoid transport damage, follow the steps in the sub-sections below.

14.1 Parking position



- ➔ Switch on the Device with the On/Off Switch (3).
- ➔ Press Settings
- ➔ Navigate to the Settings page 5/5
- ➔ Chose the option Parking position.





Fig. 14-1: Settings 5/5

The Device turns into the Parking position.

14.2 Lock the transport safety device



- ➔ Switch off the Corvis® ST with the On/Off Switch (3).
- ➔ Unplug the power cord (2) of the device.
- ➔ If necessary, disconnect the USB cable from the USB port of the computer/Laptop (1).
- ➔ Open the display cover.





Fig. 14-2: Open the display cover

- ➔ Move the Corvis® ST over the holder (1) of the transport safety device in the adjusting base.
- ➔ Lock the transport safety device (2).
Press down gently on the transport safety device and turn it counter-clockwise to the "unlocked" position (3). The transport safety device must be engaged!.



1 Holder of the transport lock

2 Spring

3 "Locked" position

Fig. 14-3: Transport safety device

- ➔ Close the display cover, [fig. 14-2, page 71](#).

14.3 Joystick lock

- ➔ Turn the Handle (1) into the right end Position..



Abb. 14-4: Joystick Arretierung

The Device is now fully prepared and ready to be packed.

14.4 Transport and Storage Information

This device can withstand the following temperature conditions for storage and transport.

Storage

Ambient temperature range	-10°C to +55°C
Relative humidity, including condensation	10% to 95%
Air pressure range	700 hPa to 1060 hPa

Transport

Ambient temperature range	-40°C to +70°C
Relative humidity range, including condensation	10% to 95%
Air pressure range	500 hPa to 1060 hPa

14.5 Transport and Storage




Note

Risk of equipment damage due to incorrect shipment or from improper storage

If you lift the appliance by the headrest, it may break off.

- ➔ Grip the device from below to lift it.
- ➔ Avoid shocks, vibrations, and contamination.
- ➔ Avoid high temperatures and humidity.
- ➔ Check the device for damage after each transport
- ➔ Do not hold the device by the joystick to carry it.
- ➔ Do not operate the appliance for approx. 3-4 hours after transportation or storage.
- ➔ Strong temperature changes from cold areas to warm rooms can cause the optical components to fog up.

15 Disposal of Used Devices



On accordance with Directive 2012/19/EU of the European Parliament and the Council of 4. July 2012, and in accordance with German law governing the circulation, return and environmentally friendly disposal of used electrical and electronic devices, such appliances must be recycled and may not be discarded as household waste.

16 Troubleshooting



Attention

Persons or equipment damage due to incorrect troubleshooting

- Do not plug in or unplug any cables while the Corvis[®] ST is switched on.
- If an error occurs which you are unable to correct by following the instructions below, label the device as "out of order" and contact our service department or an authorized dealer.

Fault	Possible cause	Help
No function when the On/Off switch is pressed	Corvis [®] ST not connected to the power supply. Power failure or power outlet is not active.	<ul style="list-style-type: none"> → Plug power cord into the outlet or into the connector on the Corvis[®] ST. → Inform the in-house electrician. → Check that the connector is plugged in properly.

17 Terms of Warranty and Servicing

Any software included in the delivery was tested by us and complies with technical standards. Note the following warranty provisions:

- Prior to and while operating the device it is important that you observe the instruction manual and safety instructions.
- The Corvis® ST carries a warranty to which you are entitled in accordance with the legal provisions.
- If any unauthorized persons interfere with the Corvis® ST, all warranty entitlements shall be void. Any inappropriate modifications or repairs can cause grave danger to the user and patient.
- Any entitlement to a warranty shall also be void if unauthorized persons interfere with the computer hardware and supplied software.
- Make transport damage claims to the shipping company during or immediately after delivery. Have the damage confirmed on the bill of lading, so that a proper claim settlement is possible.
- In general, the general terms and conditions of business and delivery apply as per the date of purchase.

17.1 Liability for proper function or damages

OCULUS will only accept responsibility for the safety, reliability and serviceability of the Corvis® ST if the unit is used in compliance with the following terms:

- Use the device in accordance with these instructions and the accompanying user manual.
- There are no user-serviceable parts either on or inside the Corvis® ST. OCULUS shall not assume any liability if assembly, extensions, adjustments, changes or repairs are carried out by unauthorized personnel; if the Corvis® ST is maintained improperly; or if it is handled incorrectly.
- If the work described above is carried out by persons authorized to do so, they must be required to supply a certification detailing the nature and scope of repairs, and, if applicable, to specify modifications to the rated data and area of work. The certificate must bear a date, a signature, specify who carried out the work, and contain company information.
- On request, and for this purpose, OCULUS will supply authorized persons with spare parts lists and additional descriptions.
- Make sure that only original OCULUS parts are used for service and maintenance.

18 Technical Data

Measure Mode

IOP + pachymetry

Tonometer

Measuring range	6 to 60 mmHg
Working distance	11 mm
Internal fixation light	red LED

Scheimpflug camera

Frame Rate	4330 frames/s
Measuring range	8.5 mm horizontally
Pachymeter measuring range	300 to 1200µm
Measuring points	576 per image
Display resolution	576 x 200 pixels
Light source	blue LED (470 nm, UV-free)

Classification according to IEC 60601 - 1

Type of protection against electrical shock	Protection class 2
Level of protection against electrical shock	Type B
Protection level of the housing against intrusion of foreign objects and liquids	IP20

Operating conditions

Temperature	+10 — +35 °C
Humidity	30 — 90 %
Air pressure	800 — 1060 hPa

Measuring head

Power supply	15 V DC; 6 A
Max. power consumption	26 W

Power adapter

Mean Well (05150285)	
Mains connection	100 - 240 V
Frequency	50/60 Hz
Max. power consumption	max. 90 W
Output	15 V DC

Computer.

Recommended computer specifications (for optional software)	Intel® Core™ i5, 500 GB SSD, 8 GB RAM, Windows® 10, Intel® HD Graphics
The IT equipment (computer, monitor, etc.) must meet the requirements of IEC 62368-1 or IEC 60950	
Recommended screen size	24"
Recommended screen resolution	1920 x 1080 Pixel (Full HD)

Other information

Dimensions W x D x H	266 x 538 x 495 to 525 mm (10.5 x 21.2 x 19.5- to 20.7 in)
Weight	approx. 14 kg (30.8 lbs)
Printer	Thermal printer
Display	TFT - LCD approx. 150 mm
Lifecycle expectancy	up to 10 years

CE in accordance with Regulation (EU) 2017/745 on Medical Devices

The unit is a class IIa product.



Conformity assessment: (EU) 2017/745 MDR, Annex IX, Chapters I and III.

19 Appendix

19.1 Correction Calculation of the Tonometrically calculated IOP

The Corvis® ST program offers the option to save the tonometrically measured IOP in the data of an examination and to correct it by Central Corneal Thickness (CCT). Different correction formulae can be applied for this. The corrected IOP value and the IOP change are also saved.

The correction formulae according to Shah, Ehlers, and the Dresdener correction formulae approximate the actual IOP related to corneal thickness from the measured value.



As the table is imported, a linear relationship will be adopted for the variates. It is sufficient to define a start value and an end value with the associated correction values.

19.2 Electromagnetic Compatibility

Medical electrical equipment is subject to special precautionary requirements with respect to EMC, and must be installed and operated according to the EMC-Instructions contained in the accompanying paperwork.

No special measures need be observed in respect of OCULUS devices and systems.

Portable and mobile RF-communications devices can interfere with electrically operated medical devices.

Produced in consideration of permissible deterioration during or caused by the EMC testing without affecting the essential performance criteria

- A short interruption of the USB connection during the examination is permissible because it will not affect the diagnosis, treatment and observation.



Warning

The use of accessories, transducers and cables not specified by OCULUS (for example as replacement parts) may result in increased emissions or decreased immunity of the Corvis® ST.

- ➔ Use only the original accessories, transducers and cables specified by OCULUS.

The use of accessories, transducers and cables specified by OCULUS with devices other than the Corvis® ST may result in increased emissions or decreased immunity of the other device.

- ➔ Do not use the accessories, transducers and cables specified by OCULUS with devices other than the Corvis® ST.
-

To be in compliance with the requirements of the IEC 60601-1-2 the following types of equipment, accessories, power adapters and cables must be used.

Order number	Description	
72100	Corvis [®] ST	
05200905	Cable with plug, EU standard	1.8 m
05200910	Cable with plug, US standard	1.8 m
05200915	Cable with plug, GB standard	1.8 m
05200920	Cable with plug, AU standard	1.8 m
05200601	USB cable	
015692000010	USB FS MED Isolator	
05150285	Power Supply	15 V DC; 6 A

19.3 Guidance and Manufacturer's Declaration - Electromagnetic Emissions and Immunity for the Corvis® ST


Guidance and manufacturer's declaration electromagnetic emissions IEC 60601-1-2		
The OCULUS Corvis® ST is intended for operation in the electromagnetic environment specified below. The user of the Corvis® ST should ensure that it is being used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Corvis® ST uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
HF-emissions CISPR 11	Class B	
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	complies	

Electromagnetic immunity, IEC 60601-1-2			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Electrical Fast transient/bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines	± 2 kV ----- ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 6100-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_{τ} : 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree	0% U_{τ} : 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Corvis [®] ST requires continued operation during power mains interruptions, it is recommended that the Corvis [®] ST be powered from an uninterruptible power supply or battery.
	0% U_{τ} : 1 period and 70% U_{τ} : 25/30 periods Single-phase: at 0 degree	0% U_{τ} : 1 period and 70% U_{τ} : 25/30 periods Single-phase: at 0 degree	
	0% U_{τ} : 250/300 periods	0% U_{τ} : 250/300 periods	

Note: U_{τ} is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity, IEC 60601-1-2

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
Conducted RF IEC 61000-4-6	3 V _{eff} 150 KHz to 80 Mhz 6 V in ISM- and amateur radio frequency bands between 150 kHz and 80 MHz 80% AM to 1 kHz	V _{eff} = 3 V	Portable and mobile RF communications equipment should be used no closer to any part of Corvis® ST, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz		$d = \left[\frac{3,5}{(V_1)} \right] \sqrt{P}$ $d = \left[\frac{3,5}{(E_1)} \right] \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{(E_1)} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interface may occur in the vicinity of equipment marked with the following symbol:			
			
Note 1:	At 80 Hz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
<p>a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Corvis® ST is used exceeds the applicable RF compliance level above, the Corvis® ST should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Corvis® ST.</p> <p>b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the Corvis® ST

The Corvis® ST is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Corvis® ST can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Corvis® ST as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 Mhz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.80	3.80	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

19.4 Description of the Connection



19.5 Data sheet Power Supply Adapter (05150285)



90W AC-DC High Reliability Medical Adaptor

GSM90B series



■ Features

- Universal AC input / Full range
- 2 pole AC inlet IEC320-C8
- Medical safety approved (2 x MOPP between primary to secondary)
- Suitable for BF application with appropriate system consideration
- Low leakage current <100uA
- No load power consumption<0.15W
- Energy efficiency level VI
- Comply with EISA 2007/DoE,NRCan, AU/NZ MEPS, EU ErP and meet CoC Version 5
- Built-in active PFC function
- High efficiency up to 91%
- Fanless design with -30~+60°C working temperature
- Class II power (without earth pin)
- Protections: Short circuit / Overload / Over voltage / Over temperature
- Fully enclosed plastic case
- LED indicator for power on
- 100% full load burn-in test
- Optional lock type DC plug
- 3 years warranty

■ Applications

- Mobile clinical workstation
- Oral irrigator
- Portable hemodialysis machine
- Breath Machine
- Medical computer monitor


■ Description

GSM90B is a highly reliable, 90W desktop style single-output green medical adaptor series. This product is equipped with a 2-pin (no FG) standard IEC320-C8 power plug, adopting the input range from 80VAC to 264VAC. The entire series supplies different output voltages between 12VDC and 48VDC that can satisfy the demands for various kinds of medical electrical devices. The circuitry design meets the international medical standards (2*MOPP), having an ultra low leakage current (<100uA), fitting the medical devices in direct electrical contact with the patients.

With the efficiency up to 91% and the extremely low no-load power consumption below 0.15W, GSM90B is compliant with USA EISA 2007/DoE, Canada NRCan, Australia and New Zealand MEPS, EU ErP, and meet Code of Conduct (CoC) Version 5. The supreme feature allows the adaptor to save the energy when it is either under the operating mode or the standby mode. The entire series utilizes the 94V-0 flame retardant plastic case, providing the double insulation that effectively prevents electrical shock. GSM90B is approved with the international medical safety certificates.

■ Model Encoding

GSM90 B 12 -P1M



File Name: GSM90B-SPEC 2016-03-16



90W AC-DC High Reliability Medical Adaptor

GSM90B series

SPECIFICATION

NO.	GSM90B12-P1M	GSM90B15-P1M	GSM90B19-P1M	GSM90B24-P1M	GSM90B48-P1M
SAFETY MODEL NO.	GSM90B12	GSM90B15	GSM90B19	GSM90B24	GSM90B48
DC VOLTAGE <small>Note.2</small>	12V	15V	19V	24V	48V
RATED CURRENT	6.67A	6A	4.74A	3.75A	1.87A
CURRENT RANGE	0 ~ 6.67A	0 ~ 6A	0 ~ 4.74A	0 ~ 3.75A	0 ~ 1.87A
RATED POWER (max.)	80W	90W	90W	90W	90W
RIPPLE & NOISE (max.) <small>Note.3</small>	120mVp-p	150mVp-p	180mVp-p	200mVp-p	240mVp-p
VOLTAGE TOLERANCE <small>Note.4</small>	±5.0%	±5.0%	±4.0%	±3.0%	±2.5%
LINE REGULATION <small>Note.5</small>	±1.0%	±1.0%	±1.0%	±1.0%	±1.0%
LOAD REGULATION	±5.0%	±5.0%	±4.0%	±3.0%	±2.5%
SETUP, RISE TIME <small>Note.6</small>	1000ms, 50ms / 230VAC 1500ms, 50ms / 115VAC at full load				
HOLD UP TIME (Typ.)	20ms / 230VAC 20ms / 115VAC at full load				
VOLTAGE RANGE <small>Note.7</small>	80 ~ 264VAC 113 ~ 370VDC				
FREQUENCY RANGE	47 ~ 63Hz				
POWER FACTOR (Typ.)	PF>0.91 / 230VAC PF>0.95 / 115VAC at full load				
EFFICIENCY (Typ.)	88%	89%	89%	90%	91%
AC CURRENT (Typ.)	1.3A / 115VAC 0.6A / 230VAC				
INRUSH CURRENT (Typ.)	30A / 115VAC 65A / 230VAC				
LEAKAGE CURRENT(max.)	Touch current < 100µA/264VAC				
OVERLOAD	110 ~ 150% rated output power Protection type : Hiccup mode, recovers automatically after fault condition is removed				
OVER VOLTAGE	105 ~ 135% rated output voltage Protection type : Shut down o/p voltage, re-power on to recover				
OVER TEMPERATURE	Shut down o/p voltage, re-power on to recover				
WORKING TEMP.	-30 ~ +60°C (Refer to "Derating Curve")				
WORKING HUMIDITY	20% ~ 90% RH non-condensing				
STORAGE TEMP., HUMIDITY	-40 ~ +85°C, 10 ~ 95% RH				
TEMP. COEFFICIENT	±0.03% / °C (0 ~ 40°C)				
VIBRATION	10 ~ 500Hz, 2G 10min./1cycle, period for 60min. each along X, Y, Z axes				
SAFETY STANDARDS	ANSI/AAMI ES60601-1 / ES60601-1-11, TUV EN60601-1 / EN60601-1-11 approved				
ISOLATION LEVEL	Primary-Secondary: 2xMOPP				
WITHSTAND VOLTAGE	I/P-O/P: 4KVAC				
ISOLATION RESISTANCE	I/P-O/P: 100M Ohms / 500VDC / 25°C / 70% RH				
EMC EMISSION	Compliance to EN55011(CISPR11) class B, EN61000-3-2,3, FCC PART 15 class B, CAN ICES-3(B)/NMB-3(B)				
EMC IMMUNITY	Compliance to EN61000-4-2,3,4,5,6,8,11, EN55024, EN60601-1-2, EN61204-3 medical level, criteria A				
MTBF	405.6K hrs min. MIL-HDBK-217F(25°C)				
DIMENSION	145*60*32mm (L*W*H)				
PACKING	0.45Kg; 30pcs/14.5Kg/1CUFT				
PLUG	See page 3 ; Other type available by customer requested				
CABLE	See page 3 ; Other type available by customer requested				


1. All parameters are specified at 230VAC input, rated load, 25°C 70% RH ambient.




90W AC-DC High Reliability Medical Adaptor

GSM90B series

Derating Curve




Static Characteristics



Mechanical Specification

Case No. GS90A Unit:mm



Plug Assignment

Standard plug: P1M

P1M	
P/N	OUTPUT
CENTER	+

Optional lock type plug: P2S
SWITCHCRAFT S761K plug equivalent

Installation Manual

Please refer to : <http://www.meanwell.com/webnet/search/InstallationSearch.html>

File Name:GSM90B-SPEC 2016-03-16

19.6 Instructions for integration into an IT-Network

The device, in conjunction with the connected computer and the device software running on it, forms a programmable electrical medical system (PEMS) according to IEC 60601-1.

It is essential that you follow the cybersecurity instructions section of "Safety Instructions" in the device instruction manual.

Please note the following information for implementing an integration of the PEMS into an IT-Network:

The purpose of integrating the PEMS into an IT-Network can be:

- Storage and retrieval of the examination data on a local network drive
- Printing
- Data export
- DICOM workflow

Required characteristics of the IT-Network into which the PEMS is to be integrated:

- Prefer wired LAN connection
- IPv4 network
- Fast-Ethernet (100 Mbit/s minimum)

Required configuration of the IT-Network into which the PEMS is to be integrated:

- Storage, printing, data export: File and Printer Sharing for Microsoft Networks (SMB 3.0 or higher - required opened port: 445)
- DICOM Storage Service Class = PACS
- DICOM Worklist Management Service Class (Modality Worklist Server)

Technical specifications of the network connection with the PEMS, including the specifications of data security:

- Refer to the cybersecurity section of "Safety Instructions" in the device instruction manual.
- Refer to device specific DICOM interface description (if applicable)

The intended flow of information between PEMS, the IT-Network and other devices in the IT-Network and the intended routing through the IT-Network

- License handling from local license server to PEMS and vice versa
- Storage and data export to local network storage and loading from local network storage
- Printout to local printer

List of hazard situations that result from the IT-Network being unable to provide the features that are required to fulfill the purpose of integrating the PEMS into to meet the IT-Network:

- Loss of data
- unsuitable data exchange
- Data corruption
- unsuitable temporal data allocation
- unexpected data reception
- unauthorized access to data



Connection of the PEMS to an IT-Network that includes other equipment could result in previously unidentified risks to patients, operators or third parties.

The responsible organization should identify, analyze, evaluate and control these risks.

Subsequent changes to the IT-Network could introduce new risks and require additional analysis.

Changes to the IT-Network include:

- changes in IT-Network configuration
 - connection of additional items to the IT-Network
 - disconnecting items from the IT-Network
 - update of equipment connected to the IT-Network
-

WWW.OCULUS.DE

OCULUS Optikgeräte GmbH

Münchholzhäuser Str. 29 , 35582 Wetzlar • GERMANY

Tel. +49 641 2005-0 • Fax +49 641 2005-255

E-mail: sales@oculus.de • www.oculus.de

G/72100/EN

Lot:

