OCULUS

Hubtische | Ophthalmic Tables

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Preface

The Ophthalmic Table was manufactured and tested under strict quality criteria.

Proper use of this equipment is indispensable for its safe operation. For this reason you should familiarise yourself thoroughly with the content of this instruction manual before operating the unit. In particular, follow the safety instructions.

Minor differences between the pictures presented herein and the equipment actually delivered to you are possible due to technical and version advances.

If you have any questions or require additional information about your equipment, call or fax us. Our Service Team will be happy to assist you.

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Table of Contents

1 Scope of Delivery	43
2 Graphic Symbols	44
2.1 On the Device / Type Plate	44
2.2 On Packaging.....	45
3 Safety Instructions	46
3.1 Pictogram Used in this Manual	46
3.2 Safety Instructions for Use	46
3.2.1 Instructions for Operating Personnel	47
3.2.2 Transport and Storage Instructions	47
3.2.3 Instructions for Setup and Connection.....	47
3.2.4 Patient environment information.....	48
3.2.5 Information about the operation of an ME system	48
3.2.6 Instructions for Operation	48
3.2.7 Further Instructions for Operation with Optional Wheels.....	49
3.2.8 Notes on Maintenance.....	49
3.2.9 Notes on Disassembly and Disposal.....	49
3.2.10 Note on Electrical Safety	49
3.2.11 Electromagnetic Compatibility (EMC) / Cables.....	50
4 Intended Use	51
4.1 Intended Purpose	51
4.2 Intended medical indication.....	51
4.3 Kontraindication.....	51
4.4 Possible side effects.....	51
4.5 Intended users	51
4.6 Patient group	51
5 Device Description	52
6 Setup and Connection	53
6.1 Setup	53
6.2 Operating Conditions.....	53
6.3 Electrical Connection.....	53
6.4 Network Cable Connection.....	54
7 Operation	55
7.1 Turn on and off.....	55

7.2	Up and down movement.....	55
8	Cleaning and Disinfection	57
8.1	Time table for cleaning and disinfection.....	57
8.2	Cleaning	57
8.3	Disinfection	58
9	Trouble shooting.....	59
9.1	Fault table.....	59
9.2	Exchange Fuse	59
10	Transport, Storage and Disassembly	61
10.1	Disassembly	61
10.2	Transport.....	61
10.3	Storage	61
11	Disposal.....	62
12	Warranty Provisions and Service	63
12.1	Warranty Provisions.....	63
12.2	Liability for Faulty Operation or Damage.....	63
13	Technical Data	64
13.1	Ophthalmic tables without electrical unit (universal tables).....	65
13.2	Ophthalmic tables with electrical unit (system tables)	68
14	Appendices	73
14.1	Electromagnetic Compatibility (EMC).....	73
14.2	Guidance and Manufacturers declaration - Electromagnetic Emmissions and Immunity	74
14.3	Description of the Connection	78

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

Ohne Y-Kabel

1 Scope of Delivery

Ophthalmic Table Kit
1. Table top
■ for Universal Table (without electrical unit)
■ for System Table (with electrical unit)
2. Table base
■ with column
■ with or without monitor mounting
3. Power cord
4. Instuction manual

Observe the technical data from page 28 onwards.

We reserve the right to make changes to deliverables as a part of technical improvements.

- ➔ If transport damage is discovered in the shipment, immediately file a complaint at once with the shipping company.
- ➔ Have the damage noted on the freight bill of lading so that an orderly handling of the complaint for damages is possible.

2 Graphic Symbols

2.1 On the Device / Type Plate


Fig. 2-1: Type plate table top - Kit part 1/2 (example)


Fig. 2-2: Type plate table base- Kit part 2/2 (example)

No.	Description	No.	Description
1	Company Logo	7	Disposal in household trash is prohibited
2	Company address	8	Medical device
3	CE Conformité Européenne	9	Max. Input
4	Read Instructions Manual	10	Power supply
5	UDI number consisting of: UDI-DI (Device-Identification), machine readable Matrix-Code and UDI-PI (Product-Identifier)	11	Fuse
		12	Serial number of the device
		13	Device number
6	Protection class	14	Device name
		15	Disposal in household trash is prohibited

Safety note

Do not push!

2.2 On Packaging

Keep dry!

This way up!

Fragile!

Transport

Permitted temperature range: transport

Lagerung

Permitted temperature range: storage

3 Safety Instructions

This section contains a compilation of the most important information of these safety-related information.

- ➔ Carefully read through the Instruction Manual.
- ➔ Keep the Instruction Manual in good condition near the device.
- ➔ Observe the legal regulations with regard to accident prevention.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer (vigilance@oculus.de) and the competent authority of the region in which the user and/or patient is established.

3.1 Pictogram Used in this Manual

Warning

Identifies a potentially dangerous situation which may cause mayhem injury.

Caution

Identifies a potentially dangerous situation which may cause minor injury or damage to property.

Notice

Denotes situations which could result in incorrect findings, user instructions and useful or other important information.

Identifies important information about the product and its use which require special attention.

3.2 Safety Instructions for Use

Warning

Personal injury or property damage due to improper operation

- ➔ Observe the following safety instructions.

Warning

Personal injury or property damage due to equipment modifications that could jeopardize safety

- It is not allowed to make any changes at the Ophthalmic Table. No modifications may be made to this device without the permission of the manufacturer.
If the ophthalmic table is changed, suitable examinations and tests must be carried out to ensure further safe use.

3.2.1 Instructions for Operating Personnel

- Make certain that the Ophthalmic Table is used exclusively by personnel that have the training and practical experience to safely and properly operate the equipment.

3.2.2 Transport and Storage Instructions

- Refer to the notes in *chap. 10, page 61*.

3.2.3 Instructions for Setup and Connection

- Do not use or store the Ophthalmic Table in damp rooms.
- Avoid droplets, condensation and spray water in the vicinity of the Ophthalmic Table and make sure that no liquids can penetrate into the Ophthalmic Table. Therefore, do not place any liquid-filled containers near or on the Ophthalmic Table.
- Only operate the Ophthalmic Table in medically used rooms after it has been set up according to VDE regulations 0100-710.
- Do not operate the equipment included in the deliverable items in explosive environments, in the presence of combustible anaesthetics or volatile solvents such as alcohol, benzene or such.
- Use only a mains cable that meets the requirements of IEC 60227-1, Type 53, min. 0.75 mm² and of IEC 60320-1.
- Position the Ophthalmic Table so that the mains plug is readily accessible. That way, you can easily disconnect it from the power supply for any repairs or maintenance work.
- Do not operate the ophthalmic table in a surgical room.
- Do not connect the electrical plug connectors with excessive force. If a connection is not possible, check whether the plug fits the jack. Make sure that the plug is fully inserted.
If you find damage to the plug connector, have the damage corrected by our Service Department.
- Be sure to provide proper grounding.

3.2.4 Patient environment information

Patient environment is the area where patients can come into contact with any part of a medical electrical equipment (ME equipment) or with another person being in contact with the ME equipment.

In the patient environment, use devices that conform to IEC 60601-1. If a multiple power socket is to be used, or if a device is to be used that does not meet the IEC 60601-1 standard, use an isolating transformer.


Fig. 3-1: Patient environment (Example with Pentacam®)

3.2.5 Information about the operation of an ME system

The Ophthalmic Table and a connected computer form a medical electrical system (ME system) according to IEC 60601-1. If you connect additional devices, such as, for example a printer, those devices become part of the ME system.

- ➔ Make sure that all devices of the ME system meet the requirements of IEC 60601-1 or IEC 60950-1.

3.2.6 Instructions for Operation

- ➔ Before initial operation: Let OCULUS or an authorized dealer train you in the operation of the Ophthalmic Table.
- ➔ Never operate a damaged Ophthalmic Table.

- Only operate the Ophthalmic Table with the original accessories supplied by us and only when the unit is in technically perfect condition.
- Only operate the equipment after you have read and understood the Instruction Manual.
- Do not touch the patient and the Ophthalmic Table simultaneously.
- Do not exceed the maximum load capacity of the Ophthalmic Table (Refer "Technical Data" on page 64).
- Do not lean any objects against the ophthalmic table.
- Do not push the Ophthalmic Table.
- Do not put any heavy objects or the device itself onto the cables.
- Make sure that the power cord does not come into contact with hot surfaces (e.g. heating).
- Do not allow the lift column to move up and down in continuous mode (one minute operation and thereafter do not move the lift column up or down for 30 minutes).

3.2.7 Further Instructions for Operation with Optional Wheels

- Do not exceed the maximum load capacity of max. 67 kg.
- Always adjust the Ophthalmic Table in the lowest position before transport.
- Always push the Ophthalmic Table with its short side first.
- After transport, secure the Ophthalmic Table by locking the brakes on all four wheels.

3.2.8 Notes on Maintenance

- To ensure satisfactory and reliable operation, we recommend having the Ophthalmic Table checked by our Service Department or by an authorized dealer every two years. If an error occurs which you cannot correct, identify the Ophthalmic Table as being non-operational and contact our Service Department.

3.2.9 Notes on Disassembly and Disposal

- When disconnecting the electrical contacts, do not pull on the cable, but rather only on the particular plugs.
- Dispose of the equipment in conformance with legal requirements.

3.2.10 Note on Electrical Safety

Attention

Personal injury or property damage due to improper safety precautions
Connecting the ophthalmic table with non-medical electrical equipment (e.g. data processing equipment) to a medical, electrical system must not result in a level of safety to the patient which is below that prescribed by IEC 60601-1. If the permissible values for shunt currents are exceeded due to the connection, then protective measures must be provided which contain a disconnect device.

- Make certain that connections to non-medical equipment are made properly.
-
- Do not use a multiple power strip to connect the Ophthalmic Table.

3.2.11 Electromagnetic Compatibility (EMC) / Cables

Attention

Personal injury or property damage due to electromagnetic interferences.

Portable and mobile RF sets may interfere with electrically operated medical devices, [chap. 14.1, page 73](#).

- Make certain that portable and mobile RF communications apparatus do not cause any interfering emissions.
- Recommendation: Maintain a minimum distance of 4 m. If the distance is shorter, you must ensure that the Ophthalmic Table functions properly.

4 Intended Use

4.1 Intended Purpose

The ophthalmic table is intended solely for accommodating ophthalmological instruments and devices in order to move them up and down for optimal positioning in front of the test person.

The Ophthalmic Table is only intended for the uses described in these Instruction of Use.

- ➔ Only operate the device with the original accessories supplied by OCULUS and when the unit is in a technically perfect condition.
- ➔ Observe the safety instructions.

4.2 Intended medical indication

There is no medical indication for the ophthalmic table itself.

However, since the ophthalmic table is only intended to accommodate medical devices, the medical indication results from the requirements of the device mounted on the ophthalmic table.

- ➔ Observe the instructions for use of the installed device.

4.3 Kontraindication

Non known.

4.4 Possible side effects

Non known.

4.5 Intended users

The OCULUS Ophthalmic Table is intended for use

- in optometrist practices,
- in clinics
- and by opticians.

It must not be used in operating rooms.

4.6 Patient group

Children from 3 years up to unlimited. No restrictions on weight, health and condition.

5 Device Description

The OCULUS Ophthalmic Table enables the instruments standing on it to be moved up and down in order to position them in the best possible way in front of the patient.

Depending on the version, the ophthalmic table is suitable for accommodating one or two devices.

Fig. 5-1: Equipment overview of the ophthalmic table

No.	Designation	Description
1	"Up" and "Down" keys	to move the ophthalmic table up and down.
2	Lift column	Depending on the version, the lift column has a different load capacity (see Technical Data chap. 13, page 64).
3	Levelling feet (4 piece)	to set the ophthalmic table horizontally.
4	Network cable connector (LAN B)	to connect the devices at the ophthalmic table to a local network
5	Network cable connector (LAN A)	to connect the devices at the ophthalmic table to a local network
6	Power supply connection	to connect the power cord.
7	Fuse	to exchange the fuses see chap. 9.2, page 59
8	On/Off Switch	The mounted devices are supplied with power via the ophthalmic table. Therefore, use the on/off switch to turn the connected instruments on or off.

6 Setup and Connection

Note

Incorrect setup and connection of the ophthalmic table may result in incorrect measurements and/or damage to the devices installed on the ophthalmic table.

- The ophthalmic table is assembled, set up and connected by OCULUS Service or an authorized dealer.
- Please observe the separate assembly instruction.

6.1 Setup

Avoid the following when setting up

- any vibrations,
 - jolts,
 - contaminants or
 - high temperatures.
- Set up the ophthalmic table on a level surface. Use the levelling feet if needed.
 - For version with wheels: Secure all wheels with the parking brake.

6.2 Operating Conditions

Temperature	+10°C - +35°C
Air humidity	30% - 90%
Air pressure	800 hPa - 1060 hPa

6.3 Electrical Connection

Warning

Personal injury due to electric shock

- Connect this device only to a mains supply with a protective earth conductor.
- Only connect this device to a mains supply with protective earth.

Warning

Personal injury caused by electric shock due to incorrect power cord

- Use only a mains cable that meets the requirements of IEC 60227-1, Type 53, min. 0.75 mm² and of IEC 60320-1.

**Warning**

Electrical safety hazard

- Do not use a multiple power strip to connect the ophthalmic table.

**Notice**

Risk of equipment damage due to incorrect connection

If you do not connect the Corvis® ST properly, and the connection is live, the unit can be damaged within a short period of time.

- Do not use excessive force when connecting the electrical plug.
- Pay attention to the specifications on the nameplate.
- If the electrical plug is damaged, contact our service department or an authorized dealer to repair the damage.

- Plug the mains plug into a mains power outlet

6.4 Network Cable Connection

There are two network cable jacks (LAN A) and (LAN B) available for connection to a network.


Fig. 6-1: Network cable jack

There are two network cables routed from the network cable sockets to the top side of the ophthalmic table. You can connect one OCULUS device to each of these.

- While making the connection, pay attention to proper connection between the network, network cable socket and the OCULUS device.

7 Operation

Caution

Personal injury or property damage due to improper operation

→ Before initial operation: Let OCULUS or an authorized dealer train you in the operation of the Ophthalmic Table.

Pinching hazard when moving the ophthalmic table

→ During movement of the ophthalmic table: Make certain that neither you nor your patient has any fingers resting on the elevating column.

Notice

Material damage due to continuous mode of the lift column

→ Do not allow the lift column to move up and down in continuous mode (one minute operation and thereafter do not move the lift column up or down for 30 minutes).

7.1 Turn on and off

→ Use the switch (1) to turn on or off the power supply:

On = I

Off = 0

Fig. 7-1: On/Off switch

7.2 Up and down movement

Press the following keys to move the ophthalmic table:

up

down

The buttons are located on the side edge of the ophthalmic table.


Fig. 7-2: „Up“ and „down“ keys

8 Cleaning and Disinfection

Cleaning and disinfection of the Ophthalmic Table is described in this chapter.

Sterilization is not required.

- ➔ Heed the product descriptions and instruction manuals of products and equipment you use to care for, clean, and disinfect the unit and/or its accessories.

8.1 Time table for cleaning and disinfection

Interval	Activity
after each examination	disinfect forehead rest
after each examination	disinfect the chin rest if paper is not used
1x monthly or as required	Cleaning

8.2 Cleaning

- ➔ If necessary, clean the ophthalmic table with a soft, lint-free cloth and, if necessary, with neutral detergent.
- ➔ Wipe the surface dry with a paper wipe.
- ➔ Do not use aggressive, chlorine containing, abrasive or sharp cleaning agents.

Warning

Risk of electric shock if the Ophthalmic Table is not completely disconnected from the mains for the cleaning.

- ➔ Switch off the electric power supply before cleaning
- ➔ Pull the power plug before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.

8.3 Disinfection

- Disinfect the forehead rest after each examination.
- If you do not use paper for the chinrest, disinfect the chinrest after each examination.

- Disinfect the chin and forehead rest.

We recommend:

Disinfection wipes suitable for medical devices, for example:

Mikrozid sensitive wipes premium

Fa. Schülke & Mayr

Softpack 48 Stück / Art. Nr. 165711

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Tel: +4940521000 / Fax: +494052100318

E-Mail@schuelke.com

www.schuelke.com

Notice

Equipment damage caused by disinfectant solution

The disinfectant solution may damage the surface of the device if it is sprayed directly on it.

- Spray the disinfectant solution onto a cleaning cloth, do not spray it directly on the device.
-

9 Trouble shooting

Caution

Personal injury or damage to equipment due to improper fault correction

- If a fault occurs which you cannot correct by following the instructions below, then identify the equipment as being out of service and notify our Service Department or your authorized dealer.

9.1 Fault table

Error	Possible cause	Remedy
Ophthalmic table does not move.	No power is being supplied to the ophthalmic table.	<ul style="list-style-type: none"> → Check whether the power supply switch is turned on. → Exchange the fuse. → Call OCULUS Service if necessary

9.2 Exchange Fuse

Warning

Risk of electric shock if the Ophthalmic Table is not completely disconnected from the mains for these jobs.

- Turn off the Ophthalmic Table, *Abb. 7-1, Seite 55*.
- Pull out the mains plug before you change the fuse. Take hold of the power plug for this purpose; do not pull on the cable.

Caution

Risk of fire if the wrong fuse is used

- Use only the following fuse:
230 V: Micro-fuse 6.3 A, Order number 05100315

You need a screwdriver to exchange the fuse.

Proceed as follows:

- ➔ Turn off the Ophthalmic Table, *Abb. 7-1, Seite 55*.
- ➔ Unplug the power cord.
- ➔ Disconnect the power cable from the Ophthalmic Table.
- ➔ Press the catches (2) together.


Fig. 9-1: Fuse holder


Fig. 9-2: Pull fuse holder out

- ➔ Pull the fuse holder (3) out.


Fig. 9-3: Replace defective fuse

- ➔ Replace the blown fuse.
A blown fuse is indicated by burnt-out filaments (4).
- ➔ Insert the fuse holder (3). Make sure it fits properly. The pin (2) must fit in the recess.
- ➔ Connect the Ophthalmic Table to the mains.
- ➔ You can now switch on the Ophthalmic Table and start the examinations.

10 Transport, Storage and Disassembly

The Ophthalmic Table must be properly dismantled before being transported or stored.

10.1 Disassembly

- Switch off the power supply.
- Unplug the power cord. Pull only on the plugs and not on the cables.
- Loosen the connections to the equipment.

10.2 Transport

Notice

Equipment damage due to improper transport

- Avoid shocks, vibration and dust.

- Transport the ophthalmic table in accordance with the transport conditions:

Ambient temperature range from	-40°C to +70°C
Relative humidity range from	10% to 95%
Air pressure range from	500 hPa to 1060 hPa

- Always move the ophthalmic table to the lowest position before transporting it.
- Always push the lift table with the short side first.
- Secure the ophthalmic table after transportation by locking all wheels.

10.3 Storage

Notice


Equipment damage due to improper storage

- Avoid contamination, high temperatures and moisture.

- Store the ophthalmic table on accordance with the storage conditions:

Ambient temperature range from	-10°C to +55°C
Relative humidity range from	10% to 95%
Air pressure range from	700 hPa to 1060 hPa

11 Disposal

In accordance with Guideline 2012/19/EG of the European Parliament and of the Council dated 27 January 2003, and additionally also the Law of the Federal Republic of Germany on the Commercialization, Recall and Environmentally Compatible Disposal of Electrical and Electronic Equipment, old electrical and electronic equipment must be sent out for recycling and may not be disposed of in household trash.

- Dispose of the Ophthalmic Table in a proper manner.
- Observe the applicable disposal guidelines in your country.

12 Warranty Provisions and Service

12.1 Warranty Provisions

The Ophthalmic Table has been manufactured with care and by using premium materials and modern production technologies. If software is included with the deliverable items, it has been tested by us and meets the technical standards. Observe the following warranty requirements:

- It is important that you follow the Instruction Manual and Safety Instructions before and/or during use of the equipment.
- You are entitled to a warranty claim to the Ophthalmic Table in accordance with legal regulations.
- If modifications are made to the Ophthalmic Table by unauthorized persons, then all warranty claims are void. This is because improper modifications and repairs may result in considerable hazards to the user and to patients.
- Any transport damage must be reported at once to the shipping company; have the transport damage noted on the bill of lading so that complaint handling and compensation of damages can proceed in an orderly manner.
- In general, our Business and Shipping Terms applicable on the date of purchase shall apply.

12.2 Liability for Faulty Operation or Damage

OCULUS shall only accept responsibility for the safety, reliability and utility of the ophthalmic table if you observe the following requirements:

- Only use the equipment in conformance with this Instruction Manual.
- Only those safety features which are to be serviced or repaired by the user are found on or in the Ophthalmic Table. If assembly work, modifications, adjustments, repairs, changes or service is performed by unauthorized personnel, or if the Ophthalmic Table is improperly maintained or handled, then any liability by OCULUS is void.
- If the above-mentioned work is performed by authorized persons, then a certification of the work shall be requested from the service technician which shall state any changes to factory defaults or to operating ranges. This certification must contain the date of performance and statement of the performing firm, with signature.
- If requested, OCULUS will provide a list of spare parts and additional descriptive material to the service technician for this purpose.
- Make certain that only original OCULUS parts are used.

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CE marking

The device is a category I product.
CE conformity according to regulation (EU) 2017/745 on Medical Devices (MDR), Annex II and III

13.1 Ophthalmic tables without electrical unit (universal tables)

	10007495 (37420)	10007496 (37430)
Table top (Kit 1)	10007482	10007484
Table base (Kit 2)	10007483	10007485
Weight	29 kg (63.1 lbs)	33 kg (73.8 lbs)
Dimensions of the desk top (W x D)	630 x 500 mm (24.8 x 19.7 in)	770 x 500 mm (30.3 x 19.7 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a	84 kg (185.1 lbs)	
evenly distributed – with wheels	67 kg (147.7 lbs)	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)	84 kg (185.1 lbs)
Max. power consumption	725 VA	725 VA
Voltage	110 – 240 V AC	110 – 240 V AC
Power supply connection	50/60 Hz	50/60 Hz
Fuses	2 x T 6.3 A	2 x T 6.3 A

a. above the lift column

	10007500 (37440)
Table top (Kit 1)	10007501
Table base (Kit 2)	10007502
Weight	39 kg (86 lbs)
Dimensions of the desk top (W x D)	980 x 500 mm (38.6 x 19.7 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a	
evenly distributed – with wheels	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)
Max. power consumption	725 VA
Voltage	110 – 240 V AC
Power supply connection	50/60 Hz
Fuses	2 x T 6.3 A

a. above the lift column

	10007584 (37450)
Table top (Kit 1)	10007585
Table base (Kit 2)	10007586
Weight	36 kg (79.3 lbs)
Dimensions of the desk top (W x D)	1050 x 570 mm (41.3 x 22.4 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a	
evenly distributed – with wheels	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)
Max. power consumption	725 VA
Spannung	110 – 240 V AC
Frequenz	50/60 Hz
Fuses	2 x T 6.3 A

a. above the lift column

	10007507 (37475)	100007510 (37460)
Table top (Kit 1)	10007508	10007508
Table base (Kit 2)	10007509	10007509
Weight	--	42 kg (93.0 lbs)
Dimensions of the desk top (W x D)	770 x 500 mm (30.3 x 19.7 in)	1050 x 500 mm (41.3 x 19.7 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg (147.7 lbs)	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)	84 kg (185.1 lbs)
Max. power consumption	725 VA	725 VA
Voltage	110 – 240 V AC	110 – 240 V AC
Power supply connection	50/60 Hz	50/60 Hz
Fuses	2 x T 6.3 A	2 x T 6.3 A

a. above the lift column

	10007513 (37485)
Table top (Kit 1)	10007514
Table base (Kit 2)	10007515
Weight	47 kg (103.6 lbs)
Dimensions of the desk top (W x D)	1250 x 500 mm (49.2 x 19.7 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a	
evenly distributed – with wheels	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)
Max. power consumption	725 VA
Spannung	110 – 240 V AC
Frequenz	50/60 Hz
Fuses	2 x T 6.3 A

a. above the lift column

	10007591 (37660)	10007699 (37665)
Table top (Kit 1)	10007592	10007600
Table base (Kit 2)	10007593	10007601
Weight	47 kg (103,6 lbs)	47 kg (103.6 lbs)
Dimensions of the desk top (W x D)	1120 x 550 mm (44.1 x 21.7 in)	1120 x 550 mm (44.1 x 21.7 in)
Lifting	655 – 905 mm (25.8 – 35.6 in)	655 – 905 mm (25.8 – 35.6 in)
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg (147.7 lbs)	67 kg (147.7 lbs)
evenly distributed – without wheels	84 kg (185.1 lbs)	84 kg (185.1 lbs)
Max. power consumption	725 VA	725 VA
Voltage	110 – 240 V AC	110 – 240 V AC
Power supply connection	50/60 Hz	50/60 Hz
Fuses	2 x T 6.3 A	2 x T 6.3 A

a. above the lift column

13.2 Ophthalmic tables with electrical unit (system tables)

	10007581 (37470)	10014624 (37471)
Table top (Kit 1)	10007582	10014663
Table base (Kit 2)	10007583	10014748
Weight	46 kg	46 kg
Dimensions of the desk top (W x D)	910 x 550 mm	910 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007575 (37480)	10015341 (37481)
Table top (Kit 1)	10007576	10015340
Table base (Kit 2)	10007577	10015339
Weight	48 kg	48 kg
Dimensions of the desk top (W x D)	1050 x 500 mm	1050 x 500 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007609 (37490)	10014642 (37491)
Table top (Kit 1)	10007610	10014667
Table base (Kit 2)	10007612	10014760
Weight	54 kg	54 kg
Dimensions of the desk top (W x D)	1250 x 550 mm	1250 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007620 (37495)	100014649 (37496)
Table top (Kit 1)	10007621	10014670
Table base (Kit 2)	10007622	10014762
Weight	58 kg	58 kg
Dimensions of the desk top (W x D)	1410 x 550 mm	1410 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007578 (37600)	100014703 (37601)
Table top (Kit 1)	10007579	10014704
Table base (Kit 2)	10007580	10014750
Weight	47 kg	47 kg
Dimensions of the desk top (W x D)	910 x 550 mm	910 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007588 (37620)	100014627 (37621)
Table top (Kit 1)	10007589	10014661
Table base (Kit 2)	10007590	10014754
Weight	49 kg	49 kg
Dimensions of the desk top (W x D)	1120 x 550 mm	1120 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007594 (37625)	100014636 (37626)
Table top (Kit 1)	10007595	10014666
Table base (Kit 2)	10007596	10014757
Weight	49 kg	49 kg
Dimensions of the desk top (W x D)	1120 x 550 mm	1120 x 550 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	230 V AC	115 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 3,15 A	2 x T 6,3 A

a. above the lift column

	10007604 (37550)
Table top (Kit 1)	10007605
Table base (Kit 2)	10007606
Weight	n.a.
Dimensions of the desk top (W x D)	1050 x 500 mm
Lifting	655 – 905 mm
Maximum load capacity ^a	
evenly distributed – with wheels	67 kg
evenly distributed – without wheels	84 kg
Max. power consumption	725 VA
Spannung	230 V AC
Frequenz	50/60 Hz
Fuses	2 x T 3,15 A

a. above the lift column

	10014654 (37670)	10014652 (37675)
Table top (Kit 1)	10014658	10014653
Table base (Kit 2)	10014764	10007485
Weight	n.a.	n.a.
Dimensions of the desk top (W x D)	1030 x 450 mm	770 x 450 mm
Lifting	655 – 905 mm	655 – 905 mm
Maximum load capacity ^a		
evenly distributed – with wheels	67 kg	67 kg
evenly distributed – without wheels	84 kg	84 kg
Max. power consumption	725 VA	725 VA
Spannung	100-240 V AC	100-240 V AC
Frequenz	50/60 Hz	50/60 Hz
Fuses	2 x T 6,3 A	2 x T 6,3 A
Classification per MPG	1	1

a. above the lift column

14 Appendices

14.1 Electromagnetic Compatibility (EMC)

Medical electric equipment is subject to special precautionary measures regarding EMC and must be installed and operated according to the EMC-instructions contained in the accompanying paperwork.

No particular measures are required for OCULUS equipment and systems.

Portable and mobile RF communications appliances can affect medical, electric equipment. The user of the Ophthalmic Table can help prevent electromagnetic interference by maintaining a minimum distance of 30 cm between portable and mobile RF communications equipment (transmitters) and the device.

Caution

Personal injury or property damage due to electromagnetic interferences

This device has been tested for immunity to interference against emitted RF fields only for certain selected frequencies and the use of wireless transmission devices in the immediate vicinity of this device, which work at different frequencies, may lead to faulty operation of the device.

Caution

The use of accessories, transducers and cables not specified by OCULUS (for example as replacement parts) may result in increased emissions or decreased immunity of the Ophthalmic Table.

→ Use only the original accessories, transducers and cables specified by OCULUS.

The use of accessories, transducers and cables specified by OCULUS with devices other than the Ophthalmic Table may result in increased emissions or decreased immunity of the other device.

→ Do not use the accessories, transducers and cables specified by OCULUS with devices other than the Ophthalmic Table.

To be in compliance with the requirements of the IEC 60601-1-2. 6.1 and 6.2 only accessories, converters and cables that have been tested and supplied by OCULUS may be used in conjunction with the device.

14.2 Guidance and Manufacturers declaration - Electromagnetic Emissions and Immunity

Guidance and manufacturer's declaration electromagnetic emissions IEC 60601-1-2, based to table 1

The OCULUS Ophthalmic Table is intended for operation in the electromagnetic environment specified below. The user of the Ophthalmic Table should ensure that it is being used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Ophthalmic Table uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
HF-emissions CISPR 11	Class B	
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	complies	

Electromagnetic immunity, IEC 60601-1-2, based on table 4

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic immunity, IEC 60601-1-2, based on table 5, 8			
Electrical Fast transient/bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines	± 2 kV ----- ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 6100-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree 0% U_T ; 1 period and 70% U_T ; 25/30 periods Single-phase: at 0 degree 0% U_T ; 250/300 periods	0% U_T ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree 0% U_T ; 1 period and 70% U_T ; 25/30 periods Single-phase: at 0 degree 0% U_T ; 250/300 periods	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ophthalmic Table requires continued operation during power mains interruptions, it is recommended that the Ophthalmic Table be powered from an uninterruptible power supply or battery.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Electromagnetic immunity, IEC 60601-1-2, based on table 4, 5

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
Conducted RF IEC 61000-4-6	3 V _{eff} 150 KHz to 80 Mhz 6 V in ISM- and amateur radio frequency bands between 150 kHz and 80 MHz 80% AM to 1 kHz	V _{eff} = 3 V	Portable and mobile RF communications equipment should be used no closer to any part of Ophthalmic Table, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{(V_1)} \right] \sqrt{P}$ $d = \left[\frac{3,5}{(E_1)} \right] \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz		$d = \left[\frac{7}{(E_1)} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interface may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1:

At 80 Hz and 800 MHz, the higher frequency range applies.

Note 2:

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Ophthalmic Table is used exceeds the applicable RF compliance level above, the Ophthalmic Table should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Ophthalmic Table.
- b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Ophthalmic Table

The Ophthalmic Table is intended for use in an electromagnetic environment in which radiated RF disturbances are uncontrolled. The customer or the user of the Ophthalmic Table can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Ophthalmic Table as recommended below, according to the maximum output power of the communications equipment.

Separation distance according to frequency of transmitter m

Rated maximum output power of transmitter W	150 KHz to 80 Mhz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.80	3.80	7.3
100	12	12	23


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.


NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

14.3 Description of the Connection

With Y-Cable

Without Y-Cable

WWW.OCULUS.DE

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Münchholzhäuser Str. 29 • 35582 Wetzlar • GERMANY

Tel. +49 641 2005-0 • Fax +49 641 2005-255

E-mail: sales@oculus.de • www.oculus.de

10031438 / Rev00
Lot:

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OCULUS Optikgeräte GmbH