

# OCULUS Centerfield®



INSTRUCTION MANUAL  
Examination of the Central Visual Field

## Notes on this instruction manual

To ensure safe operation, it is essential that you use the device correctly. For this reason you should familiarize yourself thoroughly with the contents of this instruction manual before operating the device. In particular, pay attention to the safety instructions.

- This instruction manual describes how to manage patient data, the default settings of the Centerfield® 2 program and the measuring process.
- The Centerfield® 2 user guide contains information supplementing the description of the operating concept.

Due to ongoing development your device might present minor differences compared to the information contained in this manual.

If you have any queries or would like additional information about your device, do not hesitate to call or send us an e-mail or a fax. Our service team will gladly assist you.

OCULUS Optikgeräte GmbH



OCULUS is certified according to DIN EN ISO 13485 and thus has established a high standard of quality for development, manufacture, quality assurance and service for our entire product line.

## Table of contents

1	Scope of Delivery .....	1
1.1	Software Version .....	2
2	Symbols .....	3
3	Structure of the Documentation .....	4
4	Safety instructions .....	5
4.1	About this Manual .....	5
4.1.1	Pictograms Used in this Manual .....	5
4.2	Safety Instructions for Use .....	6
5	Intended Use .....	12
6	Device Description .....	13
6.1	Overview of device components .....	13
6.2	How the Centerfield® Perimeter works .....	14
7	Operation .....	15
7.1	Installation of the software .....	15
7.1.1	Requirements .....	15
7.1.2	Installing the Software .....	15
7.2	Setup .....	16
7.3	Electrical Connection .....	16
7.4	Setup tasks for Initial Start-Up .....	17
8	Daily operation .....	18
8.1	Switching on Centerfield® 2 .....	18
8.2	Switching off Centerfield® 2 .....	18
9	Patient Data Management .....	19
9.1	Starting Patient Data Management .....	19
9.1.1	Enter a new patient .....	20
9.1.2	Select an existing patient .....	20
9.2	Starting the Centerfield® 2 program .....	21
10	The Centerfield® 2 program .....	22
11	Measurement Procedure .....	23
11.1	Preparing the Examination .....	23
11.1.1	Selecting the Examination program .....	23
11.1.2	Determining the Required Correction .....	23
11.1.3	Inserting the Corrective Lens .....	25
11.1.4	Check the Examination Conditions .....	25
11.1.5	Selecting the Examination Program .....	25
11.1.6	Preparing the Patient .....	25
11.1.7	Positioning the patient .....	26
11.1.8	Preparing the Measurement .....	27
11.1.9	Measuring the Pupil .....	27
11.2	Starting the examination .....	28

11.3	Interrupting the Examination .....	30
11.4	Ending the Examination .....	30
11.4.1	Saving the Examination Data.....	30
11.5	Re-Examination .....	31
12	Working with the Patient Data Management system.....	33
12.1	Rename Patient Data .....	33
12.2	Forwarding of examination results .....	33
12.3	Importing Patient Data.....	35
12.4	Data Backup.....	36
12.4.1	Backup Data .....	37
12.4.2	Reconstruct Data .....	37
12.4.3	Automatic Backup .....	38
13	Cleaning, Disinfection and Maintenance .....	39
13.1	Cleaning .....	39
13.2	Disinfection .....	40
13.3	Maintenance.....	41
13.3.1	Changing the Background Illumination Bulb .....	42
13.4	Changing the Chin Rest.....	43
14	Troubleshooting.....	44
15	Transport and storage.....	46
15.1	Disassembly and packing .....	46
15.2	Transport and Storage Information.....	46
16	Disposal.....	47
17	Terms of Warranty and Servicing .....	48
17.1	Terms of warranty .....	48
17.2	Assumption of liability for functions and damage .....	48
17.3	Manufacturer's address and service department:.....	49
18	Technical Data .....	50
19	Annex.....	53
19.1	Electromagnetic Compatibility.....	53
19.2	Guidance and Manufacturer's Declaration - Electromagnetic Emmissions and Immunity for the Centerfield® 2 .....	54
19.3	Description of the Connection.....	58
19.4	Data Sheet GSM60B15-P1J (05150725).....	59

# 1 Scope of Delivery

Product and accessories	Order number
Centerfield® 2-	56980
Power supply 15 V DC, 4.0 A	05150725
<ul style="list-style-type: none"> <li>■ Power cable 230 V</li> </ul> or	0520032
<ul style="list-style-type: none"> <li>■ Power cable 115 V</li> </ul>	05200210
Dust protection cover	56950 00 002
USB cable	05200570
USB cable	05200560
USB FS MED Isolator	015692000010
Replacement bulb (halogen) 12 V / 5 W for background illumination	05160060
Eye shield	44560
Hand-held button	56517
Lens holder for inserting trial lenses	085695012000
Instruction Manual	G/56980/0000/E 1219 Rev02
User Guide	BH/56980 / .../en
Software Installationn	SI/50000/.../en
Electrical safety test certificate	
Telescope box for accessories	9998027

Optional accessories	Order number
If no trial lens case is supplied	55900 XX XXX
<ul style="list-style-type: none"> <li>■ Trial lens case for perimetry</li> </ul>	10 100: +1.0 dpt 20 100: -1.0 dpt 10 300: +3.0 dpt 20 300: -3.0 dpt
Computer	70519
Laptop	59805
Floating license key with manual	77900 SI/77900/.../de

Optional accessories	Order number
■ Lifting table	37374 37377
Compact Laptop Stand	37499
■ Chin rest,electrically height adjustabler	56985
■ Chin pad (flexible plastic)	56985 01 009
■ Inkjet printer	56908
■ Transport box for Centerfield® 2	56984
■ Replacement bulbs (2 of)	56966
■ Frame for lens holder (plastic part as spare)	56950 12 003
■ CLIP Strategy	56915

- If you find transport damage upon delivery, immediately file a claim with the transport company.
- Have the damages noted on the bill of lading so that your claim for damages can be handled properly.



#### Note

We reserve the right to change the scope of delivery in line with ongoing technical developments.

## 1.1 Software Version


















This User Manual describes the following versions of the Centerfield® 2-Software together with patient data management:

- Centerfield® 2-Software: from Version 3.19r1477
- Patient Data Management: from Version 6.08



- The software version for patient data management is displayed on the "settings" page inside the patient data management.
- The software version of the Easyfield®-Program is displayed on the "settings" page inside the Centerfield® 2-Program.

## 2 Symbols

Symbols on device		Symbols packaging	
	Manufacturer		Protection class
	Keep dry		Conformité européenne
	This way up	IP XX	Type of protection
	Follow instruction for use		Article number
	Fragile		Disposal in household trash is prohibited
	Serial number		Transport Limit of temperature for transport
	Applied part Type B		Caution
			Storage Limit of temperature for storage
			Limit of humidity
	Example: UDI number, consisting UDI-DI (Device-Identification) UDI-PI (Product Identifier) machine-readable matrix code		Limit of air pressure





Fig. 2-1: Name plate (example)

## 3 Structure of the Documentation

A folder containing a set of documentation is supplied with your Centerfield® 2 perimeter:

- **Quick Guide:** In this document a checklist of the measurement procedure is supplied. The document is intended to help you familiarize with the equipment. The Quick Guide reminds you all important steps in order to obtain correct results.
- **Instruction Manual:** The design of the unit is described in detail in this document. The instruction manual also gives you general information about working with the Patient Data Management system and all safety-related instructions for use of the Centerfield® 2 perimeter.



### Caution

All safety-related instructions for use of the Centerfield® 2 perimeter are given in the Instruction Manual for the unit. It is therefore imperative that you read and understand the whole Instruction Manual before you use the Centerfield® 2 perimeter.

- 
- **User Guide:** All features of the examination and analysis software are described in the instruction manual, along with detailed information about the Patient Data Management system.
  - **Software Installation:** The introduction to the Software Installation describes how to install the Centerfield® 2 perimeter software and the associated drivers.

If you use a **Floating License Key:** information on the use of the Centerfield® 2 within networks.

## 4 Safety instructions

This chapter contains a summary of the most important safety-related information.

### 4.1 About this Manual

- ➔ Carefully read through the Instruction Manual.
- ➔ Keep the Instruction Manual in good condition near the device.
- ➔ Observe the legal regulations with regard to accident prevention.

#### 4.1.1 Pictograms Used in this Manual



##### Warning

Indicates a potentially hazardous situation that may result in irreversible injury.

---



##### Caution

Identifies a potentially dangerous situation which may cause minor injury or damage to property.

---



##### Note

Denotes situations which could result in incorrect findings, denotes user instructions and useful or other important information.

---



Identifies important information about the product and its usage, which require special attention.

---

- > This symbol denotes menu paths and screen shots. Example for starting a new examination:  
Centerfield® 2 > Examination > New patient  
which means:
  - ➔ Open the Centerfield® 2 program.
  - ➔ Select the "Examination" menu from the menu bar.
  - ➔ Select the menu item "New patient".

## 4.2 Safety Instructions for Use

---



### Caution

Personal injury or property damage due to improper operation.

- ➔ Observe the following safety instructions.

Personal injury or property damage due to equipment modifications related to safety.

- ➔ No modifications may be made to this device without the permission of the manufacturer.
- 

### Instructions to operating personnel

- ➔ The Centerfield® 2 must only be operated by people with adequate knowledge and practical experience to guarantee proper handling.

### Transport and Storage Instructions

Refer to the notes in [sec. , page 44](#).

### Instructions regarding the set-up and connection

- ➔ Do not use or store the Centerfield® 2 in damp rooms and do not store the unit in such areas.
- ➔ Keep the Centerfield® 2 away from water that may drip, surge or splash and make sure that no liquids can enter the Centerfield® 2 device. Do not place any containers with liquid either close to or on the Centerfield® 2 device.
- ➔ Germany: Only operate the Centerfield® 2 in rooms used for medical purposes if VDE 0100-710 installation procedures have been observed.
- ➔ Do not operate the devices included in the delivery in areas where explosions may occur, where there are inflammable anesthetics or volatile substances such as alcohol or petrol nearby.

- Use a power cord which meets the requirements of IEC 60227-1, type 53, minimum 0,75 m<sup>2</sup> and IEC 60320-1.
- Set up the Centerfield® 2 so that the power plug is easy to access. That way, you can easily disconnect it from the power supply for any repairs or maintenance work.
- Do not use excessive force when connecting the electrical plug. If you are unable to make a plug connection, check whether the plug fits the socket.  
If the electrical plug is damaged, contact our service department or an authorized dealer to repair the damage.
- Establish a USB connection only with the OCULUS USB FS MED-Isolator (Nr. 01 56920 00 010).
- Note that an output voltage of maximum 5.5 V DC is supplied by a device connected via USB.

#### Patient environment information

Patient environment is the area where patients can come into contact with any part of a medical electrical equipment (ME equipment) or with another person being in contact with the ME equipment.

In the patient environment, use devices that conform to IEC 60601-1. If a multiple power socket is to be used, or if a device is to be used that does not meet the IEC 60601-1 standard, use an isolating transformer




Fig. 4-1: Patient environment

#### Information about the operation of an ME system

The Centerfield® 2 and a connected computer form a medical electrical system (ME system) according to IEC 60601-1. If you connect additional devices, such as, for example a printer, those devices become part of the ME system.

- ➔ Make sure that all devices of the ME system meet the requirements of IEC 60601-1 or IEC 60950-1/IEC 62368.
- ➔ Note that an output voltage of maximum 5.5 V DC is supplied by a device connected via USB.

### Instructions for operation

- Before initial operation: Let OCVLUS or an authorized dealer provide necessary training for operation of the Centerfield® 2.
- Never operate a damaged Centerfield® 2.
- Only operate the Centerfield® 2 using original accessory parts supplied by us, and when the device is in technically correct working order.
- Do not touch the patient and the Centerfield® 2 simultaneously.
- Do not use the device if you have not understood the Instruction Manual.

### Notes on maintenance

- When cleaning, use a damp cloth and make sure that no liquid enters the Centerfield® 2.
- To ensure that it functions correctly and safely we recommend the following: Have the Centerfield® 2 checked every two years by our service department or an authorized dealer. If an error occurs which you are unable to correct, label the Centerfield® 2 as "out of order" and contact our service department or an authorized dealer.

### Instructions on disassembly and disposal

- When disconnecting electrical connections, pull on the respective plug instead of the cable itself.
- Dispose of the unit in conformance with legal requirements.

### Instructions on Electrical Safety



#### Caution

Risk of personal injury or damage to property due to an incorrect level of safety

Connecting the Centerfield® 2 with its non-medical electrical equipment (for example data processing equipment) to a medical electrical system must not result in a patient safety level below that prescribed by IEC 60601-1. If making this connection leads to the leakage current threshold being exceeded, protective measures including a circuit breaker must be in place.

- Ensure that connections with non-medical devices are made correctly.
- Only use the power adapter listed in the packing list.
- Use only a computer that meets the specifications given in this instruction manual, [sec. 18, page 50](#).
- Note that an output voltage of maximum 5.5 V DC is supplied by a device connected via USB.



### Caution

Use of a multiple socket extension cord

Risk of personal injury or material damage caused by unsafe multiple socket extension cord

If you use a multiple socket extension cord to connect the Centerfield® 2 to the power supply, you must heed the following information:

- Use an extension cord that complies with the requirements of IEC 60601-1, section 16.
- Do not place the multiple socket extension cord on the floor.
- Do not use more than one multiple socket extension cord.
- Plug only the Centerfield® 2 and the computer that is being used with the unit (if applicable) into the multiple socket extension cord.

If you are using a multiple socket extension chord it has to be supplied with a isolation transformer.

If you are using a new computer for the Centerfield® 2, you must have the electrical safety checked. Call OCULUS Service for this purpose.

### Electromagnetic Compatibility (EMC) / Cables

Risk of personal injury or damage to property due to electromagnetic interference

Portable and mobile RF communications equipment can affect medical electrical equipment, [sec. 19, page 53](#).

- Make sure that portable and mobile RF communications equipment do not cause interference.
- Recommendation: Maintain a minimum distance of 4 m. If the distance is shorter, you must ensure that the Centerfield® 2 functions correctly.

### Cybersecurity



Do not use the Centerfield® 2 with wireless technology, for example with wireless USB

To ensure cyber security in order to the usage of the device, the following security measures should be considered. Contact your computer administrator:

#### Precautions for access control of the computer

- Secure the computer with a password (for example at Windows start up).
- Choose a complex password: A good password should be at least eight characters long and are not in the dictionary. In addition to letters, it should also include numbers and special characters.
- Do not choose a name or device name for a password (for example "Centerfield").
- Change the password regularly.
- Do not note the password in an accessible location.
- Use different passwords for different users.
- Enable the screen saver and use the option for the necessity of re-entering the password when exit the screen saver.
- Choose an adequate time setting for starting the screen saver if software session is inactive (e.g. 10 minutes). Adequate time setting should consider duration of examination, number of patients, time between examinations, use of other devices in the examination room, several user, etc.
- Lock the computer if you are leaving the workstation (shortcut: 'windows logo key' + 'L').

#### Precautions if the computer is connected to a LAN or internet network

- Prefer wired connections of the computer to the network.
- If you are using Wi-Fi connections nevertheless, please ensure the usage of adequate security methods (for example WPA2/AES – Wi-Fi Protected Access / Advanced Encryption Standard – with a strong network key).
- The usage of a firewall (software or hardware) is recommended.

Recommendation: Use anti-malware tools with up to date malware definitions.



#### Note

Also observe the regulations, notes and recommendations of the *Bundesamt für Sicherheit in der Informationstechnik* for the protection of critical infrastructures.

---

## 5 Intended Use

The Centerfield® Perimeter is intended for the use described in this operating manual. It is designed for testing the visual field of the human eye. It permits kinetic and static, as well as automatic visual field examinations for this purpose.

You must use the power supply belonging to the device in order to operate it. There is no other method of connecting the device.



For this reason the Centerfield® Perimeter may only be operated by personnel instructed to do so, who, with appropriate training, knowledge and practical experience, are able to ensure proper handling of the device.

### Contraindication

not known

## 6 Device Description

### 6.1 Overview of device components



- |   |   |   |  |
|---|---|---|--|
| 1 | Centerfield® 2  | 6 | Port for connection with hand-held button      |
| 2 | Headrest  | 7 | Port for external power adapter                |
| 3 | Perimeter cone with mounts for correction lens holder | 8 | Label power adapter                            |
| 4 | Electrically height adjustable chin rest (optional)   | 9 | USB port for connection with netbook/PC/laptop |
| 5 | Power On/Off switch                                   |   |  |

Fig. 6-1: Overview Centerfield® 2

## 6.2 How the Centerfield® Perimeter works

The Centerfield® Perimeter is a bowl-type projection perimeter that is used to examine the visual field. The back surface projection principle is used in the unit. Various projectors are mounted on a movable arm. Every position of the perimeter bowl can be accessed. The perimeter bowl has a radius of 30 cm and is homogeneously illuminated to the Goldmann standard (referenced to a background luminance 10 cd/m<sup>2</sup>). The stimuli are presented accurately with exact reproducibility of the test locations - an absolute must for reliable visual field findings.

The unit is computer-controlled; connection to the computer takes place via the USB port.

The Centerfield® Perimeter meets the requirements of ISO 12866.

### Applied parts



1 Headrest  
2 Chinrest  
3 Hand-held button  
Fig. 6-2: Applied parts

## 7 Operation

Before you use the Centerfield® 2 for the first time you have to

- install the software, [sec. 7.1, page 15](#)
- set up the Centerfield® 2, [sec. 7.1, page 15](#)
- connect the Centerfield® 2, [sec. 7.1, page 15](#)
- install the firmware and drivers, [sec. 7.4, page 17](#)



### Caution

Improper set-up may result in incorrect measurements or equipment damage

- ➔ Before you use the device for the first time, the Centerfield® 2 must be installed and connected by our service personnel or by an OCULUS-authorized specialist.

## 7.1 Installation of the software



- ➔ Install the software first before you connect the unit to your PC.

### 7.1.1 Requirements

- No other software programs (screen saver, applications, etc.) must run simultaneously with the examination program in the foreground on the computer that controls the Centerfield® Perimeter.
- Energy-saving modes (BIOS or Windows) have to be deactivated ([sec. 7.1.2, page 15](#)).

### 7.1.2 Installing the Software

If you are working with a PC or laptop, you must install the Centerfield® 2 software. The Centerfield® 2 software consists of the following programs which are installed together.

- Patient Data Management
- Centerfield® 2 program
- TNT program
- ➔ Proceed as described in the [Software Installation](#)
- ➔ After installing the software, restart the PC or laptop.

## 7.2 Setup

You will find the ambient temperatures in [sec. 18, page 50](#).

- Remove the Centerfield® 2 from the packaging and keep the packaging.  
You can then ship or transport the unit in the proper manner for any servicing or repairs that may arise. You can thus avoid unnecessary damage and costs.
- Place the Centerfield® 2 on a level surface.
- Place the device so that direct light cannot affect the measurement.  
Make sure the examination is free from light reflections. To achieve this, darken the examination room.

## 7.3 Electrical Connection

You must connect the Centerfield® 2 to the power supply, and depending on the configuration, to the netbook, the laptop or PC. Connection and set-up will be demonstrated using the example of a netbook connection.



### Warning

Electrical shock due to incorrect power cord

- Use a power cord which meets the requirements of IEC 60227-1, type 53, minimum 0,75 m<sup>2</sup> and IEC 60320-1.
- 



### Caution

Electrical safety hazard

- Do not use the Centerfield® 2 adjacent to or stacked with other equipment.
  - If you have to use the Centerfield® 2 adjacent to or stacked with other equipment, verify the correct operation of the Centerfield® 2.
  - Only use the power adapter listed in the list, [sec. 18, page 50](#).
  - If you use a power strip to connect the Centerfield® 2: Use a power strip that complies with the requirements of IEC 60601-1.
  - Do not place the multiple socket extension cord on the floor.
  - Do not use more than one multiple socket extension cord.
  - Plug only the Centerfield® 2 and the computer that is being used with the unit (if applicable) into the multiple socket extension cord.
  - Use a socket with a protective earth connection which is fully operating.
-



**Note**


Risk of equipment damage due to incorrect connection

If you do not connect the Centerfield® 2 properly, and the connection is live, the unit can be damaged within a short period of time.

- ➔ Do not use excessive force when connecting the electrical plug.
- ➔ Pay attention to the specifications on the nameplate.

If the electrical plug is damaged, contact our service department or an authorized dealer to repair the damage.

- ➔ Plug in the hand-held button (1) and secure the plug into place by turning the knurled sleeve.
- ➔ Connect the Centerfield® 2 via port (2) with the PC/laptop. Connect the USB cable with the USB FS MED isolator. Connect the isolator with the computer



1 Port for hand-held button

2 Port for connection with netbook/PC/laptop

3 Port for external power supply adapter

Fig. 7-1: Connecting Centerfield® 2

- ➔ Connect the power supply adapter to the unit with the supplied power cable.

Make sure that the mains voltage is the same as the voltage specified on the rating plate of the power supply.

## 7.4 Setup tasks for Initial Start-Up

When you connect the Centerfield® 2 with a PC for the first time, you will have to perform several set up steps:

- ➔ To avoid communications problems, deactivate the power saver mode of the USB units of the operating system.
- ➔ Proceed as described in the [Software Installation](#).

## 8 Daily operation

### 8.1 Switching on Centerfield® 2

- ➔ Switch on the netbook, the PC or the laptop.
- ➔ Wait until the operating system has started completely and the Patient Data Management screen appears.
- ➔ Switch on the Centerfield® 2 at the On/Off switch.

### 8.2 Switching off Centerfield® 2


- ➔ Close the Centerfield® 2 program and the Patient Data Management.
- ➔ Shut down the Windows operating system.
- ➔ Switch off the Centerfield® 2 at the On/Off switch.

## 9 Patient Data Management


Use the Patient Data Management to input and manipulate patient data. Additional functions of the Patient Data Management system are found in [sec. 9, page 19](#) and in the [User Guide](#).

### 9.1 Starting Patient Data Management

After you have switched on the PC, it first loads the operating system.

➔ If necessary, click on the Centerfield® 2 icon: .

The user interface for the Patient Data Management will be displayed



1 "Functions" group box

2 Examination list

3 [Delete exam.] button

4 [Delete Patient] button

5 [Export] button

6 [Import] button

7 Patient list

8 Patient data

Fig. 9-1: User interface for Patient Data Management

In order to start the Centerfield® 2 program, you must first enter a new patient (9) or select a patient already present in the patient list (2).

Additional functions of the Patient Data Management system are found in the [sec. 9, page 19](#).

### 9.1.1 Enter a new patient

- ➔ Press the [New] button to enter a new patient into the Patient Data Management system.
- ➔ Enter the patient's last name, first name and date of birth completely in the patient data (8).





Fig. 9-2: Entering patients

Optionally, you can enter an ID number or additional comments for the patient.

- ➔ To save the data you entered, click [Save]. The patient you have just entered now appears in the patient list.
- ➔ Select the newly entered patient from the patient list and start the Centerfield® 2 program.

### 9.1.2 Select an existing patient

The patient list on the left-hand side of the screen displays all previously examined patients.



Lastname	Firstname	Date of birth
0_Patient	Anonymous	11.07.1976
Demo	Patient	06.10.1976
Test	Tessa	23.03.1962
Test	Theo	12.09.1965
Vv	Bb	07.08.1968

Fig. 9-3: Patient list

- ➔ Choose [Search] to quickly find the patient you are looking for in the list.
- ➔ Enter the patient's name or the first letter of the name in the "Last name" field.  
Optionally you can search for the patient using an ID number, assuming that one was assigned when the patient was first recorded.
- ➔ In the list, that appears, click the entry you were searching for to transfer the patient's name to the patient window. This also brings up a list of all previous examinations for that patient in the examination window (bottom right side).

Extended patient search: [Extended] checkbox

➔ Click on the checkbox [Extended].

The screen displays additional search parameters which reference previous examinations, for example. Proceed as for the input of a patient name.




Fig. 9-4: Extended search

## 9.2 Starting the Centerfield® 2 program

➔ Changing from patient data management > Centerfield program: After selecting a patient: Press the [Centerfield] button to start the Centerfield® 2 program.

or

➔ Double-click the selected patient name or an examination of the selected patient in order to start the Centerfield® 2 program.

## 10 The Centerfield® 2 program

You can get to the menu list from any screen of the Easyfield® program.




Fig. 10-1: Menüleiste Centerfield® 2-Programm



You will find the meaning and function of the symbols in the [User Guide](#).

### Loading previous examinations

- ➔ Select the menu item [Examination] and click [Load].  
The dialog box "Load Examination" appears.
- ➔ Make a selection by clicking the required examination.
- ➔ Confirm your selection by clicking [OK], or by double clicking.  
The Centerfield® 2 program will load the examination you have selected.

## 11 Measurement Procedure

---



### Caution

Risk of incorrect measurement due to incorrect use

- ➔ Before first use: Let OCULUS or an authorized dealer train you in the operation of the Centerfield® 2.
- 

## 11.1 Preparing the Examination

### 11.1.1 Selecting the Examination program

- ➔ Select the desired examination program from the "Program" tab panel.
- 



A description of how to write your own examination programs can be found in the [User Guide](#) of the Centerfield® 2 perimeter.

---

### 11.1.2 Determining the Required Correction

Correct measurement of the differential light sensitivity is only possible if the individual test points are focused sharply on the retina. The patient may need suitable corrective lenses for this purpose. It is possible to wear contact lenses (not colored ones); in certain conditions the patient's own glasses may be used during the examination.

To determine the required correction the exact refraction of the eye that is to be examined must be known. This can be taken either from a current refraction measurement, or from the patient's present ophthalmic lens strengths (distance vision correction).

As a patient's accommodation capacity greatly decreases with age, an age-related addition to the distance Rx is needed for patients aged approx. 40 years and older. The following are guidelines for this:

- Aged 40 – 50 years: approx. +1.00 D addition
- Aged 50 – 60 years: approx. +2.00 D addition
- Aged over 60 years: approx. +3.00 D addition

- Click in the "Correction" field. The following screen appears:




Fig. 11-1: Input of the known refraction values

- If you know the patient's refraction values: Enter the values into the fields in the „trial Lens Rx“ group box.
- Confirm by clicking [OK].
- If you don't know the patient's refraction values: Activate the „Calculate from Distance Rx“ checkbox.
- Enter the values into the fields in the „trial Lens Rx“ group box. The following screen appears:




Fig. 11-2: Input of the correction when the checkbox [Calculate from Distance Rx] is activated

- Enter the patient's previously determined refraction values into the fields in the "Current Distance Rx" group box. The values for the corrective lenses that are to be used are displayed in the fields of the "Trial Lens Rx" group box.
- Confirm by clicking [OK].

### 11.1.3 Inserting the Corrective Lens

- Place the required trial lens with the previously determined corrective power into the lens holder included in the delivery.
- Place the required trial lens with the previously determined corrective power into the holder.

### 11.1.4 Check the Examination Conditions

- Make sure that no interfering light is getting into the perimeter's viewer.
- To achieve the best possible results, darken the room somewhat.
- Make sure that the examination takes place in a quiet atmosphere and that the patient is not distracted.

### 11.1.5 Selecting the Examination Program


- Select the desired examination program on the "Program" tab panel. A description of how to write your own examination programs can be found in the User Guide for the Centerfield® 2 Perimeter.

### 11.1.6 Preparing the Patient

- Check that the forehead and chin rest and the hand held have been cleaned and disinfected before each examination.
- Make sure that the examination takes place in a quiet atmosphere and that the patient is not distracted.
- Explain the examination procedure to the patient.
- Give the patient the hand-held button for the unit and ask him to hold it in one hand.
- Ask the patient to comfortably take a seat in front of the unit. The patient should sit as upright as possible.
- Do not touch the patient and the Centerfield® 2 simultaneously.
- Cover the eye that is not being tested with the eyepatch.

### 11.1.7 Positioning the patient

- Do not touch the patient and the Centerfield® 2 simultaneously.



- |                                     |                                    |
|-------------------------------------|------------------------------------|
| 1 Chin rest                         | 3 Lens holder with corrective lens |
| 2 Height mark for the patient's eye | 4 Chin rest                        |

Fig. 11-3: Positioning aids

- Ask the patient to place his chin on the chin rest (4).
- If your perimeter is equipped with a height-adjustable chin rest, adjust the height with the cursor keys on the PC keyboard. Line up the patient's eye with the height mark (2).
- Make sure that the distance between the eye and the corrective lens, or the eye and the perimeter is no greater than 1 cm.




Fig. 11-4: Distance between eye and corrective lens

- ➔ Ask the patient to lean rest his forehead against the head rest so that he sees the fixation marks (four red dots) in the center of the perimeter bowl with the eye that is to be examined.  
Only pull out the forehead rest in exceptional cases.  
You can now see the center of the patient's eye in the camera image at the bottom left of the screen.
- ➔ Tell the patient to look towards the center of the fixation marks.



#### Note

Visual field losses in the top area could be caused by improper positioning of the patient. If the distance from the eye to the perimeter is too big (because the forehead rest has been pulled out, or the patient is not positioned properly), the patient may not have a full view into the unit.

### 11.1.8 Preparing the Measurement

- ➔ Select the eye that is to be examined in the "Eye" field of the Centerfield® 2 software.
- ➔ Click with the right mouse button in the camera image at the bottom left of the screen. The center of the pupil is automatically centered on the screen.
- ➔ Correct the position of the center of the pupil, if necessary, by pressing the appropriate arrow keys.



- ➔ If necessary, adjust the camera image settings in the Centerfield® 2 program settings so that the movements of the eye are displayed in the convenient way (mirrored or not mirrored).

### 11.1.9 Measuring the Pupil

To conclude the examination preparations, the pupil diameter must now be measured. To do this:

- ➔ Move the mouse pointer to the left edge of the pupil on the monitoring image.
- ➔ Press and hold down the left mouse button. The left edge of the pupil is marked with a green line.
- ➔ Move the mouse pointer to the right edge of the pupil and stop pressing the mouse button there.

The right edge of the pupil is also marked with a green line and the measured pupil diameter is displayed in the "Pupil" field.

## 11.2 Starting the examination

- ➔ Now instruct the patient to press the hand-held button every time a light spot is seen flashing.
- ➔ Explain to the patient the examination can be interrupted at any time by pressing and holding down the hand-held button. The examination is automatically resumed when the patient lets go of the hand-held button.
- ➔ Click on the button [Start Exam.]  
The following dialog appears so that you can check the data that you have entered:




Fig. 11-5: Display of the general data



Depending on the selected examination program, you can now select from a list box, whether the central threshold value or the peripheral threshold value is to be determined at the beginning of the examination.

- ➔ Check the data that have been entered. If you determine, for example, that you have selected the wrong eye, press the [Cancel] button.
- ➔ When all values have been entered correctly, ask the patient to once again look into the center of the four red dots.
- ➔ Press the [Start Exam.] button.

The central threshold value is determined and is displayed in the following dialog box.



Fig. 11-6: Display of the measured threshold value

- ➔ If the measured threshold deviates considerably from the normal threshold for the age group, press the [Repeat] button.



Depending on the selected examination program, you can also manually select the desired luminance class for the measurement from a list box.

- ➔ Tell the patient that the examination is about to start and press the [Start] button.

The examination program that you selected will now start to run.

Depending on the selected program and the examination area associated with it, if the patient needs a corrective lens, the following message may appear after the test has been completed in the visual field center:




Fig. 11-7: Instruction to remove the corrective lens

- ➔ First remove the corrective lens from the lens holder.
- ➔ Then press the [OK] button.
  - The threshold is determined again **without** the corrective lens and displayed.
- ➔ Press the [Start] button to continue the examination.

### 11.3 Interrupting the Examination

- If you want to interrupt the examination, press the right mouse button.

The following confirmation dialog box appears:




Fig. 11-8: Interrupt the examination

- If the examination can be resumed, press the [No] button.
- To cancel the examination completely, press the [Yes] button.

### 11.4 Ending the Examination

At the end of an examination, the following dialog box appears::




Fig. 11-9: Save the examination results

- Now decide how you want to proceed, based on the examination results.
- Tell the patient that the examination has been stopped and that the patient can relax.
- Disinfect the forehead and chin rest after each examination, [sec. 13.2, page 40](#).
- Disinfect the hand held button after each examination, [sec. 13.2, page 40](#).
- Disinfect the eyepatch after each use.

#### 11.4.1 Saving the Examination Data

If you are content with the findings, or you have already performed the desired re-examination, you can now save the examination results. To do this:


- Press the [Yes] button.

The examination data are saved and can be reloaded again later via the Patient Data Management program.

## 11.5 Re-Examination

If you want to check peculiar test points, you can conduct a re-examination. Proceed as follows:

- ➔ Press the [Re-Examination] button.
- Three additional buttons appear at the bottom right of the displayed examination results.



1 Additional buttons

Fig. 11-10: Additional buttons for a re-examination

You can determine the points for the re-examination in several ways.

- ➔ Press the [Autoselect] button. The peculiar points are then automatically selected.
- ➔ Manually select the points in the test point grid with the mouse.
- ➔ Press the [Add] button to manually add more points that are not yet present in the test point grid. Then click in the test point grid to define the additional test points.

- ➔ Press the [Add grid] button to add a predefined grid of test points. In this case, an additional dialog appears in which you can select the test point grid:




Fig. 11-11: Selection of the grid that is to be added for re-examination purposes

- ➔ Use the above options to define the points in the map that are to be tested during the re-examination.
- ➔ Finally, click on the [Finish] button. The following dialog will be displayed:




Fig. 11-12: Dialog "Continue examination"

- ➔ If you previously took the corrective lens out of the holder, re-insert it now.
- ➔ Tell the patient that you are now going to continue with the examination.
- ➔ Press the [Continue] button and, if applicable, confirm that you have re-inserted the corrective lens in the confirmation dialog box that then appears.

The examination is resumed.

After finishing the re-examination, a dialog box appears, asking whether you want to save the examination result ([Fig. 11-8, page 30](#)).

## 12 Working with the Patient Data Management system

This section describes how to work with the Patient Data Management system

- Rename it, [sec. 12.1, page 33](#)
- Export it, [sec. 12.2, page 33](#)
- Import it, [sec. 12.3, page 35](#)
- Save it, [sec. 12.4, page 36](#)



For more information on Patient Data Management, refer to the [User Guide](#).

### 12.1 Rename Patient Data


After creating of the patient data, you can edit it.

- ➔ Press the [Change] button.  
The input boxes for patient data are now enabled, and the cursor jumps to the "Last name" field.
- ➔ Change the entries in the individual boxes.
- ➔ Press the [Save] button.

### 12.2 Forwarding of examination results

For example, patient and examination data can be exported for forwarding to another clinic.

- ➔ Select the patient and also one of the examinations in the respective list as required.
- ➔ Press the [Export] button below the patients list. The following dialog will be displayed:



1 Saving destination selection

2 Button to select a folder

3 [Cancel] and [Export] buttons

Fig. 12-1: "Export patient data" dialog



The default options for import and export of data are configured in the "Settings" field, see also the [User Guide](#).

Depending on the settings, you may not have to perform all of the following steps (e.g. selection of the directory).

---

→ Select the destination type (1) you would like for the exported data.

---



Recommendation: Export the patient data using the "Single file (U12)" option.

---

- Press the [...] button (2) to select a folder.
- In the dialog that appears, select the folder or the file to which the patient data should be exported, e.g. B. TOPO.DAT for data and TOPO.BMP for the pictures.
- Confirm your selection with [OK] or [Save].
- To export the data, press the [Export] button (3).

## 12.3 Importing Patient Data

In case you keep patient data on a USB stick, you can import this data.




### Note

Loss of data due to computer viruses  
 Computer viruses can cause loss of data.

➔ Run a virus check before importing data from the USB stick.

➔ Press the [Import] button. The following dialog will be displayed:



1 Select the source of the data

2 [...] button to select a folder

3 Examinations list

4 [Import] button

5 Patient list

Fig. 12-2: "Import" dialog



The options for import and export of data are set as defaults in the "settings" field, see also the [User Guide](#).

➔ Depending on the settings you may not have to perform all the following steps (e.g. selection of the directory).

➔ Select the option (1) which contains the source data ("Folder" or "Single file (U12)").



Recommendation: Import the patient data using the "Individual file (U12)" option.

- Press the [...] button (2) to select a folder.
- In the dialog box, select the directory or the file where the patient data are located.
- Confirm your selection with [OK] or [Save].  
The patients that are located and the associated examinations are displayed in the lower part of the dialog.
- To import the data, press the [Import] button (4).  
The data will then be available in the Patient Data Management system.

## 12.4 Data Backup

You should make a backup copy of patient and examination data at regular intervals. In case of loss of data, you can reconstruct the data from a previously created backup with the help of this function. Since data backup takes several minutes, depending on the scope of the database and the data to be backed up, a backup should be carried out when the PC and the device will not be needed.



### Note

Risk of loss of data due to computer viruses

Computer viruses can cause loss of data.


- Run a virus check before making a backup to a USB flash drive.



The general rules for security backups apply to the creation of backup copies created with the help of the Patient Data Management system. Storage of backup files should always be done on a separate system (e.g. on a USB flash drive with adequate capacity).

### 12.4.1 Backup Data

- ➔ Press the [Backup] button in the "functions" field. The following dialog appears:



- |                         |                                     |
|-------------------------|-------------------------------------|
| 1 Backup data selection | 4 Backup directory and button [...] |
| 2 [Restore] button      | 5 [Save] button                     |

- 3 Display free storage space

Fig. 12-3: "Backup" dialog

- ➔ Select whether all of the data or only changed data should be backed up.



The Patient Data Management system internally tags all saved data records.

If you select the option "Changed and new data only", only the data records that were not saved during a previous backup will be backed up.

- ➔ Press the [...] button to the right of the "Backup directory" box (4).
- ➔ In the dialog that appears, select the folder to which the data should be backed up.
- ➔ Confirm your selection with [OK].
- ➔ To back up the data, press the [Save] button (5). The previously selected data will then be backed up to the corresponding folder.

### 12.4.2 Reconstruct Data

If a loss of data occurs, the data from a previous backup can be re-imported into the Patient Data Management system.

- ➔ Press the [...] button.
- ➔ In the dialog that appears, select the folder which contains the backup data.
- ➔ Confirm your selection with [OK].
- ➔ To import the data, press the [Restore] button (2). All data in the appropriate directory are copied to the Patient Data Management system.

### 12.4.3 Automatic Backup

In addition to the manually performed backup, it is also possible to automatically run a backup when exiting the Patient Data Management system. The settings required for this can be made in the "Settings" area, see [User Guide](#).

## 13 Cleaning, Disinfection and Maintenance

This section describes how to clean and disinfect the Centerfield® 2.

Sterilization is not required.

- Heed the product descriptions and instruction manuals of products and equipment you use to care for, clean, and disinfect the unit and/or its accessories.

### 13.1 Cleaning



#### Caution

Risk of electric shock if the Centerfield® 2 is not completely disconnected from the mains for these jobs.

- Turn the Centerfield® 2 off, [sec. 8.2, page 18](#).
- Pull the power plug before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.

- Do not clean the Centerfield® 2 with aggressive, chlorinecontaining, abrasive or sharp cleaning agents.

#### Required materials:

- Anti-static cleaner for plastic surfaces
- Cleaner for painted surfaces: mixture of equal parts mineral spirits and distilled water, with a few drops of household detergent if needed
- Soap solution
- Soft cloth or lens brush
- if necessary: optical brush
- Alcohol or lens cleaner

#### Cleaning intervalls

- Clean the chin rest, the head rest and the hand held button after each examination.
- Clean the device once a month or if necessary.

#### Cleaning the Painted Surfaces

- Make sure that the cleaning agent does not get into the unit.
- Wipe down the outer surfaces of the unit with a damp cloth.
- Wipe any residual particles using a mixture containing equal parts of spirits and distilled water, and add a dash of household washing up liquid. You can add a drop or two of household dish soap to this liquid.

### Cleaning the Chin Rest and Forehead Rest

- Clean the chin rest and the forehead rest with a soap solution (or with alcohol, if very dirty).
- Use a lint-free, damp cloth.

### Cleaning the Viewer Lens

- Clean the lens in the viewer with a soft cloth or an optical brush.
- If necessary, also use a little alcohol or lens cleaner.



To generally protect the unit, we recommend that after use, you cover it with the supplied cover.

---

## 13.2 Disinfection

---



### Caution

Risk of electric shock if the Centerfield® 2 is not completely disconnected from the mains for these jobs.

- Turn the Centerfield® 2 off, [sec. 8.2, page 18](#).
  - Pull the power plug before cleaning. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.
- 

Required materials (recommended) or an equivalent:

Pursept®-A Xpress disinfectant cloths,  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt | Germany  
Telefon: +49 40 52100-0  
Telefax: +49 40 52100-318  
E-Mail: [info@schuelke.com](mailto:info@schuelke.com)  
<https://www.schuelke.com/de-de/index.php>



#### Note

Equipment damage caused by disinfectant solution

The disinfectant solution may damage the surface of the device if it is sprayed directly on it.

- Spray the disinfectant solution onto a cleaning cloth, do not spray it directly on the device.

- Disinfect the forehead and chin rest after each examination; disinfect the housing as needed.
- Disinfect the hand held button after each examination.
- Disinfect the eyepatch after each use.

## 13.3 Maintenance

The Centerfield® 2 Perimeter is designed such that special maintenance is not necessary. The electronic brightness adjustment function always adjusts the unit to the specified set values. For safety reasons, we recommend that the illumination and electrical values be checked every two years.

- Please contact OCULUS Service for this.



#### Note

Damaged equipment will result in erroneous examinations

If you use a damaged unit, the examination result may be incorrect.

If an error occurs that you cannot correct:

- Identify a damaged Centerfield® 2 as being out of service.
- Report the damage to OCULUS Service department or to your authorized dealer.
- Only use an undamaged Centerfield® 2.

### 13.3.1 Changing the Background Illumination Bulb

Required materials:

- small screw driver
- halogen bulb of the type 5160060 (12 V, 5 W)

The background illumination bulb is located under a cover, below the unit's viewer (1).



1 Cover of the background illumination bulb

Fig. 13-1: Removing the cover of the background illumination bulb

- ➔ First close the program and turn off the unit in the proper manner.



#### Caution

Risk of electric shock if the Centerfield® 2 is not completely disconnected from the mains for these jobs.

- ➔ Turn the Centerfield® 2 off, [sec. 8.2, page 18](#).
- ➔ Pull the power plug before changing. When disconnecting electrical connections, pull on the respective plug and not on the cable itself.

The background illumination bulb may still be hot.

- ➔ Wait until the bulb has cooled down before you attempt to change it!

- ➔ Remove the small screw at the cover.
- ➔ Carefully tip the cover downwards out of the housing.
- ➔ Take hold of the background illumination bulb with a cloth and pull it out of the socket.
- ➔ Insert the new bulb into the socket.
- ➔ Make sure that you do not touch the glass of the bulb with your fingers.
- ➔ Put the cover back into place in the housing.
- ➔ Secure the cover with the screw.
- ➔ Connect the Centerfield® 2, [sec. 7.3, page 16](#)

## 13.4 Changing the Chin Rest

Required materials:

- ➔ hex ball driver SW 2
- ➔ new chin rest

There are two fastening screws underneath the chin rest (1).



1 Fastening screw

Fig. 13-2: Fastening of the chin rest

- ➔ Loosen the two fastening screws.
- ➔ Pull the chin rest up and remove it from the holder.
- ➔ Insert the new chin rest and fasten it into place with the two screws.

## 14 Troubleshooting



### Caution

Personal injury or equipment damage due to improper troubleshooting.

- ➔ If an error occurs which you are unable to correct by following the instructions below, label the device as "out of order" and contact our service department or an authorized dealer.

Damaged device due to improper use

- ➔ Do not plug or unplug a cable or plug while the PC or the Centerfield® 2 is running. Devices can be destroyed.

Fault	Possible Cause	Remedy
No function when the power switch is pressed or the pilot lamp on the power switch does not light up.	The Centerfield® 2 perimeter is not connected to the power supply.	Plug the power cable into the power outlet, or the inlet connector into the jack at the Centerfield® 2 perimeter.
	Power failure or power outlet is not active	Inform the in-house electrician.
	USB- or serial cable to the PC is not connected properly.	Check that the connector is plugged in properly.
No function when the power switch is pressed but the pilot lamp on the power switch is lit.	The unit has been switched off and back on again too quickly.	Wait approx. 5 seconds before turning the unit back on again.
The printer is not printing.	Connecting cable from Printer to PC is not plugged in properly. Ink cartridge empty.	Change the cartridge
Hand-held button does not respond when pressed.	The hand-held button is not properly plugged in and screwed tight in the jack at the unit.	Check the connection and plug in the cable again and screw it tight.
Camera image is too dark.	The camera brightness settings are incorrect.	Re-adjust the brightness (refer to the <a href="#">User Guide</a> ).

Fault	Possible Cause	Remedy
Background illumination not active.	Unit is in standby mode	Move the mouse, or press any key.
	The Centerfield <sup>®</sup> 2 program (examination program) has not been started.	Start the examination program <a href="#">(sec. 7.4, page 15)</a>
	Lamp is faulty.	Change the background illumination lamp <a href="#">(sec. 13.4, page 38)</a> .
After you have started the Centerfield <sup>®</sup> 2 program, the following dialog box appears: "No communication with the Centerfield!".	No power to the mains adapter.	Check whether the indicator light on the power supply is on. If not, connect the power supply to the mains.
	Connecting cable (USB cable) between the Centerfield <sup>®</sup> 2 and the PC is not plugged in properly.	Check whether the power supply cable is correctly attached to the Centerfield <sup>®</sup> 2 Perimeter.
	Software/Hardware problems.	Check whether the USB connector is properly inserted.
		Switch the Centerfield <sup>®</sup> 2 Perimeter off and restart the PC. Switch the Centerfield <sup>®</sup> Perimeter on as soon as Patient Data Management becomes active. When you start the Centerfield <sup>®</sup> 2 program, the message, "

## 15 Transport and storage

The Centerfield® 2 must be properly dismantled and packed before being transported or stored.

### 15.1 Disassembly and packing

- Select Examinations > New Patient / End.
- Exit the Patient Data Management system.
- Power down the PC/laptop.
- Shut off th device, [sec. 8.2, page 18](#).
- Disconnect the power plug from the power jack.
- Disconnect the connections to the hand-held button, to the netbook/PC/laptop.  
When disconnecting, pull on the respective plug and not on the cable itself.
- Pack up the Centerfield® 2 in the original packing.

### 15.2 Transport and Storage Information

#### Storage

Ambient temperature range	-10 – +55°C
Relative humidity, including condensation	10 – 95%
Air pressure range	700 – 1060 hPa

#### Transport

Ambient temperature range	-40 – +70°C
Relative humidity range from	10 – 95%
Air pressure range	500 – 1060 hPa

### After transport and/or storage

- ➔ Wait approx. 3-4 hours after transport before putting the Centerfield<sup>®</sup> 2 into initial operation. Extreme temperature changes from cold areas to warm rooms can cause condensation on the optical components.



#### Note

Risk of equipment damage due to incorrect shipment or from improper storage

- ➔ Avoid shocks, vibrations, and contamination.
- ➔ Avoid high temperatures and humidity.

- ➔ Transport the Centerfield<sup>®</sup> 2 in an appropriate manner.
- ➔ Store the Centerfield<sup>®</sup> 2 in accordance with the storage conditions.
- ➔ Do not store the unit near heating elements or moisture.
- ➔ Check the Centerfield<sup>®</sup> 2 for damage every time it has been transported.

## 16 Disposal



In accordance with Directive 2012/19/EC of the European Parliament and the Council and in accordance with German law governing the marketing, return and environmentally compatible disposal of used electrical and electronic devices, such appliances must be recycled and may not be discarded as household waste.

- ➔ Dispose of the Centerfield<sup>®</sup> 2 in an appropriate manner.

## 17 Terms of Warranty and Servicing

### 17.1 Terms of warranty

The Centerfield® 2 was carefully manufactured using quality materials and modern production methods. Any software included in the delivery was tested by us and complies with technical standards. Please note the following warranty provisions:

- Prior to and while operating the device it is important that you follow the user instructions, the instruction manual and safety instructions.
- The Centerfield® 2 carries a warranty to which you are entitled in accordance with the legal provisions.
- If any unauthorized persons interfere with the Centerfield® 2, all warranty entitlements shall be void. Any inappropriate modifications or repairs can cause grave danger to the user and patient.
- Any entitlement to a warranty shall also be void if unauthorized persons interfere with the PC hardware and software supplied.
- In the event of transport damage, we request that you notify the shipping company immediately and have the damage confirmed on the consignment note, to enable a proper claims settlement procedure.
- Overall, the general terms and conditions of business and delivery apply as per the date of purchase.

### 17.2 Assumption of liability for functions and damage

OCULUS will only be liable for the safety, reliability and utility of the Centerfield® 2 if you have followed the instructions below:

- ➔ Use the unit in conformance with the Instruction Manual and the included User Guide.
- There are no user-serviceable parts either on or inside the Centerfield® 2. OCULUS shall not assume any liability if assembly, extensions, adjustments, changes or repairs are carried out by unauthorized personnel, if the Centerfield® 2 is maintained improperly or if it is handled incorrectly.
- If the work described above is carried out by persons authorized to do so, they must be requested to supply documentation detailing the nature and scope of repairs, and if applicable to specify modifications to the rated data or area of work. This certificate must bear a date, a signature, specify who carried out the work, and contain company information.
- On request, and for this purpose, OCULUS will supply authorized persons with spare parts lists and additional descriptions.
- ➔ Make certain that only original OCULUS replacement parts are used for any repairs or maintenance.

### 17.3 Manufacturer's address and service department:

Our service department or authorized representatives will provide you with additional information. Address of the manufacturer and service department:

Germany:

OCULUS Optikgeräte GmbH  
Münchholzhäuser Straße 29  
35582 Wetzlar  
GERMANY  
Tel.: +49 641 2005-0  
Fax: +49 641 2005-295  
E-mail: [export@oculus.de](mailto:export@oculus.de)  
[www.oculus.de](http://www.oculus.de)



USA:

OCULUS, Inc.  
17721 59th Avenue NE  
Arlington  
WA 98223  
Tel. +1 425 670 9977  
Fax +1 425 670 0742  
E-mail: [sales@oculususa.com](mailto:sales@oculususa.com)  
[www.oculususa.com](http://www.oculususa.com)



## 18 Technical Data

### Measuring equipment

Weight (without chinrest)	11.7 kg (25.8 lbs)
Weight of the optional chin rest	1.1 kg /2.5 lbs)
Dimensions (W x D x H)	398 x 503 x 580 mm (15.7 x 19.8 x 22.8 in)
Interface	USB
Bowl radius	300 mm
Meridian	Adjustable from0°-360°
Eccentricity	up to 36° or 70° (with optional fixation shift)
Power supply	15 V DC, 4 A
Max. power consumption	30 W
Life expectancy	10 years

### Measuring parameters

Stimulus	
■ Stimulus size	Goldmann III
■ Stimulus colour	white, blue
■ Stimulus duration	200 ms/user-defined(0,2 s/0,5 s/0,8 s/adaptive)
■ Stimulus luminance range $L_s$ /steps	0 – 318 cd/m <sup>2</sup> (0 – 1 000 asb)/1 dB
■ Examination speed	adaptive / fast / normal / slow / user-defined
■ Speed (kinetic perimetry)	2°/s or user-defined
Background	
■ Luminance	10 cd/m <sup>2</sup> (31,4 asb)
■ Colour	white, yellow

### Power adapter

Power adapter	GSM60B15-P1J (05150725)
Mains connection	80 - 264 V AC 1,4 - 1 A
Frequency	47 - 63 Hz
Max. power consumption	68 W
DC output	15 V DC 4 A 60 W max.
Fuses	integrated over current shut-off

### Classification per IEC 60601 - 1 (VDE 0750)

Type of protection against electrical shock	Protection class II
Level of protection against electrical shock	Type B
Level of protection against damaging water entry	IP20

### Ambient operating requirements

Temperature	+10 - +35°C
Humidity	30 - 75%,
Air pressure	700 - 1060 hPa

### Storage requirements

Ambient temperature	-10 - +55°C
Relative humidity, including condensation	10 - 95%
Air pressure	700 - 1060 hPa

### Transport requirements

Ambient temperature	-40 - +70°C
Relative humidity, including condensation	10 - 95%
Air pressure	500 - 1060 hPa

### Computer

Use a computer which is in conformity with the IEC 60950/IEC 62368 standard.

Recommended computer specifications	Intel® Core™ i5, 500 GB HDD, 4 GB RAM, Windows® 7 Pro, Intel® HD Graphics 520
-------------------------------------	---

### CE in accordance with EC Directive 93/42/EEC for Medical Devices

The unit is a Class I product.



Conformity assessment: Directive 93/42 / EEC: annex VII

## 19 Annex

### 19.1 Electromagnetic Compatibility

Medical electrical equipment is subject to special precautionary requirements with respect to EMC, and must be installed and operated according to the EMC-Instructions contained in the accompanying paperwork.

No special measures need be observed in respect of OCULUS devices and systems.

Portable and mobile RF-communications devices can interfere with electrically operated medical devices.

**Produced in consideration of permissible deterioration during or caused by the EMC testing without affecting the essential performance criteria**

- A short interruption of the USB connection during the examination is permissible because it will not affect the diagnosis, treatment and observation.



#### Caution

The use of accessories, transducers and cables not specified by OCULUS (for example as replacement parts) may result in increased emissions or decreased immunity of the Centerfield® 2.

- ➔ Use only the original accessories, transducers and cables specified by OCULUS.

The use of accessories, transducers and cables specified by OCULUS with devices other than the Centerfield® 2 may result in increased emissions or decreased immunity of the other device.

- ➔ Do not use the accessories, transducers and cables specified by OCULUS with devices other than the Centerfield® 2.
-

To be in compliance with the requirements of the IEC 60601-1-2: 6.1 and 6.2 the following types of equipment, accessories, power adapters and cables must be used..

Order number	Description	
56980	OCULUS-Centerfield® 2	
05200320	Cable with plug, EU standard	2.5 m
05200210 (110 Volt)	Cable with plug, US standard	2.5 m
015692000010	USB FS Med Isolator	
05150725	Power supply adapter GSM60B15-P1J	see <a href="#">sec. 18</a> , <a href="#">page 50</a>
056517	Hand-held button (spi- ral cable)	< 3.0 m

## 19.2 Guidance and Manufacturer's Declaration - Electromagnetic Emmissions and Immunity for the Centerfield® 2

Guidance and manufacturer's declaration electromagnetic emissions  
IEC 60601-1-2: 2015, based to table 1

The OCULUS Centerfield® 2 is intended for operation in the electromagnetic environment specified below. The user of the Centerfield® 2 should ensure that it is being used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Centerfield® 2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
HF-emissions CISPR 11	Class B	
Harmonics emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	complies	

## Electromagnetic immunity, IEC 60601-1-2: 2015, based on table 4


Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV ± 15 kV	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30%.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

## Electromagnetic immunity, IEC 60601-1-2: 2015, based on table 5, 8

Electrical Fast transient/bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines	± 2 kV ----- ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 6100-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% $U_T$ ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree  0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single-phase: at 0 degree  0% $U_T$ ; 250/300 periods	0% $U_T$ ; 1/2 period at 0, 45, 90, 135, 180, 225, 270 and 315 degree  0% $U_T$ ; 1 period and 70% $U_T$ ; 25/30 periods Single-phase: at 0 degree  0% $U_T$ ; 250/300 periods	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Centerfield® 2 requires continued operation during power mains interruptions, it is recommended that the Centerfield® 2 be powered from an uninterruptible power supply or battery.

Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

Electromagnetic immunity, IEC 60601-1-2: 2015, based on table 4, 5

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidelines
Conducted RF IEC 61000-4-6	<p>3 V<sub>eff</sub> 150 KHz to 80 Mhz</p> <p>6 V in ISM- and amateur radio frequency bands between 150 kHz and 80 MHz</p> <p>80% AM to 1 kHz</p>	V <sub>eff</sub> = 3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of Centerfield® 2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3,5}{(V_1)} \right] \sqrt{P}$ $d = \left[ \frac{3,5}{(E_1)} \right] \sqrt{P} \quad 80\text{MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	<p>3 V/m 80 MHz to 2,7 GHz</p> <p>80% AM at 1 kHz</p>		$d = \left[ \frac{7}{(E_1)} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p>Interface may occur in the vicinity of equipment marked with the following symbol:</p> 
Note 1:	At 80 Hz and 800 MHz, the higher frequency range applies.		
Note 2:	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
<p>a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Centerfield® 2 is used exceeds the applicable RF compliance level above, the Centerfield® 2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Centerfield® 2.</p> <p>b. Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the Centerfield® 2, IEC 60601-1-2:2007, table 6

The Centerfield® 2 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Centerfield® 2 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Centerfield® 2 as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 Mhz $d= 1.2 \sqrt{P}$	80 MHz to 800 MHz $d= 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d= 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.80	3.80	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### 19.3 Description of the Connection



## 19.4 Data Sheet GSM60B15-P1J (05150725)



60W AC-DC High Reliability Medical Adaptor

**GSM60B series**


### ■ Features

- Universal AC input / Full range
- 2 pole AC inlet IEC320-C8
- Medical safety approved (2 x MOPP between primary to secondary)
- Suitable for BF application with appropriate system consideration
- Low leakage current <50uA
- No load power consumption<0.1W
- Energy efficiency level VI(Except 5~9V for Level V)
- Comply with EISA 2007/DoE,NRCan, AU/NZ MEPS, EU ErP and meet CoC Version 5
- Built-in active PFC function
- High efficiency up to 91.5%
- Fanless design with -30~+60°C working temperature
- Class II power (without earth pin)
- Protections: Short circuit / Overload / Over voltage / Over temperature
- Fully enclosed plastic case
- LED indicator for power on
- 100% full load burn-in test
- Optional lock type DC plug
- 3 years warranty

### ■ Applications


- Mobile clinical workstation
- Oral irrigator
- Portable hemodialysis machine
- Breath Machine
- Medical computer monitor

### ■ Description

GSM60B is a highly reliable, 60W desktop style single-output green medical adaptor series. This product is equipped with a 2-pin (no FG) standard IEC320-C8 power plug, adopting the input range from 80VAC to 264VAC. The entire series supplies different output voltages between 5VDC and 48VDC that can satisfy the demands for various kinds of medical electrical devices. The circuitry design meets the international medical standards (2\*MOPP), having an ultra low leakage current (<50 uA), fitting the medical devices in direct electrical contact with the patients.

With the efficiency up to 91.5% and the extremely low no-load power consumption below 0.1W, GSM60B is compliant with USA EISA 2007/DoE, Canada NRCan, Australia and New Zealand MEPS, EU ErP, and meet Code of Conduct (CoC) Version 5. The supreme feature allows the adaptor to save the energy when it is either under the operating mode or the standby mode. The entire series utilizes the 94V-0 flame retardant plastic case, providing the double insulation that effectively prevents electrical shock. GSM60B is approved with the international medical safety certificates.

### ■ Model Encoding

**GSM60B 05 - P1J**


File Name: GSM60B-SPEC 2016-03-16



60W AC-DC High Reliability Medical Adaptor

**GSM60B** series

**SPECIFICATION**

ORDER NO.		GSM60B05-P1J	GSM60B07-P1J	GSM60B09-P1J	GSM60B12-P1J	GSM60B15-P1J	GSM60B18-P1J	GSM60B24-P1J	GSM60B48-P1J	
OUTPUT	SAFETY MODEL NO.	GSM60B05	GSM60B07	GSM60B09	GSM60B12	GSM60B15	GSM60B18	GSM60B24	GSM60B48	
	DC VOLTAGE <small>Note.2</small>	5V	7.5V	9V	12V	15V	18V	24V	48V	
	RATED CURRENT	6A	6A	6A	5A	4A	3.33A	2.5A	1.25A	
	CURRENT RANGE	0 ~ 6A	0 ~ 6A	0 ~ 6A	0 ~ 5A	0 ~ 4A	0 ~ 3.33A	0 ~ 2.5A	0 ~ 1.25A	
	RATED POWER (max.)	30W	45W	54W	60W	60W	60W	60W	60W	
	RIPPLE & NOISE (max.) <small>Note.3</small>	100mVp-p	100mVp-p	100mVp-p	100mVp-p	100mVp-p	150mVp-p	180mVp-p	240mVp-p	
	VOLTAGE TOLERANCE <small>Note.4</small>	± 5.0%	± 5.0%	± 5.0%	± 3.0%	± 3.0%	± 3.0%	± 3.0%	± 2.5%	
	LINE REGULATION <small>Note.5</small>	± 1.0%	± 1.0%	± 1.0%	± 1.0%	± 1.0%	± 1.0%	± 1.0%	± 1.0%	
	LOAD REGULATION	± 5.0%	± 5.0%	± 5.0%	± 3.0%	± 3.0%	± 3.0%	± 3.0%	± 2.5%	
	SETUP, RISE TIME <small>Note.6</small>	1000ms, 30ms / 230VAC 1500ms, 30ms / 115VAC at full load								
INPUT	HOLD UP TIME (Typ.)	50ms / 230VAC 15ms / 115VAC at full load								
	VOLTAGE RANGE <small>Note.7</small>	80 ~ 264VAC 120 ~ 370VDC								
	FREQUENCY RANGE	47 ~ 63Hz								
	EFFICIENCY (Typ.)	81.5%	86%	87.5%	88%	88.5%	89%	90%	91.5%	
	AC CURRENT (Typ.)	1.4A / 115VAC 1A / 230VAC								
	INRUSH CURRENT (Typ.)	30A / 115VAC 65A / 230VAC								
PROTECTION	LEAKAGE CURRENT(max.)	Touch current < 50µA/264VAC								
	OVERLOAD	105 ~ 160% rated output power Protection type : Hiccup mode, recovers automatically after fault condition is removed								
	OVER VOLTAGE	5.25 ~ 6.75V	7.88 ~ 10.13V	9.45 ~ 12.15V	12.6 ~ 16.2V	15.75 ~ 20.25V	18.9 ~ 24.3V	25.2 ~ 32.4V	50.4 ~ 64.8V	
	OVER TEMPERATURE	Shut down o/p voltage, re-power on to recover								
ENVIRONMENT	WORKING TEMP.	-30 ~ +60°C (Refer to "Derating Curve")								
	WORKING HUMIDITY	20% ~ 90% RH non-condensing								
	STORAGE TEMP., HUMIDITY	-40 ~ +85°C, 10 ~ 95% RH								
	TEMP. COEFFICIENT	± 0.03% / °C (0 ~ 40°C)								
SAFETY & EMC <small>(Note. 8)</small>	VIBRATION	10 ~ 500Hz, 2G 10min./1cycle, period for 60min. each along X, Y, Z axes								
	SAFETY STANDARDS	ANSI/AAMI ES60601-1 / ES60601-1-11, TUV EN60601-1 / 60601-1-11 approved								
	ISOLATION LEVEL	Primary-Secondary: 2xMOPP								
	WITHSTAND VOLTAGE	I/P-O/P:4KVAC								
	ISOLATION RESISTANCE	I/P-O/P:100M Ohms / 500VDC / 25°C / 70% RH								
OTHERS	EMC EMISSION	Compliance to EN55011(CISPR11) class B, EN61000-3-2,3, FCC PART 15 class B, CAN ICES-3(B)/NMB-3(B)								
	EMC IMMUNITY	Compliance to EN61000-4-2,3,4,5,6,8,11, EN55024, EN60601-1-2, EN61204-3 medical level, criteria A								
	MTBF	720K hrs min. MIL-HDBK-217F(25°C)								
CONNECTOR	DIMENSION	125*50*31.5mm (L*W*H)								
	PACKING	0.32Kg; 40pcs/13.8Kg/1.05CUFT								
	PLUG	See page 3 ; Other type available by customer requested								
NOTE	CABLE	See page 3 ; Other type available by customer requested								
		1. All parameters are specified at 230VAC input, rated load, 25°C 70% RH ambient. 2. DC voltage: The output voltage set at point measure by plug terminal & 50% load. 3. Ripple & noise are measured at 20MHz by using a 12" twisted pair terminated with a 0.1µf & 47µf capacitor. 4. Tolerance: includes set up tolerance, line regulation, load regulation. 5. Line regulation is measured from low line to high line at rated load. 6. Length of set up time is measured at first cold start. Turning ON/OFF the power supply may lead to increase of the set up time. 7. Derating may be needed under low input voltages. Please check the derating curve for more details. 8. The power supply is considered as an independent unit, but the final equipment still need to re-confirm that the whole system complies with the EMC directives. For guidance on how to perform these EMC tests, please refer to "EMI testing of component power supplies." (as available on <a href="http://www.meanwell.com">http://www.meanwell.com</a> )								

File Name:GSM60B-SPEC 2016-03-16



60W AC-DC High Reliability Medical Adaptor

**GSM60B series**

**Derating Curve**

Ambient Temperature (°C)	Load (%)
-30	100
0	100
10	100
20	100
30	100
40	100
50	75
60	50

**Static Characteristics**

Input Voltage (VAC) 60Hz	Load (%)
80	80
90	90
100	100
110	100
120	100
130	100
140	100
150	100
160	100
170	100
180	100
190	100
200	100
210	100
220	100
230	100
240	100
250	100
264	100

**Mechanical Specification** Case No. GSM60B Unit:mm

UL2464 16AWG 1000±50mm for 5 ~ 15V  
UL1185 16AWG 1500±50mm for 18 ~ 48V

C+  
ID 2.1 x OD 5.5  
Outside ⊖ ⊕ Inside

**Plug Assignment**

Standard plug: P1J

P1J	
P/N	OUTPUT
CENTER	+

Optional lock type plug: P2S  
SWITCHCRAFT S761K plug equivalent

**Installation Manual**

Please refer to : <http://www.meanwell.com/webnet/search/InstallationSearch.html>

File Name:GSM60B-SPEC 2016-03-16





## Manufacturer and Service Address

Headquarters:

OCULUS Optikgeräte GmbH  
Münchholzhäuser Straße 29 • 35582 Wetzlar • GERMANY  
Tel. +49 641 2005-0 • Fax +49 641 2005-295  
E-Mail: [export@oculus.de](mailto:export@oculus.de) • [www.oculus.de](http://www.oculus.de)

USA:

OCULUS, Inc.  
17721 59th Avenue NE • Arlington • WA 98223  
Tel. +1 425 670 9977 • Fax +1 425 670 0742  
Email: [sales@oculususa.com](mailto:sales@oculususa.com) • <http://www.oculususa.com>

G/56980/0000/E  
Charge:  
G/56980/0000/E 1219 Rev02